
# BORT Generation Schulterorthese

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






Generation

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**REF** 215 700

Illustrationen  
Illustrations

## Sprachen Languages

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	04
<b>EN</b>	English	Instructions for use .....	09
<b>FR</b>	Français	Mode d'emploi .....	13
<b>ES</b>	Español	Instrucciones de uso .....	18
<b>IT</b>	Italiano	Istruzioni per l'uso .....	23
<b>NL</b>	Nederlands	Gebruiksaanwijzing .....	28
<b>CS</b>	Čeština	Návod k použití .....	33
<b>ET</b>	Eesti	Kasutusjuhend .....	37
<b>PL</b>	Polski	Instrukcja użytkowania .....	41
<b>RO</b>	Română	Instrucțiuni de utilizare .....	46

 PDF: [ga.bort.com](http://ga.bort.com)

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined position.

### **Indications**

Pre-operative, post-operative, post-traumatic, e.g. in case of distortion, contusion, subcapital humerus fracture, signs of paralysis, inflammation, degeneration.

### **Contraindications**

Impaired lymphatic drainage including unclear soft tissue swelling, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

### **Application risks/Important notes**

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician

- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

## Fitting instructions

- 1** Open both Velcro fasteners (1 and 2) on the arm pouch. Then, open the quick-release fastener (3) on the hanging button (4) and thus release the strap system from the arm pouch.
- 2** Put the arm pouch down on a surface in front of your body with the opening upwards. Place your arm in the arm pouch and position your elbow in the semi-circular stitched down area. Then close the Velcro fasteners (1 and 2).
- 3** Reach backwards with your free arm whilst guiding the strap system upwards over the shoulders (according to the rucksack principle).
- 4** Then, hang the strap system in front of your body with the quick-release fastener (3) into the hanging button (4) on the arm pouch. The fastener engages with resistance.
- 5** Tighten the waist (1) and shoulder strap (2) as required for the desired positioning of the arm and for optimum positioning.
- 6** **Information for technicians:** The illustrations in the fitting instructions depict the fitting on the right-hand side. When fitting on the left-hand side, the stabilisation element has its course mirrored.

For initial fitting, the stabilisation element is fixed once-only on the back between the waist and shoulder strap by the technician in the medical store. Therefore, the stabilisation element should not be placed directly centrally on the spine, but slightly offset to one side. The element also remains on the brace when it is removed, this simplifies putting it on again.

**7 Optional accessories available:** REF 991115 Hand rest with a ball  
Position the training ball in the metacarpus area and fix it using the Velcro area in the arm pouch.

## Removal

The brace is only removed by opening the Velcro fasteners on the arm pouch and the quick-release fastener. The stabilisation element guarantees easy fitting when further used.

## Material composition

Polyamide (PA), polyurethane (PUR), polyester (PES), elastane (EL), cotton (CO)  
The sewn in textile label provides the precise material composition.

## Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  
 Do not iron  Do not clean chemically

Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not

cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

### **Useful life/Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 07.2023

**MD** Medical device |  Single patient – multiple use

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