# BORT Generation Rückenorthese

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Generation

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
**REF** 215 800

**Illustrationen**  
**Illustrations**


**1**


**6**


**2**


**7a**


**3**


**7b**


**4**

**8**

**5**

## Sprachen Languages

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	04
<b>EN</b>	English	Instructions for use .....	10
<b>FR</b>	Français	Mode d'emploi .....	15
<b>ES</b>	Español	Instrucciones de uso .....	20
<b>IT</b>	Italiano	Istruzioni per l'uso .....	25
<b>NL</b>	Nederlands	Gebruiksaanwijzing .....	30
<b>CS</b>	Čeština	Návod k použití .....	35
<b>ET</b>	Eesti	Kasutusjuhend .....	40
<b>PL</b>	Polski	Instrukcja użytkowania .....	45
<b>RO</b>	Română	Instrucțiuni de utilizare .....	51

 PDF: [ga.bort.com](http://ga.bort.com)

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a brace for partial relief due to lordosis reduction of the lumbar spine.

### **Indications**

Facet syndrome, lumbosciatica, intervertebral disc protrusions/prolapse, radicular irritation syndrome, degenerative lumbar change, condition following intervertebral disc surgery, stable fractures, spondylosis, spondylolisthesis and spinal canal stenosis.

### **Contraindications**

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

### **Application risks/Important notes**

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician

- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

## Fitting instructions

### Instructions for use for patients:

- 1** When you receive your BORT Generation Back brace, the latch plate is closed. Please open it fully. Also, open both pulley handles and attach them to the side velour surfaces using Velcro.
- 2** Place the back brace around the torso so that the pelvic girdle with its curved protrusion is positioned above the iliac crest. This protrusion is central, on the inside of the pelvic girdle. Moreover, information on the inside of the back brace indicates the direction (**top = towards the head**).
- 3** Close the BORT Generation Back brace on the latch plate applying very little force, using the attached hand straps.
- 4** To close the BORT Generation Back brace fully, grasp the pulley handles and pull them both forwards simultaneously, whereby you only need to apply very little force.
- 5** Following this, attach the pulley handles to the latch plate or the velour surfaces at the side using the micro Velcro provided. In order to avoid strangulation, the brace should not be closed too tightly. Within the course

of this fitting process, the BORT Generation Back brace adapts fully to the anatomy of the body and interlocks.

**6** If the lower edge of the brace ends with the hip-bend when sitting, the brace is correctly positioned.

### **Removal**

To remove the brace, loosen the pulley handles off the velour and then open the latch plate.

### **Instructions for the technician:**

**7a** If required, you can move the pelvic girdle outwards and inwards by loosening the 4 Allen screws.

**7b** The option of placing the pelvic girdle in a tapered version in which the upper arm is positioned further inwards from the pelvic girdle and the lower arm further outwards from the pelvic girdle exists.

**Information:** Depiction without a support. 27 cm height

**8** The suspended elastic cords serve to position the brace in the pelvic region. They keep the brace interlocking on the body whilst settings are made by the technician.

**9** Finally, all 4 Allen screws should be tightened to a maximum torque of 1.5 Nm and the correct seat of the brace be checked again.

**10** If required, the elastic cords can be removed after adaptation.

## Info

The back plate is available with a height of 27 cm or 31 cm.  
A T Allen key with an across flat of 2.5 mm is provided.

## Removal

To remove the brace, loosen the pulley handles off the velour and then open the latch plate.

## Material composition

polyamide (PA), polyurethane (PU), cotton (CO), elastane (EL), polyester (PES)  
The sewn in textile label provides the precise material composition.

## Cleaning information

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer  
 Do not iron  Do not clean chemically

Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 02.2024

**MD** Medical device |  Single patient – multiple use

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