BORT Generation Rückenbandage

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Generation

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REF 215 850

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	03
EN	English	Instructions for use	07
FR	Français	Mode d'emploi	11
ES	Español	Instrucciones de uso	15
IT	Italiano	Istruzioni per l'uso	19
NL	Nederlands	Gebruiksaanwijzing	23
CS	Čeština	Návod k použití	27
ET	Eesti	Kasutusjuhend	31
PL	Polski	Instrukcja użytkowania	35
RO	Română	Instrucțiuni de utilizare	39

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace with a back pad and tension elements for stabilising the lumbar spine.

Indications

Degenerative lumbar changes, mild to moderate lumboischialgia, slipped disc with mild deficits, spondylosis, spondylolisthesis, spinal canal stenosis, radicular irritation syndrome, condition following intervertebral disc surgery.

Contraindications

Unclear soft tissue swellings, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used

- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

When you receive your BORT Generation Back Support, the front fastener is closed. Open it fully. The pad which is split into two in the back area must be positioned with the respective wider part of the pad in the direction of the head. Place both pads respectively on the outer edge of the support cover. Place the back support around the abdomen and close the BORT Generation Back Support applying very little force, using the attached handles. In order to close the BORT Generation Back Support fully, use the handles and pull them forwards simultaneously and evenly, applying very little force. Following this, attach the handles to the front fastener or the velour surfaces at the side using the micro Velcro provided. Observe here that the support is not closed too tightly. Within the course of this fitting process, the BORT Generation Back Support is fully fitted to the anatomy of the body and interlocks. If the lower part of the support ends with the hip-bend when sitting, the support is correctly positioned.

Removal

To remove the support, loosen the pulley handles and then open the front fastener.

Material composition

Polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

Cleaning information

 Hand-washing 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Close the Velcro fastening to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link:
www.bort.com/conformity

Status: 07.2023

 Medical device |  Single patient – multiple use

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