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
Stand: 05.2023

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a brace for active relief and correction of the lumbar spine and thoracic spine in the sagittal plane.

### **Indications**

Damages in the thoracic and lumbar spine region, e.g. in case of stable osteoporotic vertebral fractures and osteoporosis.

### **Contraindications**

Inflammatory or injury-caused skin changes (including swellings and redness) in the body region to be treated, circulatory problems, soft tissue swellings, neurologically caused sensory disorders and skin trophism in the body region treated (sensitivity disorders, sensory disorders with and without skin damages), ongoing immobilisation (contracture), in particular with older people.

### **Application risks/Important notes**

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material
- allow yourself time to become accustomed to the brace and increase the wearing time gradually
- if you have questions regarding the use and wearing time, please refer to your doctor or technician.

### **Fitting instructions: Instructions for use for the patient**

The brace should be worn directly over the clothes. Avoid easing the clothes in the vicinity of the brace, this may result in pressure sores.

- 1** Open the front, wide brace fastener.  
Guide the arms through the shoulder loops as accustomed for a rucksack. Then close the front, wide fastener on the stomach under slight tension. Ensure a good fit, the ventral plate should be positioned a little above the groin.
- 2** Straighten yourself up actively (“make yourself big”) and tighten both shoulder straps in sync by pulling them lightly.
- 3** If you find the correction unpleasant, you can reduce it by tilting the shoulder strap’s guide loops downwards and upwards with your thumbs as for a rucksack.

### **Removal**

To remove the brace, release the shoulder straps as described above, open the front stomach fastener and take the brace off like a rucksack.

### **Information for the technician:**

**Information:** The brace may only be adapted by expert staff!

Adaption of the brace on the patient:

**4/4.1** Then, shape the back plate to the contours of the spine. Ensure finger width distance (approx. 1.5 cm) of the course of the splint in the lumbar spine region and a little space between the back plate and the body on the

upper (cranial) edge to the body (approx. 2 cm). This supports the straightening of the patient and lordosis reduction of the lumbar spine.

**5/5.1** Place the brace on the patient, it can be positioned by the patient like a rucksack, then close the ventral plate. Ensure correct positioning of the brace, the lower edge must be at the height of the gluteal fold. The front ventral plate should be positioned at the height of the pubic bone and not restrict hip flexion.

**6** Then set the brace belt lengths. Remove the crocodile clasp ends of the belts and shorten the strap lengths accordingly. Start with the belt straps on the hip, always tighten the belts on both sides in sync. After completion, tighten the crocodile clasp ends and close the harness under slight tension. Pay attention to an appropriate length.

**7** Next, adjust the combined torso and shoulder straps and proceed as described above. Before finally shortening the straps, ensure that there is a gap of approx. 3 cm between the pull loops and the guide loops for the shoulder straps. Pad the shoulder straps in the armpit region using the enclosed strap pads.

**8** To reduce tension to the armpits and avoid strangulation, you can correct the course of the belt by changing the position of the dorsal guide loops (attached to the back plate with Velcro). The level of the costal arch serves as a reference point.

**9** In order to achieve ideal straightening of the patient, then use the pull loops in sync to tighten the torso and shoulder straps. Avoid overcorrection!

**10** To reduce tension or remove the brace, tilt the shoulder strap's guide loops downwards and upwards with your thumbs as for a rucksack.

Check the correct fit of the brace:

- Are the belt straps too tight or do they cut in?
- Is adequate space (approx. 2 cm) available at the upper (cranial) edge of the back plate so that the patient can sit up?
- Instruct the patient in putting the BORT Generation osteoporosis back brace on.

**11 Tip:** You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

### Material composition

Polyurethane (PU), polyamide (PA), cotton (CO), polyester (PES), elastane (EL), helanca (PA)  
The sewn in textile label provides the precise material composition.

### Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items.

**Information:** Before cleaning, pay attention to remove the back plate in order to avoid damages to the brace and your washing machine or tumble dryer.

### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 05.2023

 Medical device |  Single patient – multiple use

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