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

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for active relief and correction of the lumbar spine and thoracic spine in the sagittal plane.

Indications

Damages in the thoracic and lumbar spine region, e.g. in case of stable osteoporotic vertebral fractures and osteoporosis.

Contraindications

Inflammatory or injury-caused skin changes (including swellings and redness) in the body region to be treated, circulatory problems, soft tissue swellings, neurologically caused sensory disorders and skin trophism in the body region treated (sensitivity disorders, sensory disorders with and without skin damages), ongoing immobilisation (contracture), in particular with older people.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions: Information for the patient

The brace should be worn directly over the clothes. Avoid easing the clothes in the vicinity of the brace, this may result in pressure sores.

- 1** Open the front, wide brace fastener. Guide the arms through the shoulder loops as accustomed for a rucksack. Then close the front, wide fastener on the stomach under slight tension. Ensure a good fit, the ventral plate should be positioned a little above the groin.
- 2** Straighten yourself up actively (“make yourself big”) and tighten both shoulder straps in sync by pulling them lightly.
- 3** If you find the correction unpleasant, you can reduce it by tilting the shoulder strap’s guide loops downwards and upwards with your thumbs as for a rucksack.

Removal

To remove the brace, release the shoulder straps as described above, open the front stomach fastener and take the brace off like a rucksack.

Information for the technician:

Important note on use: The brace may only be adapted by expert staff.

The brace can be set with **TWO** different degrees of stability.

The following adaption describes the more stable/firmer brace setting (delivery status)

- 4** To adapt the back plate, remove the Velcro protector on the lower inner side of the brace. Then, remove the back plate from the splint bag with the attached reinforcement module.
- 5** To shape the back plate slightly, detach the clasps which connect the reinforcement module to the back plate. For this, press the side lever on the grey clasp downwards, spread it apart and push the clasp away outwards.

6 Then, shape the back plate to the contours of the spine. Ensure finger width distance (approx. 1.5 cm) of the course of the splint in the lumbar spine region and a little space between the back plate and the body on the upper (cranial) edge to the body (approx. 2 cm). This supports the straightening of the patient and lordosis reduction of the lumbar spine.

7 Then, shape the reinforcement module to the contours of the back plate just adapted. Pay attention to the position markings which define the reinforcement plate region. Pay attention to good form fit for both modules. Then connect the back plate to the reinforcement module using the 4 clasps previously removed. They are spread apart and respectively attached to the uppermost and lowest reinforcement plate openings. The clasp is properly closed when you hear a click. Finally, check the connection. Then, push the back element back into the brace and tighten the Velcro protector.

8 The guide can be configured dynamically or as stable in the shoulder height region. To realise a stable/non-elastic setting, the overlay which can be attached using Velcro must bridge the elastic restraint segment in the shoulder height region. For secure fixing of the overlay, pay attention to its guidance through the rucksack loop and fix the end on the base material.

9 Place the brace on the patient, it can be positioned by the patient like a rucksack, then close the ventral plate. The lower edge must be at the height of the gluteal fold. The front ventral plate should be positioned at the height of the pubic bone and not restrict hip flexion.

10/10.1 Then set the brace strap lengths. Remove the crocodile clasp ends of the strap and shorten the strap lengths accordingly. Start with the harnesses on the hip, always tighten the straps on both sides in sync. After completion, tighten the crocodile clasp ends and close the harness under slight tension. Pay attention to an appropriate length.

11 Next, adjust the combined torso and shoulder straps and proceed as described above. Before finally shortening the straps, ensure that there is a gap of

approx. 3 cm between the pull loops and the guide loops for the shoulder straps. Pad the shoulder straps in the armpit region using the enclosed strap pads.

12 To reduce tension to the armpits and avoid strangulation, you can correct the course of the strap by changing the position of the dorsal guide loops (attached to the back plate with Velcro). The level of the costal arch serves as a reference point.

13 In order to achieve ideal straightening of the patient, then use the pull loops in sync to tighten the torso and shoulder straps. Avoid overcorrection. To reduce tension or remove the brace, tilt the shoulder strap's guide loops downwards and upwards with your thumbs as for a rucksack [Image 3]. Check the correct fit of the brace: Are the harnesses too tight or do they cut in?

Dynamic setting of the brace

Step 1 [Image 4] Remove the Velcro protector on the lower inner side of the brace. Then, remove the back plate from the splint bag with the attached reinforcement module.

Step 2 [Image 5] Detach the connecting clasps which connect the reinforcement module to the back plate. For this, press the side lever on the grey clasp downwards, spread it apart and push the clasp away outwards.

Step 3 [Image 4] Then, push the back element back into the brace and tighten the Velcro protector.

Step 4 [Image 14] To realise dynamic reclination, the overlay which can be attached using Velcro in the shoulder height region must be removed.

Information: Store the reinforcement plate and overlay carefully in their packaging. During later care phases, it may be necessary to again make more stable settings.

15 Tip: You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

Material composition

Polyamide (PA), polyurethane (PU), cotton (CO), polyester (PES), elastane (EL), helanca (PA)
The sewn in textile label provides the precise material composition.

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items.

Information: Before cleaning, pay attention to remove the back plate in order to avoid damages to the brace and your washing machine or tumble dryer.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link:
www.bort.com/conformity

Status: 10.2024

 Medical device |  Single patient – multiple use

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