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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a support for compression of the Achilles tendon with local pressure pads consisting of flexible material. They are intended to have a massaging effect on the Achilles tendon.

Indications

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, for example in the case of Achilles tendon rupture, bursitis subachillea.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product

- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

- 1** Before putting it on, open the Velcro fasteners.
- 2** Grasp the upper edge of the BORT AchilloStabil Plus ankle support and pull it over the foot until the crescent-shaped pads are placed in the area between the Achilles tendon and the medial or lateral malleolus.
- 3** Place your foot on the ground and close first the upper and then the lower Velcro fasteners.

For additional relief of the Achilles tendon, one of the heel cushions provided can be inserted into the shoe on the side of the ailing leg. In order to avoid leg length differences between the healthy and ailing sides, place the second heel cushion provided in the shoe on the side of the healthy leg.

Removal

Open the three fasteners on the lower leg of the ankle joint support. Grasp the support on the lower edge and pull it up over the foot.

Material composition

Polyamide (PA), elastodiene/latex (LA), elastane (EL)

The sewn in textile label provides the precise material composition.

 The product contains latex and can trigger allergic reactions.

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.
Close the Velcro fastening to avoid damages to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

 Medical device |  Single patient – multiple use

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