
BORT select® ManuZip volar

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select®

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	03
EN	English	Instructions for use	07
FR	Français	Mode d'emploi	10
ES	Español	Instrucciones de uso	14
IT	Italiano	Istruzioni per l'uso	18
NL	Nederlands	Gebruiksaanwijzing	22
CS	Čeština	Návod k použití	26
ET	Eesti	Kasutusjuhend	29
PL	Polski	Instrukcja użytkowania	32
RO	Română	Instrucțiunile de utilizare	36

PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a wrist compression support for the compression of soft tissue with a local pressure pad and a reinforcement splint.

Indications

Chronic, post-traumatic or post-operative irritations in the wrist area, e.g. distortion, carpal arthrosis, tenosynovitis.

Contraindications

Circulatory disorders with various origins, lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Before putting on, bend the aluminium splint to the desired shape. Open the zip and place the support on the metacarpophalangeal joint affected and guide the thumb through the intended thumb opening at the same time. The aluminium splint runs along the surface of the palm. Finally, close the zip. The Stabilo fleece strap included can be fitted individually dosed. For this, attach the Stabilo fleece strap to the velour on the support intended for this purpose using the narrow Velcro strip. The Stabilo fleece strap can then be twisted around the wrist circumferentially.

Removal

To remove the support, please open the zip and loosen the Stabilo fleece strap. Grasp the splint in the palm area and pull the support off over the fingers. Close the Velcro fastener in order to avoid early wear to the Velcro surface.

Material composition

Polyamide (PA), elastodiene/latex (LA), viscose (CV)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer

Do not iron Do not clean chemically

Remove the aluminium splint before washing. Close the Velcro fastener in order to avoid damaging other items of clothing. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty

claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

MD Medical device |  Single patient – multiple use

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