

BORT TaloXpress

Ankle joint brace for stabilisation of the upper and lower ankle joint



- ⊕ Medial, anatomically pre-moulded brace shell and strap system for a very high degree of stabilisation of the upper and lower ankle joint
- ⊕ Supination limitation and therefore relief of the ligaments
- ⊕ Very low self-weight of the brace

BORT TaloXpress
No. 053 900

NEW

Med
In the midst
of healthy
life

4 FUNCTIONAL STRAPS

- For effective limitation of the talus prolapse and avoidance of supination movements

ANATOMICALLY PRE-MOULDED BRACE SHELL

- Special medial side profile design: ensures high material stiffness and stabilises the joint with pronation

SLIM BRACE DESIGN

- Modern and functional with a very low self-weight

MULTIFUNCTIONAL SAFETY

The combination of a medial brace shell and functional restraint guidance ensures effective stabilisation of the ankle joint. The interaction of the functional elements (brace shell and restraint) supports the ligaments and the ankle joint optimally.

BORT. Benefit at your side.

BORT TaloXpress

Ankle joint brace for stabilisation of the upper and lower ankle joint



Correct attachment - decisive for the functional stabilisation of the ankle joint



The foot is positioned in the fully opened brace. Here, the plastic shell must run along the inner side of the leg and sit flush to the foot and lower leg. The lower leg should be positioned vertical to the ground.



Closure of restraints 1 to 3 restricts the talus prolapse and anchors the brace safely on the foot. The medial brace shell has a supporting and stabilising effect. Plantar and dorsal flexion remain unrestricted, thus maintaining the natural motion sequence.



The 4th restraint (the pronation restraint) is used to set the degree of stability individually and indication-based. The positioning of the ankle joint in pronation position provides effective relief of the affected capsule/ligament.



Profile

BORT TaloXpress
No. 053 900

NEW



Med

PRODUCT-BENEFIT

- Ankle joint brace for stabilisation of the upper and lower ankle joint
- Stabilisation degree adjustable using restraints
- Very thin and slim construction
- Can be worn in shoes without any problems
- Effective stabilisation of the ankle joint and very comfortable to wear

COLOUR: ■ grey

THE MED-BENEFIT

- Can be fully opened, therefore easy to put on and simplest handling
- High stabilisation combined with low self-weight thanks to ultra-modern material
- Lightweight and slim structure for good wearing comfort

h-measure	SIZES
27 – 31 cm	small
31 – 34 cm	medium
34 – 40 cm	large

Please indicate left or right
1 cm = 0,3937 inch

MEASURING POINT



INDICATIONS

Conservative treatment of ankle distortions and ligament ruptures, post-operative protection after ligament suture/reconstruction, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, supination/pronation trauma

BORT. Benefit at your side.

BORT GmbH P.O. Box 1330 | D-71367 Weinstadt | Phone +49 (7151) 99200-0
Fax +49 (7151) 99200-50 | www.bort.com | medical@bort.com