BORT X-Walker Achillo

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Med
Mitten im
gesunden
Leben

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

REF 101 320

Illustrationen
Figures


1


4

1A


5


2

6

3

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	08
FR	Français	Mode d'emploi	12
ES	Español	Instrucciones de uso	16
IT	Italiano	Istruzioni per l'uso	20
NL	Nederlands	Gebruiksaanwijzing	24
CS	Čeština	Návod k použití	28
ET	Eesti	Kasutusjuhend	32
PL	Polski	Instrukcja użytkowania	36
RO	Română	Instrucțiuni de utilizare	40

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a lower leg foot brace for complete immobilisation of the lower leg and foot region in a defined, adjustable position.

Indications

Achilles tendon rupture, pre-operative or post-operative/post-traumatic immobilisation of the lower leg and foot region, post-operative in case of Achilles tendon joint.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- the Walker is appropriate for a maximum body weight of 130 kg
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Before going through steps 1 – 6, pull one of the two understockings provided over your foot/lower leg as required. Pay attention to it sitting crease-free.

- 1** Then, open the circumferential Velcro fasteners and remove the front cover. In order to make it easier to put on later, we recommend twisting the Velcro straps downwards, thus easing application of the brace.

1A Remove the padding liner and position the wedges provided to set the therapeutically necessary pointed position in the footrest. One wedge segment is equivalent to an elevation of approx. 10 mm. Not all wedge segments must be used. For final positioning, the individual segments can be fixed using the attached adhesive points. The wedge height can be adapted to the course of therapy. Finally, position the padding liner back in the footrest.

2 Place your foot and lower leg in the Walker. Pay attention to good heel positioning, it should be placed securely in the very back of the heel rest. The padding liner is then closed on the foot and in the lower leg region using the Velcro dots, pay attention to it sitting crease-free.

Information for the technician:

Setting the foot length:

If required, the foot length can be adapted to the patient's circumstances by moving the front binding. For this purpose, the front binding is raised and pulled upwards gently to disengage it. Following this, it can be moved and be locked in the desired position by pressing it downwards gently.

3 Then, adapt the grey front cover by attaching it centrally to the padding liner again using Velcro. The valve at the upper end should be freely accessible. If the valve is not freely accessible, pull the liner towards the knee (upwards). For the optimum functioning of the X-Walker, the front cover in the foot and lower leg region must be positioned between the padding liner and the brace rest so that both brace rests overlap.

Information: Depending on the leg circumference, it may not be possible to let the brace rests overlap in the calf region.

Close the four circumferential Velcro fasteners in the order specified below:

1. Start with the fastener on the ankle joint.
2. Close the fastener close to the ankle on the lower leg.
3. Close the strap in the metatarsus/toe region.
4. Finally, the lower leg fastener close to the knee.

If necessary, re-tighten again.

Information: The brace rest is designed flexibly in the calf area. If necessary, the circumference of the calf can be adapted by setting the calf fastener on the reverse side of the brace rest.

4 + 5 Remove the valve cap and screw the vacuum pump provided onto the valve. **Pay attention to hand-tight fitting.** In order to draw the air out of the liner and generate the vacuum, now pump several times until the pump no longer jumps back. Remove the pump from the

valve and re-attach the valve cover. You should now feel a surface stabilisation in the lower leg and foot region. If necessary, re-tighten the fastening straps again.

Tip: To ensure a good fit, we recommend that you pay attention to even distribution of the filling, which consists of a large number of small balls. By fluffing up the padding liner before putting the brace on or smoothing out the filling with your hand, you can achieve even distribution of the filling and optimise its fit.

Accessories available:

REF 100 370 Understocking

REF 100 380 Protective cover against wetness

REF 100 390 Protective cover for night-time

REF 990 175 Heel wedge set for X-Walker Achillo

Removal

6 To remove the aid, take the valve cover off and reduce the vacuum by pressing the blue button on the inside of the valve. You hear the air flowing into the liner. Then, open the harnesses and the liner, remove the front cover and the brace from the leg.

Tip: Close all the Velcro fasteners in order to avoid damaging other items, e.g. to further clothing.

Material composition

35% polyester (PES), 35% cotton (CO), 30% PUR foam (PUR)

Cleaning information

 Hand-washing  Do not bleach  Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Hand wash the Walker and the liner with lukewarm water and mild detergents. Close the Velcro fasteners in order to avoid damages, e.g. to items of clothing. Air dry it and avoid heat.

Before cleaning the liner, ensure that the valve cap is screwed on.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the responsible national body in your country under the following link: www.bort.com/md-eu-contact

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 01.2023

Medical device |  Single patient – multiple use

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