BORT Leistenbruchband

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REF 104 800

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D104800|2023-10|007 ML | Rev. 01

DE BORT Leistenbruchband

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EN BORT Inguinal Hernia Belt

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is an inguinal hernia belt with optionally one or two pads for reducing appropriate inguinal hernia.

Indications

Reducible, single-sided or double-sided inguinal hernia.

Contraindications

Unclear soft tissue swellings, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes ⚠

- The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
 - remove it during longer periods of rest, e.g. when sleeping
 - in case of numbness, loosen the medical device and remove it if necessary
 - in case of persistent complaints, consult the physician or specialist retailer
 - use the medical device in accordance with therapeutic needs
 - only use other products simultaneously after consultation with your physician
 - do not make any changes to the product
 - do not wear it on open wounds
 - do not use in case of intolerance of one of the materials used
 - no re-use. This medical aid is intended for treating one patient only
 - whilst wearing the inguinal hernia belt, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

First, position the pad on the hernial ring. Then wrap the wide truss around the body at the level of iliac crest (between the waist and hip) and fasten it by hanging it on the eyelet. Finally, guide the thigh strap between the legs around the thigh and close it with the metal clasp sewn on to the strap in the hip region. The pads are attached with Velcro and can therefore be positioned individually. This enables optimum setting on the hernial ring.

Removal

To remove, please release the thigh strap eyelets. Finally, detach the wide truss from the hook.

Material composition

polyester (PES), cotton (CO), elastodiene / latex (LA), viscose (CV)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

☑️ Hand-washing 30 °C ☒ Do not bleach ☒ Do not dry in a tumble dryer ☒ Do not iron ☒ Do not clean chemically
Hand-wash the hernia truss and pads using mild detergents. Do not use fabric conditioner. Air-dry.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2023

Medical device | 🧑 Single patient – multiple use

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