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
EN BORT SellaFlex Sport Thumb Support

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a thumb saddle joint support for compression of soft tissue with additional reinforcement elements on the thumb.

Indications

Chronic, post-traumatic or post-operative irritations in the thumb saddle joint area, e.g. in case of ulnar and radial ligament instability, early stages of rhizarthrosis.

Contraindications

Circulatory disorders with various origins, lymph drainage disorders, including unclear soft tissue swellings distal to the aid positioned, sensory disturbances of the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Important: When you receive your thumb support, the supporting module for use on the right hand is inserted. The attachment Stabilo strap module ① is attached to the support in the supporting module ②.

Place the support on the metacarpophalangeal joint affected and guide the thumb through the intended thumb opening at the same time.

1 Stabilo strap with the easyClick system for additional stabilisation of the wrist: ① and ③ attachment module, ② supporting module, ④ adjustment module.

Adapt the stabilisation degree of the Stabilo strap by moving the adjustment module ④ individually.

2 Guide the Stabilo strap circumferentially around the wrist under slight tension. Close the strap by attaching the free attachment module ③ to the adjustment module ④. Pay attention to comfortable fitting of the Stabilo strap, neither too loose nor too tight.

3 Left hand thumb support: If the thumb support with the Stabilo strap is needed for use on the left hand, hang the Stabilo strap out of the support. Remove the supporting module ② from the velour bag and slide it into the velour bag on the opposite side. The supporting module should always be located on the lower side of the arm towards the palm.

If the Stabilo strap with the easyClick system is not needed, the supporting module ② can be taken from the velour bag.

Removal

To remove the strap hang the Stabilo strap with the easyClick system out on the upper side. Then grasp the support in the palm area and pull it off over the fingers.

Material composition

Polyamide (PA, COOLMAX®), polyamide (PA), elastane (EL, LYCRA®), silicone (SI)
The sewn in textile label provides the precise material composition.

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

MD Medical device |  Single patient – multiple use

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