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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### **Intended purpose**

This medical device is a thigh support for compression of soft tissue with diagonal and elastic constraints.

### **Indications**

Acute: Swollen contusions with haematoma, muscle injuries.

Chronic: Protection against new muscle lesions, reduction of the burden on damaged muscle structures (e.g. scarring after old injuries), following muscle injuries.

### **Contraindications**

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

### **Application risks/Important notes**

Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used

- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

## Fitting instructions

A circumferential strap system is provided with the thigh support and can be attached to the thigh support as required (Step 5).

**1** Open both cross restraints and let the restraints hang downwards. Here, you see the green velour inner surface with the directional arrows and marking strips as a so-called “pull stop”.

Put the upper thigh support on when standing with a slightly bent knee. Here, the function zone (F) should be at the height of the injured region.

**Tip:** the thigh support can be used for the knee-stretching (frontal) and knee-bending (rear) muscles.

**2** Then, pull the first cross restraint downwards to the marking strip (see arrow). Here, the wide, red marking strip indicates the maximum pull (which must not be exceeded).

**3** Then, attach the cross restraints diagonally on the upper velour using Velcro.

**4** Repeat the process with the second cross restraint. Here, the restraints cross as an X. It makes no difference which cross restraint you attach first.

**5** If required, the circumferential belt fasteners provided can be attached to the velour lugs with Velcro under tensile stress. Ensure that the longer restraints are attached above and the shorter restraints below.

## Removal

To remove, loosen the circumferential strap fastener if used. Open the cross fasteners and take hold of the knitted fabric on the support to pull it over your foot.

## Material composition

polyamide (PA), polyamide / COOLMAX® (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

## Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 05.2023

**MD** Medical device |  Single patient – multiple use

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