BORT OmoBasic geschlossene Form

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
REF 121 500

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D12150012023-071007 ML 1 Rev.01

DE BORT OmoBasic geschlossene Form

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EN BORT OmoBasic Enclosed Style

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined position.

Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of distortion, luxation, subcapital humerus fracture, rotator cuff injury, scapula fracture

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- if wearing at night is prescribed, avoid negative influence on the circulatory system
 - in case of numbness, loosen the medical device and remove it if necessary
 - in case of persistent complaints, consult the physician or specialist retailer
 - use the medical device in accordance with therapeutic needs
 - only use other products simultaneously after consultation with your physician
 - do not make any changes to the product
 - do not wear it on open wounds
 - do not use in case of intolerance of one of the materials used
 - no re-use. This medical aid is intended for treating one patient only
 - whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material.

Fitting instructions

- 1 Slip the tubular knitted fabric over the lower arm affected until the thumb and fingers can be seen in the designated opening in the hand region. The tubular knitted fabric should be placed as far as possible in the axilla-shoulder region and the shoulder joint should be covered by the knitted fabric.
- 2 Then place belt 1 with its origin in the shoulder area over the neck to the opposite shoulder, guide it around the wrist on the ventral side and attach it there using Velcro. Here, the desired elbow joint position of 90° should be taken into account.
- 3 Then, guide belt 2 with its origin in the hand area around the waist on the ventral side of the side of the arm affected and attach it there.

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4 Observe that the Velcro surface not required on belt 2 is equipped with the fluffy Velcro provided in order to avoid damaging your clothing.

5 Tip: You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

Removal

To remove the OmoBasic brace, open all the Velcro fastenings and then remove the brace.

Material composition

Cotton (BW), pure ester (PUR), polyamide (PA), elastodiene/latex (LA), polyester (PES)

The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Close the Velcro fastening to avoid damage to other laundry items.
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 07.2023

MD Medical device | Single patient – multiple use

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