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
REF 121 970/121 970 F

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D121970|2022-04|007 ML 1 Rev.01

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EN BORT OmoARS with 3-Point-Strap

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a shoulder joint brace for immobilisation of the shoulder joint in a defined abduction position.

Indications

Pre-operative, post-operative, post-traumatic, e.g. in case of impingement syndrome, rotator cuff tear, endoprosthesis implantation, subacromial defilée extension (according to Neer), bursitis subacromialis

Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments on our around it as they can destroy the material

Fitting instructions

Initial fitting and adaption is normally made in the operating theatre.

Fitting and removal by the patient should take place seated at a table, hindering detrimental shoulder movement.

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Tip: The Velcro fasteners which close the arm pouch on the front side of the cushion should be open. Also open the hip strap Velcro fastener on the front part of the brace cushion.

- 1 Place the operated arm in the arm pouch so that the hand lies on the ball.
- 2 Close the arm pouch. The elbow must be well-placed in the arm pouch.
- 3 Then guide the hip strap behind the back to the front and close the Velcro fastener.
- 4 For the optimal functional position of the OmoARS (15° abduction angle), position the cushion so that the front surface is parallel to the front of your body.
Tip: If you are facing a table, the front surface of the cushion must be parallel to the edge of the table.
- 5 Then guide the shoulder strap over the uninjured shoulder to the front.
- 6 Finally, close the Velcro fastener.

7 Tip: A helper can simplify the implementation of attaching the shoulder strap. Have the shoulder strap fully removed (open the Velcro fastenings). The helper can then help you to put the shoulder strap on like a vest. Then, proceed as described above. The shoulder strap is only attached to the basic cushion using Velcro.

8 Tip: You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

If an external rotation position of 30° is desired, the optionally available external rotation wedge, REF 121 960, must be attached between the basic cushion and the arm pouch using Velcro.

Removal

To remove the OmoARS with the 3-point strap shoulder brace, open all the Velcro fasteners, loosen the Velcro in the reverse order, then remove the brace.

Material composition

PUR foam (PUR), polyamide (PA), cotton (CO), polyester (PES)
The sewn in textile label provides the precise material composition

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer
 Do not iron  Do not clean chemically
Close the Velcro fastening to avoid damage to other laundry items.
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

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Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 04.2022

Medical device |  Single patient – multiple use

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