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**REF** 145 300

4 005862 179938

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D145300I2023-05I002 ML 1 Rev.01

### **DE** BORT ImmoXpress

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### **EN** BORT ImmoXpress

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### **Intended purpose**

This medical device is a brace for immobilisation of the knee joint in a defined position (bent / stretched).

#### **Indications**

Pre-operative, post-operative, post-traumatic, joint, muscle or ligament injuries

#### **Contraindications**

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### **Application risks / Important notes** ⚠

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness or pressure sores, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

#### **Fitting instructions**

The BORT ImmoXpress consists of 3 parts: A large centrepiece and 2 side parts attached with Velcro. Before fitting the brace, the integrated aluminium splints located in the splint bags in the middle and side sections must be adapted to the patient's existing body contours. The desired flexion position in the knee joint can be readjusted by forming the aluminium splints on the centre part of the brace. The side parts which can be freely positioned enable adaption to different circumference measurements. The side parts should overlap by about 2 cm if correctly adjusted.

Place the BORT ImmoXpress around the leg from behind so that the kneecap (patella) is enclosed within the recess from behind. Thread Velcro fasteners through plastic eyelets, guide them back and attach them with Velcro under

slight tension. Fit the BORT ImmoXpress so that despite immobilisation of the knee joint no congestion is caused in the calf and foot region.

#### **Removal**

To remove the device, loosen the harnesses and take the brace off the leg.

#### **Material composition**

Polyamide (PA), PUR foam (PUR)

The sewn in textile label provides the precise material composition.

#### **Cleaning information**

☑️ Delicate wash 30 °C ☒ Do not bleach (only wash with colour or mild detergents) ☒ Do not dry in a tumble dryer ☒ Do not iron

☒ Do not clean chemically

Do not use fabric conditioner. Close the Velcro fastener in order to avoid damaging other items of clothing and remove the aluminium splints. Stretch back into shape and dry in the fresh air.

#### **Guarantee**

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### **Useful life / Lifetime of the product**

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### **Duty of notification**

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### **Disposal**

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### **Declaration of conformity**

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 05.2023

**MD** Medical device | Single patient – multiple use

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