ActiveColor Sport Epi-Sponge

active
COLOR SPORT

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

REF 1460 SP

4 005862 175695

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D1460SP12024-101003 ML I Rev. 01

DE ActiveColor Sport Epi-Sponge

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EN ActiveColor Sport Epi brace

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for relief of the muscle origins on the ulnar and/or radial epicondyle.

Indications

Overload prophylaxis, load-induced irritation, tennis or golf arm.

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/important notes

- The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
 - remove it during longer periods of rest, e. g. when sleeping
 - in case of numbness, loosen the medical device and remove it if necessary
 - in case of persistent complaints, consult the physician or specialist retailer
 - use the medical device in accordance with therapeutic needs
 - only use other products simultaneously after consultation with your physician
 - do not make any changes to the product
 - do not wear it on open wounds
 - do not use in case of intolerance of one of the materials used
 - no re-use. This medical aid is intended for treating one patient only
 - whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Guide the hand through the neoprene support. For this, loosen the circumferentially running Velcro tape as applicable. Position the Epi brace on the lower arm approx. 2 finger's width away from the elbow joint on the muscle belly. Please only close the Velcro so tightly that no congestions are caused. You should only feel pressure above the pad on the muscle with the Epi brace if you make a fist or a grasping movement.

1 In case of a tennis arm (epicondylitis humeri radialis), the pad is positioned on the upper muscle belly.

2 In case of a golf arm (epicondylitis humeri ulnaris), the pad is positioned on the lower muscle belly.

Removal

To remove the Epi brace, loosen the Velcro tape, increase the circumference and slide the brace off in the direction of the hand.

Material composition

Polyamide (PA), chloroprene rubber (CR), polyurethane (PUR), elastane (EL), polyester (PES)
The sewn in textile label provides the precise material composition.

Cleaning information

Delicate wash 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

Medical device | Single patient – multiple use

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