ActiveColor Handgelenkschiene

active
COLOR

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REF 1570 CH

4 003862 153383

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D1570CHI2020-07/001 ML

DE ActiveColor Handgelenkschiene

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EN ActiveColor Wrist Splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a compression wrist support for supporting the wrist.

Indications

Tendinitis/tenosynovitis, carpal tunnel syndrome (CTS), irritations, carpal arthrosis

Contraindications

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/important notes

The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

When delivered, the palm splint is not yet placed. To place it correctly, proceed as follows:

Open the zip and hold the support so that you can read the label and the zip points in your direction.

Using the support for the left hand: push the splint into the left splint bag, with the bend up front. Please observe that the recess must lie flush in the palm. Finally, push the rear end of the flap over the splint.

Using the support for the right hand: Push the splint into the right splint bag, with the bend up front. Please observe that the recess must lie flush in the palm. Finally, push the rear end of the flap over the splint.

Move the palm splint to the opposite splint bag if it is to be worn on the other side.

Open the zip and place the support on the metacarpophalangeal joint affected and guide the thumb through the intended thumb opening at the same time. Then, close the zip. The splint should run along the lower arm/surface of the palm.

The Stabulo fleece strap included can be fitted individually dosed. For this, attach the Stabulo fleece strap to the velour on the support intended for this purpose using the narrow Velcro strip. The Stabulo fleece strap can then be twisted around the wrist circumferentially.

Removal

To remove the support, loosen the Stabulo fleece strap. Then, open the zip, grasp the splint in the palm area and pull the support off over the fingers.

Material composition

Polyamide (PA), elastodiene/latex (LA), viscose (CV)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

Delicate wash Do not clean chemically Do not bleach
 Do not dry in a tumble dryer Do not iron
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Remove the plastic splint before washing. Close the Velcro fastening to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 11.2019

Medical device | Single patient – multiple use

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