BORT DynamicFX Rückenorthese

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Med
Mitten im
gesunden
Leben

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

REF 180 300

Illustrationen
Figures


1


6


2


7


3


8


4


9

5

10

Sprachen Languages

DE	deutsch	Gebrauchsanweisung	04
EN	english	Instructions for use	09
FR	français	Mode d'emploi	14
ES	español	Instrucciones de uso	19
IT	italiano	Istruzioni per l'uso	24
NL	nederlands	Gebruiksaanwijzing	29
CS	český	Návod k použití	34
ET	eesti	Kasutusjuhend	39
PL	polski	Instrukcja użytkowania	44
RO	românesc	Instrucțiuni de utilizare	49

 PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a brace for active relief and correction of the lumbar spine and thoracic spine.

Indications

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, e.g. stable, osteoporotic vertebral fractures, osteoporosis, cautionary treatment of secondary kyphosis (stable fractures, tumours), muscular insufficiency, pre-operative and post-operative stabilisation.

Contraindications

Inflammatory or injury-caused skin changes (including swellings and redness) in the body region to be treated, circulatory problems, soft tissue swellings, neurologically caused sensory disorders and skin trophism in the body region treated (sensitivity disorders, sensory disorders with and without skin damages), ongoing immobilisation (contracture), in particular with older people.

Application risks/Important notes 

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- in case of numbness, loosen the medical device and remove it if necessary
- should you incur pressure sores, pain or other problems whilst wearing the brace, remove it and contact your technician
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material
- allow yourself time to become accustomed to the brace and increase the wearing time gradually
- if you have questions regarding the use and wearing time, please refer to your physician or technician

Fitting instructions

Instructions for use for the patient:

The brace should be worn directly over the clothes. Pay attention to avoid creasing the clothes in the vicinity of the brace as this can result in pressure sores.

1 The BORT DynamicFX back brace which is individually compiled for you and adapted by the orthopaedic technician is pulled over the shoulders like a rucksack. Here, the basic support is open. Please observe that the shoulder straps are not twisted.

2 Close the basic support using the highest tension possible. In order to enable closing the basic support with as little application of force as possible, please use the hand straps attached to the basic support's latch plates. When closing, observe that the basic support is attached above the groin and does not cause unpleasant pressure when sitting. Both closure components from the basic support must be attached to each other parallel using Velcro.

3 Then, pull both side restraints forwards evenly and attach the restraints to the basic support using Velcro.

4 Finally, attach the chest strap to your BORT DynamicFX back brace.

Special notes for the patient:

Please have the fit and function of your BORT DynamicFX Back brace checked regularly by your orthopaedic technician. In particular when starting therapy with the BORT DynamicFX back brace, it is particularly important to check the length of the shoulder straps and the base plate and to re-adjust them as necessary. It is particularly important that you practice putting the frame support brace on and removing it with your physician/physiotherapist/orthopaedic technician in order to gain security in handling your aid.

Wearing time:

The wearing time for your brace should be determined within a treatment plan by your treating physician, depending on the current indication and your personal situation. Please observe that the therapy with the BORT DynamicFX back brace is a long-term therapeutic measure. Sustainable therapeutic success will only come after a while.

Information for the technician:

Information: The BORT DynamicFX back brace may only be adapted by trained staff!

Determination of the waist circumference:

measure the waist circumference and select the appropriate size from the basic support sizes available.

Selection of the base plate:

The length of the spine is measured along its course via the spinous processes of the spine from the intergluteal fold to the head (cranially). The base plate is decisively selected according to the affected segment/s as well as the patient's individual body size! The base plate is available in 4 different lengths.

Function of the spring tension elements:

The dynamics of the brace can be selected individually. The spring strength determines the extent of straightening and therefore postural improvement. The following spring strengths are available: 30 Newton, 40 Newton, 50 Newton, 60 Newton and 80 Newton. Primarily the spring strengths of 80N to 60N are used in the post-operative application area, the treating physician determines the strength to be used. Conservatively treated patients with muscular insufficiency, osteoporosis and also in case of therapy for secondary kyphosis can profit from spring strengths from 30N to 50N, always in relation to the patient's individual constitution.

Strap lock:

If in both the conservative and post-operative application area (temporarily) dynamics are not applied via spring tension elements, the strap lock in both the individual straps' belt connectors can be used. The installation of the strap lock in all three tension straps is used by the BORT DynamicFX back brace to fix the lumbar and thoracic spine when the pantyhose, and shoulder straps are closed correctly.

Inserting/exchanging the spring tension elements:

All spring tension elements are colour-marked, dependent on their spring constant:

30 Newton ➔ yellow, 40 Newton ➔ white, 50 Newton ➔ blue, 60 Newton ➔ red, 80 Newton ➔ black

5 To exchange the spring elements, release the spring protection bag press stud and push it upwards.

6 Open the Allen screws on both strap connectors with the Allen key provided and exchange the spring tension elements as required. Then, close the Allen screw again.

7 Now, pull the spring protection bag downwards again and close the press stud again on the base plate. All 3 BORT DynamicFX back brace straps should be equipped with the same spring strength.

Putting the pantyhose on:

8 Before putting the basic support on, separate it from the base plate by releasing the Velcro fastener. Position the basic support with the lower edge of the back section at the height of the gluteal fold. The patient then closes the basic support at the front.

Adjustment of the base plate:

9 Then, adapt the base plate to the patient's spine, adapted to the body. The base plate is cold-mouldable without tools. Please observe optimum straightening of the spine for patients with a strong inclination stance. Then, attach the base plate to the basic support using Velcro. Ensure that the base plate is fully flush fitted to the base support at the lower edge.

Adjustment of the reclination and shoulder straps:

10 Shorten the reclination and shoulder straps in the following order:

Step 1: Shorten the middle strap to that the shoulder plate with the upper material edge ends approx. 5 to 10 cm beneath cervical vertebra C7.

Step 2: Next, shorten both shoulder straps evenly. To shorten the shoulder straps, please remove the crocodile clasps which are respectively attached to the ends of the straps. After shortening, re-attach the crocodile fasteners with Velcro and close the brace straps again.

Ensure that the spring elements are not tense when the body is reclined. If the patient assumes a neutral body stance, light dosage of the intensity of the tension of the shoulder straps increases wearing comfort and nonetheless retains the therapeutic function.

Adjustment and attachment of the chest strap:

To shorten the chest strap, please remove the crocodile clasps which are respectively attached to the ends of the straps. After shortening, re-attach the crocodile fasteners with Velcro and re-attach to the shoulder straps. Shorten the chest strap so that both shoulder straps don't run in the axilla area. Slight tension on the chest strap moves the course of the shoulder straps forwards, hindering friction under the axilla. The chest strap is easy to open and close using the quick-release fastener in the middle. If required, attach the shoulder pads provided below the shoulder straps.

Instruct the patients in putting it on and using the BORT DynamicFX back brace.

Removal

To remove the brace, release the shoulder straps as described above, open the front stomach fastener and take the brace off like a rucksack.

Material composition

Basic support: Polyamide (PA), elastane (EL)

Back plate/belt system: Polyamide (PA), polyurethane (PU)

The sewn in textile label provides the precise material composition.

Cleaning information

☞ Hand-washing ☒ Do not clean chemically ☞ Do not bleach ☒ Do not dry in a tumble dryer
☞ Do not iron

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastening to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 11.2019

Medical device | ☞ Single patient – multiple use

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