# BORT Generation Hallux-Valgus-Schiene

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


REF 215 930

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D215930I2022-04I007 ML I Rev. 01

4 005862 134449

### DE BORT Generation Hallux-Valgus-Schiene

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### EN BORT Generation Hallux valgus splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a hallux valgus correction brace for correction of an x-malposition of the large toe. An adjustable fastener enables stepless, individual setting of the correction pressure.

#### Indications

Hallux valgus after hallux surgery to ensure the success of treatment, post-operative capsule relief.

#### Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAVK), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/Important notes ▲

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system

- in case of numbness, loosen the medical device and remove it if necessary

- in case of persistent complaints, consult the physician or specialist retailer

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– use the medical device in accordance with therapeutic needs  
– only use other products simultaneously after consultation with your physician

– do not make any changes to the product

– do not wear it on open wounds

– do not use in case of intolerance of one of the materials used

– no re-use – this medical aid is intended for treating one patient

– whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

– The brace can be worn both during the daytime and at night

– The brace can be worn in footwear

#### Fitting instructions

1 Open the BORT Generation Hallux Valgus Splint fasteners fully. The relief pad is attached to the brace with Velcro and can be freely positioned as needed.

2 Select the stabilisation element and slide it fully into the velour side pocket intended for this purpose. Static and dynamic application possibilities, configuration-dependent: 1. Spiral spring (flexible), 2. Splint (rigid).

3 In the next step, pull the support over the large toe (hallux). Here, the relief pad should be precisely on the same level as the metatarsophalangeal joint.

4 Close the BORT Generation Hallux Valgus Splint over the instep of the foot with an even pull. Pay attention to pleasant wearing comfort.

5 Pull the laterally positioned corrective support towards the heel with a well-dosed pull, here a visible correction of the hallux valgus (large toe position discrepancy) should take place.

The wearing time can vary from patient to patient and is determined by the existing anatomical circumstances. Remove the brace at regular intervals and check the skin appearance for redness or pressure sores.

## Removal

To remove the BORT Generation Hallux Valgus splint, open the Velcro fastener of the brace casing and the lateral corrective tape.

**Tip:** Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

## Material composition

Polyamide (PA), elastane (EL).

The sewn in textile label provides the precise material composition.

## Cleaning information

Wash delicate 30 °C. Do not bleach. Do not dry in a tumble dryer. Do not iron. Do not clean chemically. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Remove the stabilisation element. Close the Velcro fasteners in order to avoid damaging other items of clothing.

## Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of

application risks, instructions and unauthorised modifications to the product.

## Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 04.2022

Medical device | Single patient – multiple use

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