# BORT Stack'sche Schiene

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**REF** 470 100

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D47010012023-071004 ML I Rev.01

### DE BORT Stack'sche Schiene

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### EN BORT Stack Finger Splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a finger splint for immobilisation of the DIP joint.

#### Indications

Fracture, luxation, ligament injuries, extensor tendon rupture.

#### Contraindications

Unclear soft tissue swelling, impaired sensation and circulation in the specific body region, skin diseases in the specific part of the body.

#### Application risks/important notes

- The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
- in case of numbness, remove the medical device if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only

#### Fitting instructions

Place your stretched finger on a flat surface (e.g. a table) and put the splint over the finger affected so that the nail can be seen in the semi-oval opening. The fingertip should be placed flush at the end of the brace. Please observe that the rear end of the brace must not reach beyond the PIP joint. To attach the splint, you can use a skin-friendly, self-adhesive fixing tape (not included in the scope of delivery) and wrap it circumferentially around the rear end of the splint. If the splint is too long, it can be processed using a funnel grinding machine at low rotational speed. (Only trained expert staff should carry out this processing.)

#### Removal

To remove the splint, remove the fastener used and pull the splint off cautiously.

#### Material composition

Beige: 100% polyethylene (PE)  
Transparent: 100% (ABS)

#### Cleaning information

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Do not use fabric conditioner.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 07.2023

Medical device | Single patient – multiple use

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