# BORT Valco Hallux-Valgus-Schiene

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
**REF** 930 010

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D93001012022-041004 ML I Rev.01

### DE BORT Valco Hallux-Valgus-Schiene

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94% Polyamid (PA), 3% Polyester (PES), 3% Elasthan (EL).

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### EN BORT Valco Hallux Valgus splint

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

#### Intended purpose

This medical device is a hallux valgus correction brace for correction of an x-malposition of the large toe. An adjustable fastener enables stepless, individual setting of the correction pressure.

#### Indications

Hallux valgus after hallux surgery to ensure the success of treatment, post-operative capsule relief.

#### Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated.

#### Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material
- the BORT Valco Hallux valgus splint is mostly worn at night

#### Fitting instructions

**Important:** The splint cannot be worn in shoes!

Before application: Loosen the closure on the large toe (no. 1). Open the fasteners on the dorsum of the foot completely (nos. 2 and 3)

**1** Slide the BORT Valco Hallux valgus splint between the large toe and the second toe. The padded cup surrounds the inside of the ball of the large toe and is flush with the tip of the ball of the large toe. The guide strap (no. 1) is closed on the metatarsophalangeal joint and can be individually adjusted as required using the serrated splint.

**2** Gently pull on the non-elastic strap on the outside of the foot (no. 2) to close it so that any adjustment on the metatarsophalangeal joint can be felt and seen.

**3** To hold the brace firmly in place, close the elastic strap on the inside of the foot (no. 3). If applicable, retighten the large toe pad (no. 1) and the non-elastic strap (no. 2).

#### Removal

To remove the Hallux valgus splint, please fully open the non-elastic restraint (no. 2) fully and loosen the guide strap (no. 1). Finally, remove the Hallux valgus splint in a forwards direction over the large toe.

**Tip:** Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

#### Material composition

94% polyamide (PA), 3% polyester (PES), 3% elastane (EL)

#### Cleaning information

Hand-washing Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

#### Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 04.2022

Medical device | Single patient – multiple use

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