BORT Fersentlastungsschuh

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REF 930 160

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D930160I2024-10I007 ML I Rev.01

DE BORT Fersentlastungsschuh

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EN BORT Heel Relief Shoe

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a therapeutic shoe/relief shoe for relieving the heel.

Indications

Post-operative, wounds/injuries in the back foot area and on the heel.

Contraindications

No contraindications are known if used correctly.

Application risks/Important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
 - you may not drive a vehicle whilst wearing the BORT heel relief shoe
 - if wearing at night is prescribed, avoid negative influence on the circulatory system
 - in case of numbness, loosen the medical device and remove it if necessary
 - in case of persistent complaints, consult the physician or specialist retailer
 - use the medical device in accordance with therapeutic needs
 - only use other products simultaneously after consultation with your physician
 - do not make any changes to the product
 - do not wear it on open wounds
 - do not use in case of intolerance of one of the materials used
 - no re-use. This medical aid is intended for treating one patient only
 - whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Important:

Protect the open injury on the foot with a bandage or similar. Always wear a non-slip shoe with the same heel height on the healthy foot. Do not wear the shoe directly on the skin.

1. Open all the Velcro straps on the shoe fully, then open the shoe uppers' flaps in the front area of the shoe and above the heel.
2. Check the shoe for foreign bodies and remove them before putting it on.
3. Put your foot into the opened shoe.
4. Close the shoe uppers' flap in order to avoid any strangulation. Then, the flap over the heel.
5. Finally, close the Velcro tapes. They can be shortened, if required.

Important: a relief must be noticeable in the heel area.

Removal

Open all the Velcro straps on the shoe fully, then open the shoe uppers' flaps in the front area of the shoe and above the heel. Remove your foot from the heel relief shoe.

TIP: Close the Velcro fastener in order to avoid early wear to the Velcro surface.

Material composition

Upper: 50% polyurethane (PUR), 35% polyester (PES), 15% polyamide (PA)
Sole: 100% ethylene vinyl acetate (PUR)

The sewn in textile label provides the precise material composition.

Cleaning information

Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer Do not iron Do not clean chemically
Do not use fabric conditioner.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

Medical device | Single patient – multiple use

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