BORT Verbandschuh offene Form

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REF 930 180

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D930180|2024-10|008 ML | Rev. 01

DE BORT Verbandschuh offene Form

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EN BORT Therapeutic Shoe, open design

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose
This medical device is a ready-made therapeutic shoe for holding the foot provided with a wound dressing and possible insoles.

Indications
Post-operative and other wounds, even when treated with extensive and possibly padded bandaging, where the healing process takes place over an extended period of time and weight can be placed partially on the foot whilst walking, even outdoors. An additional insole is not required.

Contraindications
No contraindications are known if used correctly.

Application risks / Important notes
This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
– remove the medical device prior to radiological examinations
– you must not drive any vehicle whilst wearing the BORT therapeutic shoe
– if wearing at night is prescribed, avoid negative influence on the circulatory system
– in case of numbness, loosen the medical device and remove it if necessary
– in case of persistent complaints, consult the physician or specialist retailer
– use the medical device in accordance with therapeutic needs
– only use other products simultaneously after consultation with your physician
– do not make any changes to the product
– do not wear it on open wounds
– do not use in case of intolerance of one of the materials used
– no re-use. This medical aid is intended for treating one patient only
– whilst wearing the therapeutic shoe, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material

Fitting instructions
Important: Protect the open injury on the foot with a bandage or similar. Always wear a non-slip shoe with the same heel height on the healthy foot. Do not wear the shoe directly on the skin.
– Fully open all the Velcro shoe straps, then open the shoe tongue.
– Check the shoe for foreign bodies and remove them before putting the shoe on.
– Place your foot into the opened shoe when putting it on.
– Pay attention to safe / correct positioning of the foot in the shoe, i.e. the heel must be positioned as far back as possible in the shoe.
– Close the wide tongue of the shoe on the opposite Velcro side first. Then, the Velcro straps. If required, the Velcro strap can be shortened individually

in the back foot area. It is not necessary to adhere to a specific sequence for closing the fasteners.

Removal
Fully open all the Velcro shoe straps, then open the shoe tongue. It is not necessary to adhere to a specific sequence for opening the fasteners. Remove your foot from the therapeutic shoe.
Tip: Close the Velcro fastener in order to avoid early wear to the Velcro surface.

Material composition
40 % polyester (PES), 25 % thermoplastic elastomers (TPE), 15 % nylon, 15 % ethylene vinyl acetate (EVA), 5 % polyurethane (PU)

Cleaning information
 Hand-washing 30 °C Do not bleach Do not dry in a tumble dryer
 Do not iron Do not clean chemically
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

Guarantee
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life / Lifetime of the product
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).
You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2024

Medical device | Single patient – multiple use

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