BORT AktiVen®

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EN

BORT AktiVen® Medical compression stockings/tights with a closed toe

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a pair of medical compression stockings for the legs.

Indications

Varicosis, primary and secondary varicosis, varicose veins during pregnancy, supporting sclerotisation therapy, following surgical vein operations, thromboembolism, thrombophlebitis (superficial) as well as condition after healed phlebitis, deep leg vein thrombosis, condition after thrombosis, postthrombotic syndrome, thrombosis prophylaxis for mobile patients, chronic venous insufficiency (CVI), stages I – III CVI according to Widmer or C1 – C6 according to CEAP, ulcer prevention, deep vein insufficiency, angiodynplasia, oedemata, lymphoedemata, oedemata during pregnancy, post-traumatic oedemata, post-operative oedemata, cyclically idiopathic oedemata, lipoedemata from stage II, lip-lymphoedemata, phlebo-lymphoedemata, congestion as a result of immobility (arthrogenic congestion syndrome, pareses and partial pareses of the extremity). Further indications: Post-burn condition, scar treatment.

Contraindications

Absolute contraindications:

Advanced peripheral arterial occlusive disease, decompensated heart failure, septic phlebitis, phlegmasia coerulea dolens.

Relative contraindications:

Distinct weeping dermatoses, intolerance of compression stocking material, severe paraesthesia of the extremities, advanced peripheral neuropathy (e.g. diabetes mellitus), primary chronic polyarthritis.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, remove the medical device
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient only
- whilst wearing the medical device, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

1 Pull the stocking over the lower arm and take hold of the inside of the heel.

2 Hold the heel firmly, pull the stockings forwards and turn them inside out up to the heel.

Tip: Use the BORT AktiVen® Special Gloves for the following steps.

3 Take hold of the overlapping material in the heel area and slip over your foot until the knitted heel fabric is over the heel.

Tip: For the stockings, carry out steps 1–5 first with one leg and then repeat for the second leg.

4 Pull the stockings up over the heel gradually and thus simultaneously back with the right side out.

5 Then slip the stockings up step by step towards the knee, using the palm of your hand to shape them to the leg. At the same time, never pull the upper edge!

Tip: Turn the gripping edge (if one exists) inside out before putting them on.

6 Then repeat steps 1 to 4 for the stockings with the second leg. Turn both sides back the right way out and

slip them up both legs. Then, adjust the positioning of the fabric with the palm of your hand.

7 Finally pull the body belt up to the waist whilst observing that the gusset is vertical in the crotch.

Attention: Your compression stockings must never be twisted. Put the stockings gradually in place from the bottom to the top.

These tips make them easier to put on:

- Put the compression stockings on in the morning, directly after getting up.
- Rough, unfiled fingernails, sharp-edged rings and toenails, cornea and faulty shoe lining can damage the knitted fabric.
- Please always use BORT AktiVen® Special Gloves for putting on the compression stockings when getting dressed.

Removal

To remove the compression stockings/tights, please pull them from the top downwards over the foot.

Material composition

Polyamide (PA), elastane (EL)

The sewn in textile label provides the precise material composition.

Cleaning information

Before washing, all fasteners (if existing) must be closed! Please always wash stockings with a grip top inside-out. Wash the compression stockings separately first time (the dye may run). Your compression stockings should ideally be hand-washed or washed using a delicate wash cycle (40 °C) daily. We recommend the use of BORT special detergents for BORT AktiVen® medical compression stockings. **Tip:** a laundry net provides additional protection for the knitted fabric. If hand-washing, rinse thoroughly, do not wring. You can reduce the drying time by placing the compression stockings on a thick terry towel after washing them, rolling them up firmly and firmly pressing out excess moisture. Then hang the stockings up loosely. Do not leave them around in the towel and neither dry them on a radiator nor in the sun. The compression stockings must not be cleaned chemically. Please observe the care instructions on the sewn in textile label.

 Delicate wash 40 °C  Do not bleach
 Dry in a tumble dryer  Iron
 Do not clean chemically

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Wearing time and useful life

If not prescribed differently by your physician, wear your compression stockings daily from the morning to the evening. As a result of daily wearing and washing, the necessary medical compression pressure and elasticity of the knitted fabric may be reduced. The recommended duration of use is a maximum of 6 months. It is also dependent on the daily strain and correct handling when putting them on and removing them. The date specification on the hourglass symbol on the packaging label incorporates the maximum duration of use of 6 months. Products made-to-measure are determined for immediate use and also designed for a maximum duration of use of 6 months. After this, in case of subsequent care, the body dimensions must be re-checked by expert medical staff. In case of significant change to relevant body dimensions, new measurement and new care may also already be founded prior to expiry of the six-month duration of use.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the responsible national body in your country under the following link:
www.bort.com/md-eu-contact.

EN

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link:

www.bort.com/conformity

Status: 03.2024

[MD] Medical device |  Single patient – multiple use

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