Hammerzeh-Korrektur-Pad
Hammer Toe Correction Pad / Pad de correction d'orteil marteau
Cuscinetto correttivo per dita a martello / Correctiepad voor hamerteen
Korekční polštářek na kladívkové prsty / Haamervarvast korrigeeriv padjand
REF 137 120 large
ab Schuhgröße 40
Shoe size over 40 / A partir de pointure 40 / Dal numero 40
Vanaf schoenmaat 40 / Od velikosti bot 40 / Alates jalanumbrist 40
1 Stück / piece / pièce / pezzo / stuk / kus / tükk

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EN PediSoft Hammer Toe Correction Pad

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a device for the correction of dislocation of one or more hammer toes.

Indications

Hammer toe dislocation.

Contraindications

Do not use in case of open wounds. In case of diabetes, inflammations or vascular diseases, consult a physician prior to use!

Application risks/Important notes

- The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
 - remove it during longer periods of rest, e.g. when sleeping
 - in case of numbness, loosen the medical device and remove it if necessary
 - in case of persistent complaints, consult the physician or specialist retailer
 - use the medical device in accordance with therapeutic needs
 - only use other products simultaneously after consultation with your physician
 - do not make any changes to the product
 - do not wear it on open wounds
 - do not use in case of intolerance of one of the materials used
 - no re-use – this medical aid is intended for treating one patient
 - whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Fix affected toes on the side covered with textile using the fleece tie, tighten it on the reverse side and attach it to the Velcro. Protruding fleece tie can be shortened. For an optimum fit, wear socks or tights over the hammer toe correction pad if possible.

Removal

Loosen the fleece tie on the lower side. Then, pull the pad away to the front over the toe/s.

Material composition

Sealing tape: 91 % polyamide (PA), 9 % elastane (EL)
Reinforcement: 100 % polyethylene (PE)
Micro Velcro: 100 % polypropylene (PP)

Cleaning information

- ☞ Hand-washing ☒ Do not clean chemically ☞ Do not bleach ☒ Do not dry in a tumble dryer ☞ Do not iron
- Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2019

Medical device | ☞ Single patient – multiple use

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