BORT OmoBasic CH

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
REF 121 850 CH

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D121850CHI2020-05I002 ML

DE BORT OmoBasic CH

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EN BORT OmoBasic CH

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is an open shoulder-arm support for immobilising the shoulder joint in a defined position.

Indications

Pre-operative, post-operative, post-traumatic, e.g. for rotator cuff tear, endoprosthesis implantation, distortion, simple humerus fractures, supraspinatus syndrome.

Contraindications

Lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid congestion of the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient
- whilst wearing the shoulder-arm support, please neither use any creams nor ointments locally in the region of the attached aid as they can destroy the material

Fitting instructions

- Guide the arm affected through the closed, sewn strap until it is positioned beneath the elbow joint. Ensure that the strap is not twisted. The sewn arm strap position should be positioned towards the body.
- Then guide the strap behind the back and let it fall downwards vertically over the healthy shoulder.
- Then attach the guide loop provided to the black velour pad on the shoulder with Velcro. The guide loop should point downwards.
- Place the end of the strap around the wrist of the arm affected and then guide the Velcro strip through the previously attached guide loop. Close the Velcro fastener. The upper Velcro strip can be moved to change the arm position.

5 When you have set the correct position, attach the wide Velcro tape provided with the label OmoBasic to the support. For this, attach one end of the wide Velcro tape to the black velour beneath the elbow and the other end to the black velour on the wrist using Velcro. The wide Velcro tape should run horizontally from the elbow to the wrist.

6 If required, the product "Orthosan thorax strap for OmoBasic, art. no. 121 860 CH", which is available separately, can be attached to the support. For this, guide the strap circumferentially beneath the chest and also bind the upper arm with it. Finally, close the Velcro fastener in order to attach the thorax strap.

7 **Tip:** You can find our video, which shows the attachment process, both on our homepage and via the product-specific QR code.

Removal

To remove the support, remove the thorax strap first if used. Remove the wide Velcro tape with the OmoBasic label. Open the lower Velcro fastener to open the wrist flap. Then remove the strap from the lower arm cautiously and pull the support off.

Material composition

OmoBasic CH: Cotton (BW), polyamide (PA), polyester (PES)
Thorax strap for OmoBasic: Cotton (BW), polyamide (PA), polyester (PES)
The sewn in textile label provides the precise material composition.

Cleaning information

Delicate wash Do not clean chemically Do not bleach
 Do not dry in a tumble dryer Do not iron
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Remove the guide loop before washing. Close the Velcro fastenings to avoid damage to other laundry items.

Guarantee


The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

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Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

You can find the current declaration of conformity under the following link:
www.bort.com/conformity

Status: 11.2019

Medical device | Single patient – multiple use

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