

Orthopaedics | Surgery | Sports medicine | Rehabilitation

Playing
is possible!



**BORT-Bandages:
Also for children**



BORT GmbH

P.O. Box 1330 · D-71367 Weinstadt
Phone **49 7151 99200-0
Fax **49 7151 99200-50
medical@bort.com · www.bort.com

No game is without risk.



**BORT Immob-Splint
with patella recess**
No. 145 000

Body size
80 - 145 cm

3-component splint to immobilize knee joint in 0° or 20° flexion. Medial and lateral splint can be freely positioned. Inner face lined with terry cloth.

Indications: for complete immobilization of knee joint.



BORT TaloStabil® Eco
No. 054 600



Circumference
of ankle:
17 - 19 cm

Three-dimensional form-knitted dual-tension bandage with two silicone pads for massage. Applicable left and right.

Indications: distortion, ligament instability, after plaster cast.



BORT TaloFX®
No. 054 800



Circumference
of ankle:
17 - 19 cm

Soft orthosis with DynamicForceControl technology (DFC) for effective stabilization of the ankle with bilateral, soft padded cups. Limitation of talus advance by means of crossed straps. Not bulky, can be worn in shoes.

Indications: conservative treatment of pronation and supination traumas.



**BORT Silicone length
difference compensation**
No. 930 200

Shoe size:
33 - 35
Heights:
3mm, 5mm or 10mm

Heel cushion made of silicone for compensation of leg length differences, for application in shoe. Heights: 3 mm, 5 mm and 10 mm.

Indications: Compensation of leg length differences.

In this game, the winner is certain:
your health!



BORT bandages: Also for children

When children play, injuries can occur. This is why BORT now offers bandages for thumbs, wrists, shoulders, knees and ankles in smaller sizes.

As with all our medical products, the “BORT Program **for Children**” was developed in close cooperation with medical practitioners, therapists, orthopedic technicians and designers. Based on this cooperation, we are able to constantly improve established products and develop new high-grade product innovations. BORT bandages **for children** comply with stringent quality checks to ensure that the little patients can concentrate on a greater goal: recover as soon as possible.



This mark and the text
“For children”, from now on, designates
BORT bandages that are in children’s sizes!

Look for this, because children’s injuries demand
optimal support.




BORT bandages for children

The natural urge for action sometimes leads children to hazards that adults have long forgotten. No difference whether you happen to be a Cowboy or an Indian: "Gee, this is sore!"



BORT Cervical support
No. 127 260



Cervical support made of soft foam material with recess for larynx and Velcro® closure for wearing at day and night.

Indications: whiplash injury of the cervical spine, aching wryneck, after operations.



BORT DigiSoft®
Finger orthosis
No. 112 770




Finger splint with malleable alu splint and Velcro® closures at metacarpus and fingers, which can be placed in variable positions.

Indications: injuries and surgeries, immobilization after plaster cast, redressement.



BORT SellaTex®
No. 112 050




Sturdy bandage with embraced thumb, volar and dorsal wrist splint as well as dorsal thumb splint made of breathable material.

Indications: after thumb fracture, ulnar/radial thumb ligament injury, protection against wrong movements, tendopathy.



BORT ManuStabil® short
No. 112 730




Sturdy bandage with three Velcro® fasteners, volar and dorsal splint and circular Stabilo® strap made of breathable material. **Can also be worn at night.**

Indications: after removing of cast, stabilizing wrists after non-bone injuries, radial nerve palsy (wrist-drop).

The small heroes now need big help: BORT bandages take the load off, protect, alleviate the pain and promote recovery. Why? Because the next playmate is already waiting.



BORT StabiloFix®
No. 104 600



Bandage used as support for **upright position** for posture correction and as **clavicle bandage** in case of clavicular injuries. Made of soft, skin-friendly material.

Indications: for postural disturbance, clavicle fracture.



BORT Patella Bandage
in case of
Osgood-Schlatter
No. 114 510




Knee bandage with silicone pad for pressure relief at the patellar tendon base, and one lateral splint each side.

Indications: Morbus-Osgood-Schlatter, articular effusion, swellings.



BORT StabiloGen® Eco
No. 114 520




Patella fixation bandage with silicone pad and one splint on each side.

Indications: Patella instability, articular effusions, swellings, post-operative irritations, ligament instability.



BORT Stabilo®
Knee bandage with
articulated joint
No. 150 140



Knee bandage with silicone pad, interlocking hinge with extension limit at 0° and circular fastening straps for thigh and shank.

Indications: lateral ligament instability of knee joint.