

Catalogue 2025

Orthopaedics



BORT. Benefit at your side.

BORT is your performant partner

Convincing with a comprehensive selection, high quality and many years' experience

We have maintained a trustworthy partnership with our specialist shop customers for more than 40 years, convincing with our quality, selection and experience.

A full range of medical devices from the orthopaedic and phlebology areas with provable therapeutic benefits and top manufacturing quality as well as a known good value for money makes BORT a performant partner.



ORTHOPAEDICS

Our range of orthopaedic aids is structured with various product lines and is attuned to the individual needs of various target groups.

Innovative product development and production, state-of-the-art logistics

New and innovative product solutions for competent patient care arise repeatedly from knowledge within the company, intense expert exchange with the doctors, therapists and our specialist trade partners.

BORT AktiVen*: You can also trust proven BORT quality in the area of compression therapy

BORT AktiVen® Medical compression stockings for phlebology are manufactured using the circular knitting technique whilst adhering to stringent quality standards and are also available made-to-measure



PHLEBOLOGY



BORT AktiVen® INDIVIDUAL: for fast, individual patient treatment made-to-measure

INDIVIDUAL benefit Each stocking is adapted to the individual leg dimensions

- Comfort benefit Soft haptics for top wearing comfort
- Perfect fit Made-to-measure for you



Therapeutic success with diversity, quality and competence.





BORT select* The first choice for demanding people.

SENSIL, SENSIL[®] is a registered trademark of NILIT.

BORT Sport: The sporty side of BORT.

COOLMAX[®]





BORT Generation: For better mobility for the aged.



BORT Med: In the midst of healthy life





5

BORT **select**^{*} reflects the pursuit of excellence in quality at BORT, where special product features with added value meet the highest user requirements. Here, exclusive design and perfect functionality meet with therapeutic effectiveness, providing wearing comfort at its finest.

Stabilisation and relief for moving moments: the Sport product line combines specific functional requirements with modern, sporty design. Attractive equipment features, the use of high-tech fibre and an innovative fastener system make the supports an ideal accompaniment for all kinds of bodily activity.

Especially easy to handle, particularly comfortable to wear: our special offer for people with restricted freedom of movement or gripping ability is called Generation. The products help to retain independence in everyday life and can considerably improve the quality of living.

From young to old, from the head to the foot - the Med product line encompasses our diverse range of supports and braces. The broad delivery programme enables the care of virtually all joint-specific indications in accordance with the therapy.

An overview of everything which is new

Year by year we actively develop new product solutions that meet current treatment needs.





For reliable immobilisation

Modular brace system, wrist stabilisation, downgradeable.

Page 67

Pain relief in case of overstress

Wrist brace with a medial window for immobilisation and stabilisation in case of Tendovaginitis de Quervain.

Page 77



For correction of the patella glide path

Sports knee support with a grip top, also available in Plus and Double Plus sizes.

Page 15, 145

For ruptures of the muscle fibre and muscular injuries

Lower leg active support with a proprioceptive effect.

Page 14, 142





OUR PROGRAMME

Product ranges	8
KIDS	27
Overview of indications	37
Neck Shoulder	38 41
51001021	41
Overview of indications Elbow	51 52
	52
Overview of indications	63
Fingers Thumb	64 68
Overview of indications Hand	81 82
Overview of indications Thoracic Spine	101 102
Back	103
Overview of indications	121
Torso	122
Groin Hip	128 131
Overview of indications Knee	133 134
Thigh/Calf	142
Overview of indications	157
Ankle	158
Foot	171
Shoes Insoles	172 175
	100
Reha Products BORT Special Detergent	180 183
Overview OTC Products	184
EasyLife® products for rehabilitation ActiveColor®	185 186
CLIMACare® warmers	186
PediSoft® foot care	198



The first choice for demanding people.

Select®

THE SELECT-BENEFIT

- SOFTflex Plus knit fabric in feel-good quality with cooling luxury SENSIL fibre
- State-of-the-art knitting technology for top wearing comfort and a perfect fit
- Attractive product design with high-quality packaging
- 🖲 Made in Germany



INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **epicondylitis** (acute and chronic), tendomyopathy

See page 52



INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

See page 82



INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

See page 134



INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

10 SELECT __





INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

See page 104



INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, lumbosciatica, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

SELECT 11



INDICATIONS

Moderate lumboischialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade I spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history), i.e. after inter-vertebral disc surgery with mild to moderate prolapse, facet joint arthrosis See page 103



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, tendomyopathies

See page 158



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

See page 159



INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea



The sporty side of BORT.



THE SPORT-BENEFIT

- SOFTflex 3D knitted fabric with a velvety soft surface structure
- 📀 For more elasticity, stability and very comfortable to wear
- High-tech fibre COOLMAX®
- Innovative: Stabilo[®] Strap with easyClick system
- Modern, sporty design with a reflecting label
- 📀 Made in Germany

BORT EpiBasic Sport REF 122 600 SP



INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **epicondylitis** (acute and chronic), tendomyopathy

BORT EpiContur[®] Sport REF 022 500 SP



INDICATIONS

Epicondylitis, chronic or recurrent epicondylopathy, tendopathy, overstraining caused by work and sport, 'tennis/golf elbow'

See page 52

SPORT 13

BORT SellaFlex Sport Thumb Support REF 112 030 SP



INDICATIONS

Chronic, post-traumatic or post-operative irritations in the area of saddle joint of the thumb, in the case of ulnar and radial ligament instability, early stages of rhizarthrosis BORT ManuBasic[®] Sport



COOLMAX

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

See page 68

See page 82



INDICATIONS

Post-operative, post-traumatic, or degenerative irritations in the area of saddle joint of the thumb, e.g arthrosis of the CMC-I joint (thumb saddle joint), post-operatively following arthroplasty of the CMC-I joint, instability of the CMC-I joint

See page 68

COOLMAX® is a trademark of The Lycra Company



BORT Corset Sport

NDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine, back pain (pain in the region of the thoracic spine)

14 SPORT ____



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), ischialgia, dorsalgia

BORT StabiloBasic Lady Sport Back Support with Pad REF 104 680 SP



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), ischialgia, dorsalgia

See page 106



INDICATIONS

For soft tissue compression on the lower leg in case of acute muscle fiber tears and/or acute muscle strains

See page 106



INDICATIONS

Preventive: For highly strained muscle groups, Acute: Swollen contusions with haematoma, muscle injury, Chronic: following muscle injury, protection against new muscle lesions, reduction of the burden on damaged muscle structures (e.g. scarring after old injuries)

SPORT **15**



INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy



INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

See page 136



INDICATIONS

Patellofemoral pain syndrome, high riding patella, patellar lateralisation, condition after lateral release surgery, recurrent or chronic, patellar instability, meniscopathy See page 136



INDICATIONS

Patellar chondropathy, femoropatellar arthrosis, patellar tendonitis

See page 145

See page 143

COOLMAX® is a trademark of The Lycra Company

16 SPORT _____



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

See page 159



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, arthrosis



INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, in the case of Haglund syndrome heel irritation, achillodynia, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

See page 160

For better mobility for the aged.

-10

DAS GENERATION-PLUS

- Easy handling for more independence
- Safe application due to clear marking
- Support without limits
- Comfortable to wear



INDICATIONS

Post-operative, post-traumatic, rhizarthrosis, degenerative, inflammatory, i. e. arthritis, ulnar/radial ligament injury of the thumb

BORT Generation Thumb-Wrist Brace REF 215 200



INDICATIONS

Pre-operative, post-traumatic, postoperative, i. e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, i. e. rhizarthrosis, in the case of degenerative diseases such as carpal tunnel syndrome



INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e.g. distorsion, carpal arthrosis, carpal tunnel syndrome (CTS), rheumatism, tenosynovitis

See page 87



INDICATIONS

Pre-operative, post-operative, posttraumatic, carpal arthrosis, inflammatory-rheumatic, degenerative diseases, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis



INDICATIONS

Pre-operative, post-operative, posttraumatic carpal arthrosis, inflammatory-rheumatic, degenerative diseases, e.g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

_ GENERATION 19

BORT Generation Schulterorthese REF 215 700



INDICATIONS

Pre-operative, post-operative, posttraumatic, e. g. sprain, contusion, **subcapital humerus fractures**, signs of paralysis, inflammatory, **degenerative**

See page 41



INDICATIONS

Tendency to fall for various reasons, danger of falling

See page 131

BORT Generation Dynamic Osteoporosis Brace REF 215 880



INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, damage in the thoracic and lumbar spine region, stable, osteoporotic vertebral fractures, osteoporosis BORT Generation Osteoporosis Brace REF 215 870



INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, damage in the thoracic and lumbar spine region, stable, osteoporotic vertebral fractures, osteoporosis

See page 118

20 GENERATION



INDICATIONS

Degenerative lumbar change, mild to moderate lumboischialgia, slipped disc with mild deficits, spondylosis, spondylolisthesis, spinal canal stenoses, radicular irritation syndrome, condition after inter-vertebral disc surgery



INDICATIONS

Facet syndrome, lumbosciatica, invertebral disc protrusions/prolapse, radicular irritation syndrome, degenerative lumbar change, condition after inter-vertebral disc surgery, stable fractures, spondylosis, spondylolisthesis, spinal canal stenoses

See page 113



INDICATIONS

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability See page 113



INDICATIONS

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability

_ GENERATION 21



INDICATIONS

Degenerative diseases, mild to moderate lateral ligament instability of the knee joint, **gonarthrosis**, meniscus injuries



INDICATIONS


Partial relief of the medial or lateral compartment in case of level II-III gonarthrosis, Bone bruise (posttraumatic bone marrow edema), post-operative care following tibial plateau fracture

See page 149



INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis See page 150



INDICATIONS

Hallux valgus, after hallux surgery to ensure the success of treatment, post-operative capsule relief

See page 177

In the midst of healthy life.



THE MED-BENEFIT

- Full cover of indications
- Very practical
- Comfortable to wear
- Easy to use



INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, post-operative



INDICATIONS

Severe painful conditions, minor whiplash injuries, cervical syndrome, post-operative

See page 39

See page 40



INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **epicondylitis** (acute and chronic), tendomyopathy



INDICATIONS

Post-operative, post-traumatic, degenerative, inflammatory conditions, such as in case of arthritis, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, rhizarthrosis

See page 53

See page 70

*Length information only denotes approximate measures and may vary by size.

24 MED ____



INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis BORT ManuStabil® Short II2 730

INDICATIONS

Pre-operative, post-operative, posttraumatic, e. g. sprain degenerative, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

See page

See page 83



INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

_ MED 25



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

See page 107



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

26 MED ____





INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis See page 163



INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea



INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament rupturesankle fractures (Weber A), chronic ligament instability, shortening of time spent in plaster, ligament injuries, distortion of upper ankle joint, fibular ligament rupture, pronation and/or supination traumas



New leeway: BORT orthopaedic aids



THE KIDS-RANGE

 BORT supports and braces for kids Best protection, so that small joints can grow healthy. Cids

28 BORT. Benefit at your side.



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- Soft, open-pored foam
- With recess for larynx
- Adjustable Velcro fastening



PRODUCT-BENEFIT

- Open shoulder/arm adduction brace for immobilisation of the shoulder joint
- Stabilising strapping system ensures secure fit and immobilisation of the arm in the adduction position
- Soft, tension-resistant straps

COLOUR: blue

THE KIDS-BENEFIT

- Can be worn day and night
- Particularly rounded edges
- Skin-friendly 100% cotton cover

COLOUR: blue



THE KIDS-BENEFIT

- Flexible hand and arm adjustment
- Easy access to the surgical area to allow for wound care
- Skin-friendly

INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, cervical syndrome, post-operative

INDICATIONS

Pre-operative, post-operative, posttraumatic, sprain, dislocation, subcapital humerus fractures, scapular fracture, rotator cuff injury

SIZES

Circumference of neck

Children (x-small)_up to 33 cm

Adult sizes see page 38

Chin heights 1 cm = 0.3937 inch 5,0 cm / 6,0 cm

Length Olecranon to SIZE metacarpal knuckle point V



right and left symmetrical 1 cm = 0.3937 inch

BORT Clavicle Support REF 121 700 KIDS



- Backpack support
- Soft padded straps
- Three Velcro fastenings, individually adjustable, ensure an ideal blockage

BORT StabiloFix® REF 104 600 KIDS



- Corset for correcting the posture of the thoracic spine
- Soft, non-elastic strap material

COLOUR: blue

THE KIDS-BENEFIT

- Conforms ideally to the body's contours, especially in the critical armpit area
- Comfortable to wear
- Skin-friendly
- Can be worn day and night

COLOUR: blue

THE KIDS-BENEFIT

- Fastening in the chest area makes it possible to put on/remove the item independently
- Conforms perfectly to the body contours, especially in the critical armpit area
- Fitted with tucks to enable straps to be shortened easily
- l atex-free

INDICATIONS

Treatment of clavicle fractures

INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine

SIZE	Underbust circumference	SIZE	Underbust circumference
Ĩ	Children up to 65 cm	Ĭ	Children 50 – 70 cm
<u>∖</u> } *+*{ {	Adult sizes see page 40	$= \int \int_{-\infty}^{\infty} \frac{1}{1+1} \int_{-\infty}^{\infty} \int_{-\infty}^{\infty} \frac{1}{1+1} \int_{-\infty}^{\infty} \int_{-\infty}^{\infty} \frac{1}{1+1} \int_$	Adult sizes see page 102
1 cm = 0.3	937 inch	1 cm = 0.3	

30 BORT. Benefit at your side.



PRODUCT-BENEFIT

- Finger brace for immobilisation of interphalangeal joints
- Adjustable aluminium splint with soft padding



PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal ioint and saddle ioint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Adjustable thumb element
- Length: approx. 17 cm (6.7 inches)

COLOUR: grey

THE KIDS-BENEFIT

- Immobilisation of either one or two fingers
- Dorsal fitting possible after remodelling
- Velcro fasteners, on metacarpals and fingers can be positioned as required

COLOUR: grey

THE KIDS-BENEFIT

- Slip-on aid
- Comfortable to wear
- Stabilo® Strap


INDICATIONS

Post-operative, post-traumatic, i. e. redressement, positioning, immobilisation of the interphalangeal joint

INDICATIONS

SIZE

Pre-operative, post-traumatic, postoperative, i. e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes



right and left symmetrical 1 cm = 0.3937 inch Children ____ 13 – 15 cm Adult sizes see page 79

Circumference of wrist

please indicate left or right 1 cm = 0.3937 inch

BORT ManuStabil® short REF 112 730 KIDS



PRODUCT-BENEFIT

- Rigid wrist brace with palm and dorsal reinforcement splints, individually shaped
- Short design
- Length: approx. 17 cm (6.7 inches)

COLOUR: 🔳 grey

THE KIDS-BENEFIT

- Stabilo[®] Strap for effective stabilisation of the wrist
- Breathable material provides optimal comfort
- Can also be worn at night

Osgood-Schlatter Knee Support REF 114 510 KIDS



PRODUCT-BENEFIT

- Knee support for local decompression of tibial tuberosity
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Anatomically contoured knit
- Multi-zone silicone pad
- One lateral spiral stay on each side
- Compression class II

COLOUR: blue

THE KIDS-BENEFIT

- Velvety soft surface structure, creating very pleasant wearing comfort
- Highly elastic bending zone, also protects the popliteal fossa when worn for longer periods
- Adjustable setting using V-shaped strap correction

INDICATIONS

Pre-operative, post-operative, posttraumatic, i. e. sprain, peripheral nerve paralysis, tenosynovitis


INDICATIONS

SIZE

Facilitation of the movement of the patella, **Osgood-Schlatter disease**, tendinosis of the tibial tuberosity ligament, articular effusion, irritations

d-measure f-measure

Children 23 - 26 cm 28 - 32 cm



please indicate left or right 1 cm = 0.3937 inch Adult sizes see page 141

right and left symmetrical 1 cm = 0.3937 inch

32 BORT. Benefit at your side.



PRODUCT-BENEFIT

- Classical knee support for compression of soft tissues with pad for stabilisation and relief of the knee joint
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- One lateral spiral stay on each side
- Viscoelastic pad
- Compression class II

COLOUR: blue

THE KIDS-BENEFIT

- Velvety soft surface structure
- Highly elastic bending zone, also protects the popliteal fossa when worn for longer periods

BORT Stabilo® Knee Support with Articulated Joint REF 150 140 KIDS



PRODUCT-BENEFIT

- Knee support with side joint splints for guidance of the knee joint
- Viscoelastic pad
- Circular fastening straps on the thigh and lower leg
- Anatomically contoured knit
- Compression class II

COLOUR: silver

THE KIDS-BENEFIT

- Extension limit at 0°
- Individually mouldable joint splints

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. ligament instability, meniscopathy

INDICATIONS

Mild to moderate lateral ligament instability of the knee joint

SIZE d-measure f-measure

Children 23 – 26 cm 28 – 32 cm

Adult sizes see page 138

right and left symmetrical 1 cm = 0.3937 inch

d-measure	f-measure

Children 23 – 26 cm 28 – 32 cm

Adult sizes see page 147

right and left symmetrical 1 cm = 0.3937 inch

SIZE

BORT Immobilisation Splint with Patella Recess REF 145 000 KIDS



PRODUCT-BENEFIT

- Knee brace for immobilisation of the knee joint with a 0° and/or 20° flexion
- 3-piece with open patella
- Median and lateral aluminium splints, which can be freely positioned
- Length: 38 cm (14.9 inches)

COLOUR: 🗖 black

THE KIDS-BENEFIT

- The support adapts to body shape thanks to a 3-piece design
- Lined with terry cloth on the inside for excellent comfort
- The aluminium splints can be removed for washing

INDICATIONS

Pre-operative, post-operative, posttraumatic, joint, muscle or ligament injuries

BORT Helix S Spiraldynamik® Lower Ankle Support [REF] 053 200 KIDS



RODUCT-BENEFIT

- Ankle bandage to support torsion of the back of the foot and to raise the instep
- Thin, elastic straps
- 3-D contoured knit
- Two silicone pads
- Extra soft, snag proof material edge at the calf

COLOUR: 🗖 black

THE KIDS-BENEFIT

- Stabilisation of the ankle and activation of foot muscles through compression
- Functional supportive restraint in the case of classical foot deformities (i.e. splayfoot)
- Optimal comfort even during intense activity
- Separate slip-on aid included

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of lower ankle joint and the sole of the foot, ligament instability, foot problems due to static anomalies and defects in the shape of the foot due to dysfunction of the lower ankle joint e. g. skew and splayfoot

SPLINT LENGTH


Fits height

Children ______ 80 – 145 cm Please indicate 0° or 20° flexion Adult sizes see page 154

right and left symmetrical 1 cm = 0.3937 inch

SIZE

Circumference of ankle



please indicate left or right 1 cm = 0.3937 inch

34 BORT. Benefit at your side.



PRODUCT-BENEFIT

- Active ankle support for compression of soft tissues
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II

COLOUR: blue

BORT MalleoStabil® Brace REF 100 500 KIDS

PRODUCT-BENEFIT

 Combined Air/Gel brace and foam pad to stabilise the ankle, while limiting pronation and supination

COLOUR: grey



- With memory foam padding only



- With Air/Gel padding only

THE KIDS-BENEFIT

- No impairment when walking
- Easy to put on thanks to circular Velcro fastening

THE KIDS-BENEFIT

- Velvety soft surface structure
- Pressure relief at the base of the ankle

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of lower ankle joint and the sole of the foot, sprain, ligament weakness, after removal of plaster cast, ligament injuries

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament ruptures, ligament injuries, supination/pronation trauma, fibular ligament rupture

SIZE	Circumference of ankle

b-measure	Children	 17	-	19	ст

Adult sizes see page 162

right and left symmetrical 1 cm = 0.3937 inch

SIZE	Fits height
Size 1/Children	up to 155 cm
Adult sizes see page 169	

right and left symmetrical 1 cm = 0.3937 inch

BORT Heel Lift

made of Silicone

REF 930 200 KIDS

- Silicone heel lift to compensate dif-

- Heights: 3 mm, 5 mm, 8 mm, 10 mm

(0.1, 0.2, 0.3 and 0.4 inches)

ferences in leg-length

- Pack: Single item

- Made of durable silicone



PRODUCT-BENEFIT

- Soft ankle brace for stabilisation of the lower and upper ankle joints
- Bilateral shells padded with soft memory foam
- Cross restraint featuring Dynamic-Force-Control (DFC) technology
- Elastic material that carries away moisture
- Brace can be opened fully from the rear for easy access

COLOUR: 🗖 black

THE KIDS-BENEFIT

- Slim design

- Can be worn without difficulty with sports and walking shoes
- Limitation of talus prolapse with optimal stabilisation
- For positioning the ankle in pronation or supination

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries, supination / pronation trauma

THE KIDS-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Compensation of differences in leg-length, achilles tendon disorders (Achillodynia)

SIZE	Circumfere	nce of ankle		
b-measure	Children	17 – 19 cm		
Adult sizes see page 167				

please indicate left or right 1 cm = 0.3937 inch

SIZE	Shoe size UK
Children	1-2
Adult sizes see page 179	
right and left symmetrical	

1 cm = 0.3937 inch

Keep your posture: Stabilisation and relief of cervical spine and shoulder.



 O V E	R١	V I [W	C) F	IN	DI	C A	AT I	0	N S	37	!

Overview of indications for selected products	Cervical Support 127 260/ 360/660	Clavicle Support 121 700	OmoFX 121 100	OmoFX Vario Generation 121 020 Shoulder Brace 12 215 700 11	Generation Shoulder Brace 215 700	Omo- Basic [®] 121 350/ 121 500	Omo- Xpress 120 900 Omo-	0mo- Stabil® 121 200	OmoSAT 121 980	OmoARS 121 900	OmoTex Traction 120 800	Omo Control 120 700
	Stabilo [®] Nec 127 560						Xpress light 120 920					
Page	P. 38/39	P. 40	P. 46	P. 43	P. 41	P. 45	P. 44	P. 43	P. 42	P. 46	P. 49	P. 48
Position Outward rotation			.06	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~			-			1.35°an		
Abduction				- AL			1			700 / 200		
Adduction			1.0	<u>Í</u> ,			20°-0°				<u>></u>	-
Traction			Abduction	Abduction		4	Adduction			Outward rotation / Abduction	Traction	Reduction
Cervical syndrome	•											
Clavicle fracture		•										
Arthrosis					•							
Impingement syndrome											•	
Supraspinatus syndrome			•									
Contusion						•	•					
Post-operative			•	•	•	•	•	•	•	•		
Dislocation						•	•		•			•
Subcapital humerus fractures					•							
Signs of paralysis					•				•			•
Injury to the rotator cuff								•				
Slight to moderate instability												•



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of the cervical spine
- Stiff, open-pored foam, for daytime wear
- With recess for larynx
- Adjustable Velcro fastening

COLOUR: beige

THE MED-BENEFIT

- Particularly rounded edges
- Skin-friendly 100% cotton cover

BORT Cervical Support REF 127 260/360/660



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- Soft, open-pored foam
- With recess for larynx
- Adjustable Velcro fastening
- For children too, see page 28

THE MED-BENEFIT

- Can be worn day and night
- Particularly rounded edges
- Skin-friendly 100% cotton cover

INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, cervical syndrome, post-operative

SIZES	Circumference of neck	SIZES	Circumference of neck
	small up to 38 cm medium 38 - 43 cm large 43 - 48 cm	small	127 260, blue) _ up to 33 cm up to 38 cm 38 - 43 cm 43 - 48 cm

Chin heights 5,5 cm / 7,5 cm / 9,0 cm / 11,0 cm Children 5,0 cm / 6,0 cm



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- Integrated reinforcement splint
- Stiff, open-pored foam
- With recess for larynx
- Adjustable Velcro fastening

COLOUR: beige

BORT Cervical Support Eco REF 127 460



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- With recess for larynx
- Soft, open-pored foam
- Adjustable Velcro fastening
- Length: 56 cm (22.05 inches)

COLOURS: beige blue/black

THE MED-BENEFIT

- Functional stabilisation
- Particularly rounded edges
- Skin-friendly 100% cotton cover

THE ECO-BENEFIT

- It can be worn day and night
- Particularly rounded edges
- Skin-friendly 100% cotton cover

INDICATIONS

Severe painful conditions, minor whiplash injuries, cervical syndrome, torticollis (rheumatic), osteochondrosis/ spondylarthrosis of the cervical spine, post-operative

INDICATIONS

Pain in the cervical spine area, minor whiplash injuries, sprain, post-operative

SIZES

Circumference of neck



small	up to 38 cm
medium	38 – 43 cm
large	43 – 48 cm

Chin heights 5,5 cm / 7,5 cm / 9,0 cm / 11,0 cm

SIZES		Ch	in height
*	22.1 inch circumfe	th only (56 ies), to fit n rence from 2.6 to 17.7	eck 32 to



PRODUCT-BENEFIT

- Anatomical cervical support for stabilisation of cervical spine
- Integrated reinforcement splint
- With recess for larynx
- Stiff open-pored foam
- Adjustable Velcro fastening
- Length: 56 cm (22.05 inches)

COLOURS: beige blue/black

THE ECO-BENEFIT

- Functional stabilisation
- Particularly rounded edges
- Skin-friendly 100% cotton cover



PRODUCT-BENEFIT

- Backpack support
- Soft padded straps
- Three Velcro fastenings, individually adjustable, ensure an ideal blockage
- For children too, see page 29

COLOURS:
white (adults)
blue (children)

THE MED-BENEFIT

- Conforms ideally to the body's contours, especially in the critical armpit area
- Skin-friendly
- Can be worn day and night

INDICATIONS

Severe painful conditions, minor whiplash injuries, cervical syndrome, post-operative

SIZES		Chin height
	Size 0 Size 1 Size 2 Size 3	5,5 cm 7,5 cm 9,0 cm 11.0 cm

One length only (56 cm/22.1 inches), to fit neck circumference from 32 to 45 cm (12.6 to 17.7 inches).

INDICATIONS

Treatment of clavicle fractures

SIZES	Under	bust ci	ircumference
	Children (only blue) small medium large x-large		<mark>up to 65 cm</mark> 65 – 80 cm 80 – 95 cm 95 – 110 cm 10 – 125 cm

1 cm = 0.3937 inch

1 cm = 0.3937 inch

BORT Generation Shoulder Brace REF 215 700



PRODUCT-BENEFIT

- Shoulder/arm adduction brace with pouch for immobilisation of the shoulder joint
- Waist and shoulder strap for safe fixing of the shoulder and arm



COLOUR: blue

ACCESSORIES: REF 991 115 BORT Hand rest with a ball



BORT Generation: For better mobility for the aged Product range with special handling advantages

THE GENERATION-BENEFIT

- Simple fitting process for the shoulder and waist strap
- Positioning aid for the back prevents twisting of the straps
- Quick-release fastener for simple and safe attachment of the arm pocket to the strap system



INDICATIONS

Pre-operative, post-operative, post-traumatic, e.g. sprain, contusion, **subcapital** humerus fractures, signs of paralysis, inflammatory, **degenerative**

SIZES	Fits height
Size 1	up to 160 cm
Size 2	160 – 180 cm
Size 3	over 180 cm



PRODUCT-BENEFIT

- Shoulder/arm adduction brace with pouch for immobilisation of the shoulder joint
- Hand rest with ball for training forearm muscles

THE MED-BENEFIT

- Simple fitting process for the shoulder and waist strap
- Quick-release fastener for simple and safe attachment of the arm pocket to the strap system



COLOUR: grey

INDICATIONS

Pre-operative, post-traumatic, post-operative, sprain, contusion, subcapital humerus fractures, signs of paralysis

SIZES	Fits height
Size 1	up to 160 cm 160 – 180 cm over 180 cm

BORT OmoStabil® Desault's Bandage [REF] 121 200

PRODUCT-BENEFIT

- Shoulder/arm adduction brace as a so-called **Desault's Bandage** for the immobilisation of the shoulder joint in a defined position
- Cap-shaped recess in the shoulder area
- Fastens around the torso area
- Perforated composite velour material

COLOUR: 🗖 blue



THE MED-BENEFIT

 Surface fixation of the shoulder and arm allowing free movement of the hand

BORT OmoFX Vario REF 121 020



PRODUCT-BENEFIT

- Shoulder/arm abduction brace for the immobilisation and relief of the shoulder joint
- Freely positionable straps for fixing the arm
- Ball-shaped hand rest for training the forearm muscles
- Skin-friendly Velour cover



THE MED-BENEFIT

COLOUR: grey

- Reversible cushion: positioning in 30° or 45° abduction
- High therapeutic safety due to the anti-slip anatomic cut of the cushion
- Comfortable to wear

INDICATIONS

Pre-operative, post-traumatic, postoperative, contusion, dislocation, subcapital humerus fractures, injury to the rotator cuff

INDICATIONS

Pre-operative, post-operative, posttraumatic, for example in the case of impingement syndrome, rotator cuff injury, implantation of shoulder prosthesis

SIZES

Circumference of waist

 Size 0
 up to 75 cm


 Size 1
 75 – 95 cm

 Size 2
 95 – 125 cm

 Size 3
 125 – 145 cm

right and left symmetrical 1 cm = 0.3937 inch SIZES

Circumference of torso





PRODUCT-BENEFIT

- Shoulder/arm adduction brace for immobilisation of shoulder joint
- Stable fixation of the shoulder and arm
- Open design for easy wound inspection
- Skin-friendly material

COLOUR: 🔳 grey



THE MED-BENEFIT

- Hand strap can be fully opened
- Easy to handle by rotating buckle, which allows an easy early-functional post-treatment
- Comfortable to wear

BORT OmoXpress light REF 120 920



PRODUCT-BENEFIT

- Shoulder/arm adduction brace for immobilisation of shoulder joint
- Stable fixation of the shoulder and arm
- Open design for easy wound inspection
- Skin-friendly material

COLOUR: grey



THE MED-BENEFIT

- Very easy adaptation and safe handling due to the Velcro fastening
- Freely positionable arm slings
- Slight, visually hardly noticeable care, is comfortable to wear

INDICATIONS

Pre-operative, post-traumatic, post-operative, sprain, contusion, dislocation, humerus fractures

INDICATIONS

Pre-operative, post-traumatic, post-operative, sprain, contusion, dislocation, humerus fractures

SIZES

Underbust circumference

)={	Size O	_ up to 80 cm
$\langle \rangle$	Size 1	_ 80 – 120 cm
$\left \right\rangle^{*}$	Size 2	_120 – 160 cm
$I \Lambda = \Lambda \Lambda$		

right and left symmetrical 1 cm = 0.3937 inch SIZES

Underbust circumference

	Size 0	up to 80 cm
$\langle \rangle$	Size 1	80 – 120 cm
$\left \right\rangle^{*} + \left \left\langle \right\rangle$	Size 2	120 – 160 cm



PRODUCT-BENEFIT

- Open shoulder/arm adduction brace for immobilisation of the shoulder joint
- Stabilising strapping system ensures secure fit and immobilisati on of the arm in the adduction position
- Soft, tension-resistant straps
- For children too, see page 28

COLOURS: 🗖 blue 🗌 white



THE BASIC-BENEFIT

- Flexible hand and arm adjustment
- Easy access to the surgical area to allow for wound care
- Skin-friendly

BORT OmoBasic Enclosed Style REF 121 500



PRODUCT-BENEFIT

- Closed shoulder/arm adduction brace for immobilisation of shoulder join
- Highly elastic hose material
- Elastic support band at the base of the thumb
- Stabilising strapping system ensures secure fit and immobilisation of the arm in the adduction position

COLOURS: blue white



THE BASIC-BENEFIT

- Stabilisation of the shoulder while allowing free movement of the thumb and fingers
- Elastic border in the shoulder area
- Skin-friendly

INDICATIONS

Pre-operative, post-operative, posttraumatic, sprain, dislocation, subcapital humerus fractures, scapular fracture, rotator cuff injury

SIZES	Length Olecranor metacarpal knuckle poin	
	Children (blue) 22 – 25 Size 0/S 25 – 28 Size 1/M 28 – 31 Size 2/L 31 – 34 Size 3/XI 34 – 38	cm cm cm

right and left symmetrical 1 cm = 0.3937 inch

INDICATIONS

Pre-operative, post-traumatic, postoperative, sprain, dislocation, subcapital humerus fractures, scapular fracture, rotator cuff injury

SIZES		Arm circumference	Arm length
Am long	0/S	26 - 32 cm	48 – 52 cm
	1/M	28 - 34 cm	50 – 56 cm
	2/L	32 - 40 cm	54 – 60 cm
	3/XL	38 - 48 cm	58 – 68 cm



PRODUCT-BENEFIT

- Shoulder/arm abduction cushion for immobilisation and relief of the shoulder joint
- Abduction: ca. 30°
- Stabilisation straps for the shoulder and elbow
- Skin-friendly Velour cover



COLOUR: grey

THE MED-BENEFIT

- Easy to handle
- Comfortable arm rest

INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of periarticular irritations and disorders in the glenohumeral joint area, acromioclavicular joint dislocation (Tossy Classification Grades I + II), supraspinatus syndrome

SIZES	Fits height
-------	-------------

Size 1 _	u	p to	165	cm
Size 2 _	(ver	165	ст


please indicate left or right 1 cm = 0.3937 inch

BORT OmoARS REF 121 900



PRODUCT-BENEFIT

- Arm support cushion for immobilisation and relief of the shoulder joint with outward rotation at 15° or 30° and 15° abduction
- Strap system with Velcro fastening
- Hand rest with ball for training forearm muscles



COLOUR: grey

THE MED-BENEFIT

- Safe stability and positioning through anatomical cushion
- Easy to handle
- Skin-friendly Velour cover

ACCESSORIES: REF 121 960 Wedge for positioning of outward rotation at 30°

INDICATIONS

Pre-operative, post-traumatic, post-operative, for example in the case of impingement syndrome, rotator cuff rupture, implantation of shoulder prosthesis, subacromial defilee dilation (acc. to Neer), subacromial bursitis

SIZES	Fits height
Size 1 Size 2	up to 165 cm over 165 cm
JIZE Z	

_SHOULDER 47

BORT OmoARS with 3-Point-Strap [REF] 121 970



PRODUCT-BENEFIT

- Arm support cushion for immobilisation and relief of shoulder joint in outward rotation at 15° or 30° and 15° abduction
- 3-point strap system with Velcro fastening
- Hand rest with ball for training the forearm muscles



COLOUR: grey

THE MED-BENEFIT

- 3-point strap system ensures optimal fit and stability
- Easy to handle
- Skin-friendly Velour cover

ACCESSORIES: REF 121 960 Wedge for positioning of outward rotation at 30°

INDICATIONS

Pre-operative, post-traumatic, post-operative, for example in the case of impingement syndrome, rotator cuff rupture, implantation of shoulder prosthesis, subacromial defilee dilation (acc. to Neer), subacromial bursitis

SIZES	Fits height
Size 1	up to 165 cm
Size 2	over 165 cm

right and left symmetrical 1 cm = 0.3937 inch

BORT OmoARS short REF 121 030



PRODUCT-BENEFIT

- Arm support cushion for immobilisation and relief of the shoulder joint in outward rotation at 15° and in abduction at 15°
- Strap system with Velcro fastening
- Hand rest with ball for training the forearm muscles



COLOUR: grey

THE MED-BENEFIT

- Safe stability and positioning through anatomical cushion
- Skin-friendly Velour cover

SIZES	Length metacarpal kn	Olecranon to uckle point V
	Size 1 Size 2 Size 3	up to 34 cm 35 – 39 cm over 39 cm



PRODUCT-BENEFIT

- Shoulder joint brace for physiological guidance and stabilisation of the shoulder joint
- Correction of internal rotation and return of the humeral head
- Relief of the capsules, tendons, muscles and nerves
- Enables free pendular motion of the arm

COLOUR: grey



THE MED-BENEFIT

- Comfortable materials for maximum wearing comfort
- Its construction allows single-handed application
- Easy handling due to innovative BOA® Fit system

Arm and shoulder pads can also be ordered separately in (S), (M) and (L)

REF 120 500 Arm pad **REF** 120 600 Shoulder pad

INDICATIONS

Slight to moderate instability of the shoulder joint, **conservative in case of dislocation of the shoulder**, inflammatory conditions, degenerative

SIZES		Arm cuff (Module A) lower arm circumference (5 cm under the elbow joint)	Shoulder cuff (Module B) chest circumference (under the armpit, measured in cm)
J.	Size 2	(S) 20 – 25 cm	(S) 75 – 90 cm (M) 90 – 105 cm (L) 105 – 125 cm
J.	Size 5	(M) 25 – 30 cm	(S) 75 – 90 cm (M) 90 – 105 cm (L) 105 – 125 cm
	Size 8	(L) 30 – 35 cm	(S) 75 – 90 cm (M) 90 – 105 cm (L) 105 – 125 cm

please indicate left or right 1 cm = 0.3937 inch

BORT OmoTex REF 121 600



PRODUCT-BENEFIT

- Shoulder Joint compression support
- Additional Velcro straps to restrict outward rotation
- Silicone pad, which can be freely positioned
- Anatomical, muscle-activating contoured knit

COLOUR: 🗖 black

THE MED-BENEFIT

- Improved guidance of the glenohumeral joint
- Length of the straps can be adjusted individually
- Comfortable to wear even under clothing

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the shoulder area, arthrosis

SIZES Circumference of upper arm

jël	Size 1	 24 – 28 cm
6)	Size 2 🔔	 28 – 32 cm
apper ann	Size 3 🔔	 32 – 36 cm
6	Size 4 🔔	 36 – 40 cm

right and left symmetrical 1 cm = 0.3937 inch

BORT OmoTex Traction REF 120 800



PRODUCT-BENEFIT

- Modular shoulder joint support with separately useable weight band for subacromial decompression
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Four 250 g weights for individual traction weighting
- Improved adaptability due to the flexibly adjustable belt guidance plate

COLOUR: 🗖 black



THE MED-BENEFIT

- Improved joint guidance and relief of pressure in the humeral head region
- Weight band removable from the support and adjustable with Velcro tape in accordance with the course of therapy
- Comfortable to wear even under clothing

INDICATIONS

Impingement syndrome, chronic, post-traumatic or post-operative irritations of soft tissues in the shoulder area, periarthropathia humero-scapularis (PHS)

SIZES Circumference of upper arm

)≋{ Size 1 24 -	- 28 cm
Size 2 28 -	- 32 cm
(Size 3 32 -	- 36 cm
Size 4 36 -	- 40 cm

Remain mobile: Relief and Stabilisation for the elbow



_____OVERVIEW OF INDICATIONS 51

Overview of indications for selected products	select [*] EpiPlus [®] 122 700	EpiBasic Sport 122 600 SP	EpiBasic 122 600 EpiBasic Short 122 600 K	Kubi Tal® 055 100	EpiContur® Sport 022 500 SP	Stabilo [®] Epi Brace 122 580	Elbow Support 104 200	KubiFX 055 300 KubiFX light 055 350	KubiFX long 055 400 KubiFX long light 055 450	Elbow ROM Brace 055 500
Page	P. 52	P. 52	P. 53	P. 54	P. 57	P. 59	P. 55	P. 60	P. 61	P. 59
Post-traumatic			•					•	•	•
Post-operative							•	•	•	•
Dislocation										•
Sulcus-ulnaris syndrome				•				•		
Tendopathies						•			•	
Epicondylitis	•	•	•		•	•			•	
Irritations	•	•		•	•					•
Tendomyopathy			•							
Sprain							•			
Arthrosis							•			

Elbow



PRODUCT-BENEFIT

- High quality elbow brace with additional Stabilo[®] Strap
- Two structured silicone pads for intermittent massage
- Positioning aid in the antecubital fossa
- Lengthwise-elastic stretch zone in the elbow area ensures an ideal fit at all bending angles
- Compression class II
- Length: 20 cm (7.9 inches)

COLOURS: silver black

THE SELECT-BENEFIT

- SOFTflex Plus knit fabric:

- SENSIL® Breeze: Luxury fibre with a cooling effect
- · Soft-velvety surface structure
- Extra-soft, snag proof material edge
- Very comfortable to wear, ergonomic and slip-free

BORT EpiBasic Sport REF 122 600 SP



PRODUCT-BENEFIT

- Elbow compression support with additional Stabilo[®] Strap
- SOFTflex 3D knitted fabric with a velvety soft surface structure
- Two structured silicone pads
- Positioning aid in the antecubital fossa
- Flexible stretching zone in the elbow area
- Length: 20 cm (7.9 inches)

COLOUR: I black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body
- Easy fixing and length setting of the Stabilo[®] Strap using the easyClick system
- Latex-free knitted fabric

Circumference of forearm 5 cm (2 inches) below elbow joint

 The pressure on the pads can be individually adjusted through Stabilo[®] Strap

INDICATIONS

SI7ES

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **epicondylitis** (acute and chronic), tendomyopathy

01220	
x-small small medium large x-large	20 – 22 cm 22 – 24 cm 24 – 26 cm 26 – 28 cm 28 – 30 cm 30 – 32 cm 32 – 34 cm

BORT Epi Basic REF 122 600



PRODUCT-BENEFIT

- Elbow compression support
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Two structured silicone pads for intermittent massage
- Positioning aid in the antecubital fossa

BORT EpiBasic Short REF 122 600 K



PRODUCT-BENEFIT

- Elbow compression support
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Two structured silicone pads for intermittent massage
- Positioning aid in the antecubital fossa

COLOURS: beige silver black

THE BASIC-BENEFIT

- Velvety soft surface structure
- Special elastic bending zone, also protects the antecubital fossa when worn for longer periods

COLOUR: silver

THE BASIC-BENEFIT

- Velvetv soft surface structure
- Special elastic bending zone, also protects the antecubital fossa when worn for longer periods
- Reduced upper arm length

*Length information only denotes approximate measures and may vary by size.

INDICATIONS

01750

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, epicondylitis (acute and chronic), tendomyopathy 0. (((-

21ZE2	CIrcumference of forearm 5 cm (2 inches) below elbow joint
x-small small medium large x-large xx-large	20 – 22 cm 22 – 24 cm 24 – 26 cm 26 – 28 cm 28 – 30 cm 30 – 32 cm 32 – 34 cm
"xxx-larde (d	only silver) 34 – 36 cm

*xxx-large only for REF 122 600



PRODUCT-BENEFIT

- Extensive padded elbow support
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Extensive silicone padding with recess at olecranon for dimensional stability during flexion
- Compression class II
- COLOURS: beige grey

THE MED-BENEFIT

- Velvety soft surface structure
- Extensive protection of entire elbow area
- Seamless comfort zone in the antecubital fossa

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, for example in the case of compression of the nervus ulnaris at elbow (sulcus ulnaris syndrome/cubital tunnel syndrome), bursitis

SIZES	Circumferenc 5 cm (2 inches) belo	
x-small small medium large x-large		$\begin{array}{c} 20-22 \text{ cm} \\ 22-24 \text{ cm} \\ 24-26 \text{ cm} \\ 26-28 \text{ cm} \\ 28-30 \text{ cm} \\ 30-32 \text{ cm} \\ 32-34 \text{ cm} \end{array}$



PRODUCT-BENEFIT

- Elbow compression support
- Ulnar and radial silicone pads
- Compression class II

COLOURS: beige mineral grey graphite

THE MED-BENEFIT

- Elastic comfort zone in the antecubital fossa for pleasant wearing comfort
- Soft surface structure



PRODUCT-BENEFIT

- Elbow compression support
- Anatomically contoured knit for optimal fit
- Compression class II

COLOUR: beige

THE MED-BENEFIT

- Comfortable to wear
- Not bulky

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, **epicondylitis** (acute and chronic), tendomyopathy


Circumference of forearm SIZES 5 cm (2 inches) below elbow joint

x-small _____ 22 – 24 cm small _____ 24 – 26 cm medium _____ 26 – 28 cm large _____ 28 – 30 cm

x-large _____ 30 – 32 cm

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, articular effusion, sprain, arthrosis, contusion



right and left symmetrical 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Elbow compression support with lateral splints
- Rounded knitted fabric, seamless
- Elastic in longitudinal and lateral directions
- Compression class II

COLOUR: silver

THE MED-BENEFIT

- Comfortable to wear
- Lateral splints for extra support

BORT Dual-Tension Elbow Support **REF** 054 200



PRODUCT-BENEFIT

- Circular knitted elbow compression support
- Elastic in longitudinal and lateral directions
- Compression class II

COLOUR: beige

- THE MED-BENEFIT
- Comfortable to wear thanks to seamless knitted fabric

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of the elbow joint, articular effusion, sprain, arthrosis, contusion

SIZES	
-------	--

Circumference of elbow

	small medium large x-large	 26 – 28 –	26 cm 28 cm 30 cm 32 cm
للمسا	x largo	 00	02 0111

right and left symmetrical 1 cm = 0.3937 inch

SIZES	Circu	umference of elbo	W
	x-small small medium large x-large	22 - 24 c 24 - 26 c 26 - 28 c 28 - 30 c 30 - 32 c	m m m

Elbow

BORT EpiContur[®] Sport REF 022 500 SP

PRODUCT-BENEFIT

- Epicondylitis soft brace for relieving the radial or ulnar muscle origin
- Standard size adjustable to individual circumferences using Velcro fastening
- COLOUR: black/green

THE SPORT-BENEFIT

- Fitted cut brace design enables slipfree fit for sport activities
- Brace pressure individually metered
- Specially shaped friction pad, freely positionable
- Modern, sporty design (with a reflecting label)

INDICATIONS

Epicondylitis, chronic or recurrent epicondylopathy, tendopathy, overstraining caused by work and sport, 'tennis/golf elbow'

> Circumference of forearm 5 cm (2 inches) below elbow joint

____ 21 – 35 cm

SIZE

One size only _



PRODUCT-BENEFIT

- Epicondylitis soft brace for relieving the radial or ulnar muscle origin
- Standard size adjustable to individual circumferences using Velcro fastening
- Freely positionable friction pad

COLOURS: beige grey

THE MED-BENEFIT

- Brace pressure individually metered
- Fitted cut brace design enables slipfree fit

BORT EpiContur® with 2 Pads REF 122 500



PRODUCT-BENEFIT

- Soft epicondylitis brace to relieve the muscle origins at the ulnar and/or radial epicondyle
- Ulnar and radial friction pad, which can be freely positioned

COLOURS: beige grey

THE MED-BENEFIT

- Brace pressure is individually adjustable
- Fitted cut brace design enables slipfree fit

INDICATIONS

Epicondylitis, chronic or recurrent epicondylopathy, tendopathies, overstraining, 'tennis/golf elbow'

SIZE	Circumference of forearm 5 cm (2 inches) below elbow joint	SIZES	Circumference of forearm 5 cm (2 inches) below elbow joint
	One size only 21 – 35 cm		Size 0 21 – 25 cm Size 1 25 – 30 cm ∋ Size 2 30 – 35 cm

right and left symmetrical 1 cm = 0.3937 inch

BORT Stabilo® Epicondylitis Strap with Ulnar Relief [REF] 122 580

PRODUCT-BENEFIT

- Epicondylitis brace to relieve the muscle origins at the ulnar and / or radial epicondyle
- Special anatomical design allows to individually meter strap pressure

COLOUR: silver

THE MED-BENEFIT

- Elastic Velcro fastening, which results in no pressure on the elbow
- Simple design
- Padded inside

INDICATIONS

Epicondylitis, chronic or recurrent epicondylopathy, tendopathies, overstraining, tennis elbow, golf elbow



PRODUCT-BENEFIT

- Elbow brace for mobilisation with simultaneous immobilisation of the proximal radio ulnar joint
- Limitation of extension: from 0° 90°
- Limitation of flexion from 0° 110°
- Fixation of the joint at 0°, 10°, 20°, 30°, 40°, 50°, 60°, 70°, 80°, 90°
- Avoidance of supination and pronation in the lower arm with a hand rest

COLOUR: 🔳 gre🖗

THE MED-BENEFIT

- Comfortable to wear due to extensive contact surfaces on the upper arm and lower arm
- Individual length setting for the upper arm and lower arm
- Extra lightweight brace construction

INDICATIONS


All indications that require early functional mobilization with movement limitation of the elbow joint and immobilization of the proximal radius-ulnar joint is required, such as. b., conservative/postoperative/degenerative care, mobilisation after surgery and trauma involving the proximal radioulnar joint, inflammations, ize radial head fractures, elbow luxation, distal humerus fractures, i.e after prosthesis implant

	Circumference of forearm
SIZES	5 cm (2 inches) below elbow joint

	,	
Size 1 Size 2 Size 3 Size 4 Size 5	 $\begin{array}{r} 22-24\\ 24-26\\ 26-28\\ 28-31\\ 31-34 \end{array}$	cm cm cm

right and left symmetrical 1 cm = 0.3937 inch

SIZES Circumference of upper arm



please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Brace for immobilisation of the elbow joint
- Two removable aluminium splints
- Non-elastic material
- Inner side lined with terry cloth
- Length: ca. 39 cm (15.4 inches)

BORT KubiFX Elbow Brace REF 055 300



PRODUCT-BENEFIT

- Brace for immobilisation of elbow joint
- Two removable aluminium splints
- Non-elastic material
- Inner side lined with terry cloth
- Freely adjustable elbow pad
- Four reversible straps
- Length: ca. 37 cm (14.6 inches)

COLOUR: 🔳 grey

THE MED-BENEFIT

- Secure support
- Individually mouldable from 0° to 90°
- Extra light and open frame construction

COLOUR: grey

THE MED-BENEFIT

- Reliable support
- Individually shaped from 0° to 90°

INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of proximal sulcus ulnaris syndrome, insertional tendinopathy, post-bursectomy conditions

SIZE

One size only

SIZE

One size only

right and left symmetrical 1 cm = 0.3937 inch

BORT KubiFX Long Elbow Brace REF 055 400



PRODUCT-BENEFIT

- Brace for immobilisation of the elbow joint and the proximal radioulnar joint
- Volar splint extending to palm area
- Two removable aluminium splints
- Non-elastic material
- Inner side lined with terry cloth
- Freely positionable elbow pad
- Five reversible straps
- Length: ca. 53 cm (20.9 inches)

COLOUR: grey

- THE MED-BENEFIT
- Reliable support
- Individually shaped from 0° to 90°

BORT KubiFX Long Light Elbow Brace REF 055 450



PRODUCT-BENEFIT

- Brace for immobilisation of the elbow joint and the proximal radioulnar joint
- Inner side lined with terry cloth
- Hand rest which can be used on both sides
- Freely adjustable elbow pad
- Standard size due to stepless setting of the lower arm length (olecranon to metacarpal knuckle point V) from 28 cm (11.0 inches) to 38 cm (15.0 inches)

COLOUR: grey

THE MED-BENEFIT

- Reliable support
- Individually shaped from 0° to 90°
- Extra light and open frame construction

INDICATIONS

Pre-operative, post-operative, post-traumatic, for example in the case of proximal sulcus ulnaris syndrome, insertional tendinopathy, epicondylitis, post-bursectomy conditions

	Length Olecranon to
SIZES	metacarpal knuckle point V

SIZE

Size 1 _____ 28 – 31 cm Size 2 _____ 31 – 34 cm Size 3 _____ 34 – 38 cm

please indicate left or right 1 cm = 0.3937 inch One size only

Fit for everyday commitments: Pain reduction, more stability and functionality for the fingers and thumb.




Overview of indications for the selected products	MobiDig [®] Flexion 112 560	MobiDig [®] Extension 112 510	DigiSoft [®] Vario Finger Brace 112 790	SOFT Thumb Splint Long 112 720	SOFT Thumb Splint Plus 112 780	Generation Thumb Splint 215 300	Rhizo Flex® 105 100	Sella- Xpress Short 112 840	Sella- Xpress Plus 112 860	Thumb Splint Sport 112 180 SP	SellaTex [®] SellaTex [®] Plus – Classic Quervain 112 050 C 112 160 SellaTex [®] 112 050 C 112 050	SellaTex [®] Classic 112 050 C SellaTex [®] 112 050	Generation Thumb-Wrist Brace 215 200
Page	P. 65	P. 65	P. 67	P. 70	P. 71	P. 73	P. 73	P. 71	P. 72	P. 68	P. 77	P. 78/79	P. 76
Immobili- Metacarpophalangeal joint	t			•	•	•		•	•		•	•	•
sation Saddle joint				•	•	•	•	•	•		•	•	•
Interphalangeal joint					•				•		•		
Flexion DIP/PIP joints	•												
Extension PIP joints		•											
Positioning			•										
Ulnar/radial lesion of ligaments						•					•	•	•
Rhizarthrosis				•	•	•	•	•	•	•		•	•
Tendovaginitis de Quervain											•	•	
Tendopathy					•				•			•	•
Arthrosis						•	•			•	•	•	•
Additional wrist immobilisation											•	•	•

Fingers Thumb



PRODUCT-BENEFIT

- Plastic finger splint for hyperextension of the last damaged phalanx



- THE MED-BENEFIT
- Conservative treatment
- Beige colour: Can be thermoplastically remoulded
- Machining possible



PRODUCT-BENEFIT

- Plastic finger splint for treating buttonhole deformity (Boutonnière Deformity) on the long fingers of the hand
- Immobilisation of the PIP joint

COLOUR: beige

THE MED-BENEFIT

- Conservative treatment
- Thermoplastically mouldable
- Machining possible

Fixing tape, e.g. REF 470 400 is required for functionality.

INDICATIONS

Fracture, dislocation, ligament injuries, extensor tendon rupture

INDICATIONS

Buttonhole deformity (Boutonnière Deformity), immobilisation of the PIP joints

SIZES		F	inger width at DIP joint
- A A	Size 1 Size 2 Size 3 Size 4 Size 5 Size 5.5 Size 6 Size 7		1,4-1,6 cm 1,5-1,8 cm 1,7-2,1 cm 1,8-2,2 cm 1,9-2,3 cm 2,3-2,6 cm 2,0-2,8 cm 2,6-2,9 cm

Size 1 1,9 cm Size 2 2,1 cm Size 3 2,3 cm Size 4 2,4 cm Size 5 2,5 cm Size 6 2,6 cm Size 7 2,8 cm	SIZES		diameter
		Size 2 Size 3 Size 4 Size 5 Size 6	2,1 cm 2,3 cm 2,3 cm 2,4 cm 2,5 cm 2,6 cm

Innor colint



PRODUCT-BENEFIT

 Dynamic Quengel splint for mobilisation of interphalangeal joints

- Tension can be increased or reduced

by carefully bending the spring

- Ergonomic, durable padding

COLOUR: White

THE MED-BENEFIT

BORT MobiDig® Extension Finger Splint REF 112 510

PRODUCT-BENEFIT

- Dynamic Quengel splint for mobilisation of interphalangeal joints
- Proximal padding in the palm area
- Ergonomic, durable padding

COLOUR: grey

THE MED-BENEFIT

- Linear spring force
- Tension can be increased or reduced by carefully bending the spring
- Plus sizes

INDICATIONS

Post-operative, post-traumatic flexion of PIP (L₁) or DIP (L₂) joints

INDICATIONS

Post-operative, post-traumatic extension of PIP joints

SIZES Length	Finger SIZES Length ¹⁾ width ²⁾
small 3,5 – 5,0 cm medium 5,0 – 5,8 cm large 5,8 – 6,5 cm x-large over 6,5 cm	x-small plus 5,0- 6,3 1,5-1,7 small 6,3- 7,6 1,5-1,7
L ₁ = Flexion PIP joint or L ₂ = Flexion DIP joint	¹⁾ Distal fold of palm to distal finger fold ²⁾ at PIP joint
right and left symmetrical	right and left symmetrical

1 cm = 0.3937 inch

right and left symmetrica 1 cm = 0.3937 inch

510



PRODUCT-BENEFIT

- Dynamic Quengel splint for mobilisation of interphalangeal joints
- Two lateral spiral springs



PRODUCT-BENEFIT

- Finger brace for immobilisation of interphalangeal joints
- Adjustable aluminium splint with soft padding
- Also for children, see page 30

COLOUR: White

- THE MED-BENEFIT
- Velcro strip for optimal fixing
- Soft padding
- Progressive spring force

COLOUR: grey

THE MED-BENEFIT

- Immobilisation of either one or two fingers
- Dorsal fitting possible after remodelling
- Velcro fasteners, on metacarpals and fingers can be positioned as required

INDICATIONS

Post-operative, post-traumatic extension of PIP joints

INDICATIONS

SI7FS

Post-operative, post-traumatic, i. e. for redressement, positioning, Immobilisation of the interphalangeal joint

Ζ	ES	
	Ζ	ZES

Lenath¹⁾

	small medium large	 6,3 – 7,0 cm 7,0 – 7,6 cm 7,6 – 8,0 cm
¹⁾ Length (L)	= Distal fo finger fo	m to distal

right and left symmetrical 1 cm = 0.3937 inch

	Children Size 1 Size 2 Size 3	11 - 14 cm 13 - 16 cm 16 - 19 cm 19 - 22 cm
--	--	--

Hand length

Fingers

BORT DigiSoft® Vario-Finger Brace REF 112 790

PRODUCT-BENEFIT

- Finger brace for immobilisation of the interphalangeal joints
- Aluminium splint, adjustable with soft padding
- Velcro fasteners, variably positionable
- Can be opened fully

COLOUR: grey

THE MED-BENEFIT

Modular brace system:
 Wrist stabilisation, downgradeable



INDICATIONS

Post-operative, post-traumatic, redressement, positioning, immobilisation of the interphalangeal joint

SIZES

Metacarpus circumference

Size 2 19-22 c			17 – 19 cm 19 – 22 cm 22 – 26 cm
----------------	--	--	--



PRODUCT-BENEEIT

- Flexible wrist support to stabilise the metacarpophalangeal joint and saddle joint of the thumb using integrated reinforcing elements and additional Stabilo® Strap
- SOFTflex 3D knitted fabric with a velvety soft surface structure very comfortable to wear
- Anatomically contoured knit
- Supportive in load situations

COLOUR: black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body
- Soft, snag proof material edge
- Latex-free knitted fabric
- Easy fixation using the easyClick-System

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the area of saddle ioint of the thumb, in the case of ulnar and radial ligament instability, early stages of rhizarthrosis

SIZES	
-------	--

Circumference of wrist

r I	k-small small nedium arge k-large	up to 13 cm 13 – 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm
--------	---	---

right and left symmetrical 1 cm = 0.3937 inch



PRODUCT-BENEELT

- Thumb brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Supportive in load situations

COLOUR: black/green

THE SPORT-BENEFIT

- Breathable material provides optimal comfort
- Can be fully opened
- Can also be worn at night
- Modern, sporty design
- (with a reflecting label)

INDICATIONS

Post-operative, post-traumatic, or degenerative irritations in the area of saddle joint of the thumb, e.g arthrosis of the CMC-I joint (thumb saddle joint), post-operatively following arthroplasty of the CMC-I joint, instability of the CMC-I joint

SIZES	Circumference of wrist
H Con	x-small 13 – 15 cm small 15 – 17 cm medium 17 – 19 cm large 19 – 21 cm

please indicate left or right 1 cm = 0.3937 inch

BORT SellaFlex Thumb Support REF 112 030



PRODUCT-BENEFIT

- Flexible wrist support to stabilise the metacarpophalangeal joint and saddle joint of the thumb
- Anatomically contoured knit
- Additional stabilising elements at the thumb

BORT Thumb Hand Support REF 112 420



PRODUCT-BENEFIT

- Compression wrist support with thumb inlay to support the metacarpophalangeal joint and saddle joint of the thumb and to stabilise the wrist
- New, slim anatomically pre-moulded plastic stabilising element
- Anatomically contoured knit
- Compression class II

COLOURS: 📕 beige 🔳 black

THE MED-BENEFIT

- Soft, snag proof material edge
- Flat design, not bulky
- Very comfortable to wear

COLOURS: beige black

THE MED-BENEFIT

- Stabilo[®] Strap for additional stabilisation of the wrist
- Special lightweight and slim structure for good wearing comfort

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the area of saddle joint of the thumb, in the case of ulnar and radial ligament instability, early stages of rhizarthrosis

INDICATI	ONS
----------	-----

Chronic, post-traumatic, post-operative, irritations in the area of saddle joint of the thumb, i.e. late stage post-fracture, arthritis, ulnar and radial ligament lesions

SIZES	Circur	mference of wrist
H	x-small small medium large x-large	up to 13 cm 13 – 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm

right and left symmetrical 1 cm = 0.3937 inch

right and left symmetrical 1 cm = 0.3937 inch

Circumference of wrist

1/3	small	 up to 15	cm
H	medium	 15 – 17	ст
	large	 17 – 19	cm
	x-large	 19 – 21	cm



PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint of the thumb
- Soft inner padding
- Integrated, anatomically pre-moulded stabilising element
- Length: approx. 11 cm (4.3 inches)

COLOURS: beige grey black

THE MED-BENEFIT

- Free movement of the fingers
- Breathable material provides optimal comfort



PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Soft inner padding
- Integrated, anatomically pre-moulded stabilising element extending to the wrist
- Length: approx. 14 cm (5.5 inches)

COLOURS: beige grey

THE MED-BENEFIT

- Free movement of the fingers
- Breathable material provides optimal comfort

INDICATIONS

Post-operative, post-traumatic, degenerative, inflammatory conditions, such as in case of arthritis, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, rhizarthrosis

SIZES

Circumference of wrist

13	x-small small	_	13 – 15 cm 15 – 17 cm	
~	medium large		17 – 19 cm 19 – 21 cm	

right and left symmetrical 1 cm = 0.3937 inch

INDICATIONS

Post-operative, post-traumatic, degenerative, inflammatory conditions, such as in case of arthritis, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, **rhizarthrosis**

SIZES

Circumference of wrist

	x-small small medium large		13 – 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm
--	-------------------------------------	--	--



PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint of the thumb, saddle joint and distal interphalangeal joint of the thumb
- Soft inner padding
- Long integrated, anatomically pre-moulded stabilising element extending to the wrist area

COLOUR: grey

THE MED-BENEFIT

- Extended thumb element with Velcro fastening
- Free movement of the fingers
- Breathable material provides optimal comfort



PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Integrated aluminium stabilisation element
- Extremely thin lightweight structure
- Latex-free

COLOUR: grey

THE MED-BENEFIT

- Customisable, cold-mouldable
- Extra soft padding provides optimal comfort
- Skin-friendly velour

INDICATIONS

Post-operative, post-traumatic, i.e. following tendon reconstruction, immobilisation in the case of degenerative processes such as de Quervain's Disease, **rhizarthrosis**

INDICATIONS

Post-traumatic, post-operative, i. e. following ligament lesions, sprain, immobilisation in the case of rheumatic-inflammatory processes i. e. rhizarthrosis

Circumference of wrist

____ 15 – 17 cm

x-small _____ 13 – 15 cm

medium _____ 17 – 19 cm

SIZES

Circumference of wrist

x-small _____ 13 – 15 cm small _____ 15 – 17 cm medium _____ 17 – 19 cm large _____ 19 – 21 cm

right and left symmetrical 1 cm = 0.3937 inch large ____ 19 - 21 cm please indicate left or right

small

1 cm = 0.3937 inch

SIZES



PRODUCT-BENEFIT

- Brace for partial immobilisation of the thumb saddle joint and metacarpophalangeal joint
- Integrated aluminium stabilisation element
- Extremely lightweight structure
- Soft padding
- Special guide on the thumb tube
- Latex-free

COLOUR: 🔳 grey

THE MED-BENEFIT

- Long splint spanning the wrist joint
- Individual, cold-mouldable interlocking system
- Skin-friendly velour, comfortable to wear



PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb plus the distal interphalangeal joint
- Integrated aluminium stabilisation element
- Extremely lightweight structure
- Soft padding
- Latex-free

COLOUR: grey

THE MED-BENEFIT

- Stabilisation including the distal interphalangeal joint of the thumb
 Customisable, cold-mouldable
- Customisable, cold-mouldable
- Extra soft padding provides optimal comfort
- Skin-friendly velour

INDICATIONS

Post-traumatic, post-operative, after thumb fractures, arthritis, arthrosis of the thumb saddle and/or metacarpophalangeal joints, ulnar/radial ligament injury to the thumb, late stage after sinew reconstruction, tendopathy, for controlled mobilisation in the case of **rhizarthrosis**, ski thumb

SIZES

Circumference of wrist

H CO	x-small small medium large	 13 – 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm
	0.	

please indicate left or right 1 cm = 0.3937 inch

INDICATIONS

Post-traumatic, post-operative, i. e. following tendon reconstruction, immobilisation in the case of degenerative processes such as de Quervain's Disease, for inflammatory and rheumatic diseases, such as **rhizarthrosis**, tendopathy

SIZES

Circumference of wrist

H J	x-small small medium large		13 – 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm
-----	-------------------------------------	--	--

please indicate left or right 1 cm = 0.3937 inch

BORT Generation Thumb Splint REF 215 300



PRODUCT-BENEFIT

- Soft brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Soft inner padding
- Removable, anatomically pre-moulded stabilising elements extending to the wrist
- Length: approx. 14 cm (5.5 inches)

COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- Variably adjustable thumb strap with a reversal loop for an extra safe hold
- Suitable for use together with a walking aid/wheeled walker
- Attractive, functional design
- Can also be worn at night

INDICATIONS

Post-operative, post-traumatic, rhizarthrosis, degenerative, inflammatory, i. e. arthritis, ulnar/ radial ligament injury of the tumb



PRODUCT-BENEFIT

- Thumb ring brace for the stabilisation and relief of the thumb saddle joint
- Positioning the thumb in functional position
- Integrated, cold adjustable stabilisation element
- Removable fastening strap
- Optimal wearing comfort
- Easy to put on

COLOUR: grey

THE MED-BENEFIT

- Extremely slim and functional brace design
- Free movement of the fingers and the wrist
- Washable, highly durable

INDICATIONS

Post-operative, post-traumatic or degenerative irritations in the area of saddle joint of the thumb, e.g arthrosis of the CMC-I joint (thumb saddle joint), post-operatively following arthroplasty of the CMC-I joint, instability of the CMC-I joint

SIZES

Circumference of wrist

small _____ 13 – 16 cm medium _____ 16 – 19 cm large _____ 19 – 22 cm

please indicate left or right 1 cm = 0.3937 inch

SIZES

Metacarpus circumference

P	Size 0	15 – 17 cm
	Size 1	17 – 19 cm
	Size 2	19 – 22 cm
	Size 3	22 – 26 cm



PRODUCT-BENEFIT

- Thermoplastic brace for immobilisation of the metacarpophalangeal joint of the thumb
- Inner side is completely fitted with soft foam
- Circumferential Velcro fastening
- Extensive perforation

COLOUR: grey

THE MED-BENEFIT

- Free movement of the wrist (short design)
- Thermoplastically mouldable (warm air)
- Individually adjustable
- Comfortable to wear



PRODUCT-BENEFIT

- Thermoplastic brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Inner side is completely fitted with soft foam
- Circumferential Velcro fastening
- Extensive perforation

COLOUR: grey

THE MED-BENEFIT

- Thermoplastically mouldable (warm air)
- Individually adjustable
- Comfortable to wear

INDICATIONS

Post-traumatic, post-operative, sprain, contusion

INDICATIONS

SIZES

Post-traumatic, post-operative, i. e. sprain, contusion, immobilisation in the case of rheumatic-inflammatory processes (rhizarthrosis), ligament lesions (skiers thumb)

> small medium

large

x-large

Circumference of wrist x-small ____ up to 15 cm small ____ 15 - 17 cm

> ___ 17 – 19 cm ___ 19 – 21 cm

21 – 23 cm

SIZES

Circumference of wrist

	x-small small medium large x-large		up to 15 cm 15 - 17 cm 17 - 19 cm 19 - 21 cm 21 - 23 cm
--	--	--	---

please indicate left or right 1 cm = 0.3937 inch

BORT SellaFix® N Black Edition REF 112 040 B



PRODUCT-BENEFIT

- Thermoplastic brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Innovative fastening system for optimum adaptation to the patient's circumstances
- Slim structure

COLOUR: 🗖 black

THE MED-BENEFIT

- Individually adjustable
- Thermoplastically mouldable (warm air)
- Comfortable to wear due to the use of innovative materials
- Washable

BORT SellaDur® REF 105 400



PRODUCT-BENEFIT

- Brace for immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Special raised edge in the critical thumb forefinger area (commissure)
- Elastic sleeve with downward traction

COLOUR: anthracite

THE MED-BENEFIT

- Extremely durable brace
- Extensive perforation for comfortable wear
- Washable

INDICATIONS

Post-traumatic, post-operative, i.e. sprain, contusion, immobilisation in the case of rheumatic-inflammatory processes (i.e. rhizarthrosis), ligament lesions (skiers thumb)

A	x-small small medium	 18 – 20 cm 20 – 22 cm 22 – 24 cm
	large	 24 – 26 cm
	x-large	 26 – 28 cm


please indicate left or right 1 cm = 0.3937 inch

INDICATIONS

Post-traumatic, post-operative, i. e. sprain, contusion, immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis, ligament lesions (skiers thumb)

SIZES

Circumference of wrist





PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Soft inner padding
- Splints individually adjustable and removable

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- The specially padded zone in the palm area enables use combined with walking aids/a wheeled walker
- Variably adjustable thumb wrap with a reversal loop for an extra safe hold

COLOUR: 🔳 blue



BORT Generation: For better mobility for the aged Product range with special handling advantages

INDICATIONS

Pre-operative, post-traumatic, post-operative, i. e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, i. e. rhizarthrosis, in the case of degenerative diseases such as carpal tunnel syndrome

SIZES

Circumference of wrist

H A	x-small small medium large	up to 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm
	x-large	 _ 21 – 23 cm

THUMB 77

BORT SellaTex® Plus – Quervain REF 112 160



PRODUCT-BENEFIT

- Wrist brace for stabilisation of the metacarpophalangeal joint, saddle joint and distal interphalangeal joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Can be fully opened
- Length approx. 20 cm (7.9 inches)

COLOUR: grey

THE MED-BENEFIT

- Large opening in the region of the problematic medial area in case of Tendovaginitis de Quervain
- Easy downgrade option for the thumb splint within the course of therapy
- Long thumb element can be individually positioned, can be shaped, fully removable
- Pleasant wearing comfort



INDICATIONS

Tendovaginitis de Quervain, pre-operative, post-traumatic, post-operative, i.e. after thumb ligament reconstruction, immobilisation in the case of rheumaticinflammatory processes, in the case of degenerative diseases, i. e. rhizarthrosis, carpal tunnel syndrome

SIZES

Circumference of wrist

1/3	x-small	up to 15 cm
-HB	small	 15 – 17 cm
~	medium	 17 – 19 cm
	large	19 – 21 cm
	x-large	 21 – 23 cm



PRODUCT-BENEFIT

- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Dorsal thumb splint for additional thumb stabilisation
- Length: approx. 20 cm (7.9 inches)
- Also for children, see page 30

COLOUR: grey

THE MED-BENEFIT

- Stabilo[®] Strap for additional stabilisation of the wrist
- Extremely flat design
- Breathable material
- Comfortable to wear

INDICATIONS

Pre-operative, post-traumatic, post-operative, i.e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis, in the case of degenerative diseases such as de Quervain's Disease, carpal tunnel syndrome

SIZES	Circumference of wrist
Children x-small small medium large x-large	13 – 15 cm up to 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm 21 – 23 cm



- Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb
- Palm and dorsal splints to stabilise wrist and metacarpus
- Can be fully opened
- Length: approx. 19 cm (7.5 inches)

THE MED-BENEFIT

- Adjustable thumb element, which can be individually positioned, completely removable
- Comfortable to wear
- Stabilo[®] Strap

COLOUR: 🔳 grey

INDICATIONS

Pre-operative, post-traumatic, post-operative, i.e. after thumb ligament reconstruction, immobilisation in the case of rheumatic-inflammatory processes, rhizarthrosis, in the case of degenerative diseases such as carpal tunnel syndrome

SIZES		Circumference of wrist
x-sm smal med large x-lar	um	up to 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm 21 – 23 cm

All in hand: More functional stability for the wrist.



Overview of indications for selected products	select* ManuZip® volar 112 220	ManuBasic [®] Sport 112 920 SP	ManuBasic [®] 112 920 ManuBasic [®] Plus 112 940	ManuBasic* 112 920Wrist Support with Alu-ManuStabil* 102 200ManuBasic* Plus103 300Classic 103 300112 940Generation Wrist Support102 200 C 215 000	ManuStabije 102 200 ManuStabije Classic 102 200 C ManuStabije Short 112 730	SOFT Volar CTS Splint 112 750	Generation Wrist Brace 215 100 Generation Wrist Brace 215 150	ManuCarpal [®] Combi 102 400	Dorsal Intrinsic Plus Splint (D.I.P.S) 102 500	ManuDisc [®] 102 300
Page	P. 82	P. 82	P. 83/84	P. 86 / 87	P. 90/91	P. 92	P. 88	P. 94	P. 95	P. 93
Irritations	•	•	•	•						
Inflammatory rheumatic							•			
Degenerative				•			•			
Sprain	•	•	•	•		•				
Carpal arthrosis	•	•	•	•			•			•
Carpal tunnel syndrome (CTS)				•	•	•	•	•		•
Tenosynovitis	•	•	•	•	•		•	•		
Disc injury										•
(Post-operative) immobilisation					•	•		•		•
Peripheral nerve paralysis					•			•		
Metacarpal and metacarpo- phalangeal joint fractures									٠	
Ligament injuries of the long fingers									•	

_____ OVERVIEW OF INDICATIONS 81

Hand



PRODUCT-BENEFIT

- High quality compression wrist support to stabilise the wrist
- Aluminium palm splint, anatomically premoulded
- Stabilo® Strap, freely positionable
- Viscoelastic silicone pad for intermittent massage
- Anatomically contoured knit

COLOURS: silver black

THE SELECT-BENEFIT

- Zip in the region of the back of the hand, therefore very comfortable to put on and off
- Extra-soft, snag proof material edge
- Integrated slip-on aid

BORT ManuBasic[®] Sport REF 112 920 SP



PRODUCT-BENEFIT

- Compressing wrist support with additional Stabilo[®] Strap
- SOFTflex 3D knitted fabric with a velvety soft surface structure
- Anatomically contoured knit
- Plastic palm splint, anatomically preformed, thermoplastically mouldable
- Viscoelastic silicone pad

COLOUR: I black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body, for dryness and comfort
- Ideal fit for sports activities
- Latex-free knitted fabric
- Extra elastic in the wrist and forearm area, making it very easy to put on
- Easy fixing and length setting of the Stabilo[®] Strap using the easyClick-System

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

SIZES	Circumference of wrist	SIZES	Circumference of wrist
n	small 13 - 15 cm medium 15 - 17 cm large 17 - 19 cm x-large 19 - 21 cm	H	x-small up to 13 cm small 13 - 15 cm medium 15 - 17 cm large 17 - 19 cm x-large 19 - 21 cm

please indicate left or right 1 cm = 0.3937 inch



- Compression wrist support to stabilise the wrist
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Plastic palm splint, anatomically preformed, thermoplastically mouldable
- Stabilo® Strap, freely positionable
- Silicone pad with hollow at the ulnar styloid process

COLOURS: beige silver

THE BASIC-BENEFIT

- Extra elastic in the wrist and forearm area, making it easier to put on
- Especially soft material edge in the forearm area
- Plus-size system results in a good fit

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

SIZES		c1-measure	c2-measure	c1-measure c2-measure
C2 C11 D28 S S r r	small small plus medium medium plus	13 - 15 13 - 15 15 - 17	18 - 20 cm 20 - 22 cm	x-large 19 - 21 22 - 24 cm x-large plus 19 - 21 24 - 26 cm

please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Compressing wrist support with additional stabilisation of the metacarpophalangeal joint and saddle joint of the thumb
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Anatomically pre-moulded, thermoplastic palm splint
- Silicone pad with hollow at the ulnar styloid process

COLOUR: 🗖 black

THE BASIC-BENEFIT

- Velvety soft surface structure
- Extra-soft, snag proof material edge
- Stabilo[®] Strap for additional stabilisation of the wrist



PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- Palm splint, anatomically pre-formed, thermoplastically mouldable
- Stabilo[®] Strap
- Silicone pad
- Anatomically contoured knit

COLOUR: Dlack

THE ECO-BENEFIT

- Zip in the region of the back of the hand, therefore very comfortable to put on and off
- Integrated slip-on aid

INDICATIONS

Chronic post-traumatic or post-operative irritation in the wrist and thumb saddle joint areas, sprain, carpal arthrosis, tenosynovitis, in the case of ulnar and radial ligament instability, early stages of rhizarthrosis

SILES	S	IZES
-------	---	------

Circumference of wrist

13	small	 13 – 15 cm
	medium	 15 – 17 cm
<u>/</u>	large	 17 – 19 cm
	x-large	 19–21 cm

please indicate left or right 1 cm = 0.3937 inch

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e.g. sprain, carpal arthrosis, tenosynovitis

SIZES	Circ	
1/3	small	

Circumference of wrist

- H	small medium	 13 – 15 cm 15 – 17 cm
	large	 17 – 19 cm
	x-large	 19 – 21 cm



- Compression wrist support to stabilise the wrist
- Palm splint, anatomically pre-formed, thermoplastically mouldable
- Stabilo® Strap
- Silicone pad with hollow at the ulnar styloid process

BORT Ganglion Support REF 112 820

PRODUCT-BENEFIT

- Compression wrist support to stabilise the wrist
- Palm splint, anatomically pre-formed, thermoplastically mouldable
- Stabilo[®] Strap, for additonal stabilisation of the wrist
- Silicone pad with hollow at the ulnar styloid process

COLOURS: beige mineral grey graphite

THE MED-BENEFIT

 Seamless expansion zone, perfectly surrounding the base of the thumb and optimally adjusting to anatomical form

COLOUR: beige

THE MED-BENEFIT

- Ganglion pad, freely positionable

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, tenosynovitis

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, intracapsular ganglion

SIZES

Circumference of wrist

small _____ 13 – 15 cm medium _____ 15 – 17 cm large _____ 17 – 19 cm x-large _____ 19 – 21 cm

please indicate left or right 1 cm = 0.3937 inch

SIZES

Circumference of wrist

|--|



PRODUCT-BENEFIT

- Classical wrist support
- Individually adjustable aluminium palm splint, in a durable pocket
- Length: 16 cm (6.3 inches)

COLOURS: beige grey

THE MED-BENEFIT

 Integrated slip-on aid and a Velcro strip on the back of the hand make it easy to apply/take off, even if gripping ability is restricted

BORT Wrist Support with Aluminium Splint and Additional Strap REF 103 301

- Version as REF 103 300

COLOURS: beige grey

 Additional stabilisation thanks to Stabilo[®] Strap

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, carpal tunnel syndrome (CTS), rheumatism, tenosynovitis

SIZES

Circumference of wrist

n	x-small small medium large	 up to 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm
	x-large	 13 - 21 cm $21 - 23$ cm

_HAND 87



PRODUCT-BENEFIT

- Padded wrist support with palm reinforcement
- Elastic basicmaterial
- Palm splint, anatomically pre-formed, thermoplastically mouldable

COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- The specially padded zone in the palm area enables use combine with a walking aid/wheeled walker
- Attractive, functional design



BORT Generation: For better mobility for the aged Product range with special handling advantages

INDICATIONS

Chronic, post-traumatic or post-operative irritations in the wrist area, e.g. distorsion, carpal arthrosis, carpal tunnel syndrome (CTS), rheumatism, tenosynovitis

SIZES

Circumference of wrist

H

, small ______ 13 – 16 cm medium ______ 16 – 19 cm large ______ 19 – 23 cm



PRODUCT-BENEFIT

- Padded wrist brace with palm and dorsal reinforcement
- Splints individually adjustable and removable
- Soft inner padding

COLOUR: blue

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- The specially padded zone in the palm area enables use combine with a walking aid / wheeled walker
- Attractive, functional design
- Can also be worn at night



PRODUCT-BENEELT

- Padded wrist brace with palm and dorsal reinforcement
- Splints individually adjustable and removable
- Breathable material provides optimum comfort

COLOUR: blue

THE GENERATION-BENEFIT

- Can be put on with one hand due to innovative strapping fastening
- The specially padded zone in the palm area enables use combine with a walking aid / wheeled walker
- Attractive, functional design
- Can also be worn at night


INDICATIONS

Pre-operative, post-operative, post-traumatic, carpal arthrosis, inflammatoryrheumatic, degenerative diseases, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

SIZES

SIZES	Circumference of wrist		
R CON	x-small up to 15 cm small 15 - 17 cm medium 17 - 19 cm large 19 - 21 cm x-large 21 - 23 cm		

please indicate left or right 1 cm = 0.3937 inch



Circumference of wrist

BORT Arm and Wrist Support [REF] 103 360

- Classical arm and wrist support to stabilise the wrist
- Individually adjustable aluminium palm splint in durable pocket
- Velcro fastening on back of hand
- Length: 22 cm (8.7 inches)

COLOURS: beige grey

THE MED-BENEFIT

- Longer aluminium splint resulting in better stabilisation
- Easy to put on/remove, even if gripping ability is restricted, thanks to integrated slip-on aid and dorsal Velcro fastening



PRODUCT-BENEFIT

- Classical arm and wrist support to stabilise the wrist
- Individually adjustable aluminium palm splint in durable pocket
- Velcro fastening on back of hand
- Length: 22 cm (8.7 inches)

COLOURS: beige grey


THE MED-BENEFIT

- Additional stabilisation thanks to the Stabilo[®] Strap
- Longer aluminium splint resulting in better stabilisation
- Easy to put on/remove, even if gripping ability is restricted, thanks to integrated slip-on aid and dorsal Velcro fastening

INDICATIONS

1 cm = 0.3937 inch

Chronic, post-traumatic or post-operative irritations in the wrist area, e. g. sprain, carpal arthrosis, carpal tunnel syndrome (CTS), tenosynovitis





PRODUCT-BENEFIT

- Rigid wrist brace with palm and dorsal reinforcement splints
- Breathable material
- Right-hand and left-hand versions
- Long design
- Length: approx 22 cm (8.7 inches)

COLOUR: grey

THE MED-BENEFIT

- Easy to shape individually
- Excellent stabilisation of the wrist thanks to extended length and special fastening technology
- Can also be worn at night



PRODUCT-BENEFIT

- Rigid wrist brace with palm and dorsal reinforcement splints
- Breathable material
- Can be fully to opened
- Can be used on the left or right
- Long design
- Length: approx 22 cm (8.7 inches)


COLOUR: grey

THE MED-BENEFIT

- Individually mouldable
- Circular strap for effective stabilisation of the wrist
- Can also be worn at night

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis



please indicate left or right 1 cm = 0.3937 inch right and left symmetrical 1 cm = 0.3937 inch



- Rigid wrist brace with palm and dorsal reinforcement splints, individually shaped
- Short design
- Length: approx. 18 cm (7.1 inches)
- Also for children, see page 31

COLOURS: beige grey

THE MED-BENEFIT

- Stabilo[®] Strap for effective stabilisation of the wrist
- Breathable material provides optimal comfort
- Can also be worn at night



PRODUCT-BENEFIT

- Wrist brace with palm and dorsal reinforcing splint
- Breathable, skin-friendly material
- Can be fully to opened
- Extra-long design
- Length: approx 25 cm (9.8 inches)

COLOUR: grey

THE MED-BENEFIT

- Easy to shape individually
- Stabilisation of the wrist and forearm
- Easy to put on, even following surgery
- Can also be worn at night

INDICATIONS

Pre-operative, post-operative, posttraumatic, degenerative, carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

SIZES	Circumference of wrist		
	Children (only grey) x-small small medium large x-large	13 - 15 cm up to 15 cm 15 - 17 cm 17 - 19 cm 19 - 21 cm 21 - 23 cm	


please indicate left or right 1 cm = 0.3937 inch

INDICATIONS

Post-operative, post-traumatic, degenerative, e. g. carpal tunnel syndrome (CTS), tenosynovitis

SIZES

Circumference of wrist





PRODUCT-BENEFIT

- Soft palm rest splint for immobilisation of wrist
- Extra soft padded splint and fastening straps
- Latex-free

COLOUR: grey

THE MED-BENEFIT

- Easily mouldable
- Lightweight and comfortable to wear
- Can also be worn at night

BORT SOFT Dorsal CTS Splint REF 112 760

PRODUCT-BENEFIT

- Soft dorsal rest splint for immobilising the wrist
- Palm splint
- Extra soft padded splint and fastening straps
- Three point Velcro fastening
- Latex-free

COLOUR: grey

THE MED-BENEFIT


- Maintained gripping ability
- Easily mouldable
- Lightweight and comfortable to wear
- Can also be worn at night

INDICATIONS

Post-operative, post-traumatic, e. g. carpal tunnel syndrome (CTS), sprain

SIZES	

Circumference of wrist



please indicate left or right 1 cm = 0.3937 inch SIZES

Circumference of wrist



BORT ManuCarpal® Brace REF 103 200



PRODUCT-BENEFIT

- Thermoplastically adjustable planthand brace
- Full padding
- Large contact surface with bilateral support of the forearm
- Extra lightweight design
- Washable

COLOUR: 🗖 black

THE MED-BENEFIT

- Thermoplastically mouldable at 90°C (warm air)
- Easy for the patient to position and use
- Can be fully opened
- Comfortable to wear

BORT ManuDisc® acc. to Dr. Schütz REF 102 300

PRODUCT-BENEFIT

- Wrist brace for multidirectional stabilisation of wrist
- Ulnar support and bilateral forearm rest
- Individually adjustable
- Soft padding on all inside surfaces

COLOUR: grey

THE MED-BENEFIT

- Lightweight design with T-shaped palm splint
- Can be fully opened
- Comfortable to wear

INDICATIONS

Pre-operative, post-operative, posttraumatic, degenerative, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

SIZES

Circumference of wrist

 small medium large	 15 – 17 cm 17 – 19 cm 19 – 21 cm
x-large	 21 – 23 cm


please indicate left or right 1 cm = 0.3937 inch

INDICATIONS

Pre-operative, post-operative, posttraumatic, e. g. carpal tunnel syndrome (CTS), disc injury, carpal arthrosis

SIZES

Circumference of wrist





PRODUCT-BENEFIT

- Hand brace with removable finger rest
- Individually adjustable finger rest
- Can be fully opened
- Can be used on the left or right

COLOUR: grey

THE MED-BENEFIT

- Modular orthotic system
- Additional relief in the carpal tunnel when using the finger rest
- Circular strap for effective stabilisation of the wrist



- Functionally shaped relief pad consisting of memory foam
- Relief of the carpal tunnel and the sensitive scar region (condition after surgery)
- Comfortable to wear

INDICATIONS

Pre-operative, post-operative, post-traumatic, degenerative, e. g. carpal tunnel syndrome (CTS), tenosynovitis, peripheral nerve paralysis

SIZES

Circumference of wrist

1/3	x-small	 up to 15 cm
-H -B	small	 15 – 17 cm
~	medium	 17 – 19 cm
	large	 19 – 21 cm
	x-large	 $21-23\ \text{cm}$

BORT Dorsal Intrinsic Plus Splint (D.I.P.S.) REF 102 500



PRODUCT-BENEFIT

- Wrist splint for immobilisation with finger fixation
- Ideal pre-forming of the splint in Intrinsic Plus position
 - · 25° wrist extension
 - · 70° metacarpophalangeal flexion
- Two elastic Velcro fasteners in the lower arm region, respectively a non-elastic Velcro fastening in the metacarpophalangeal joint area of the hand as well as the phalanx area
- Continuously softly padded on both sides, ideally comfortable to wear
- Easy post-operative use due to anatomical pre-forming

THE MED-BENEFIT

- Immobilising dorsal attachment
- Post-operative fast and safe attachment
- Avoidance of contractures of the long fingers
- Early mobilisation of the long fingers possible via fist closure
- Stretching the long fingers in the central and interphalangeal joints with Velcro for wearing at night-time

COLOUR: 🗖 black

INDICATIONS

Pre-operative, post-operative, post-traumatic, e. g. metacarpus and proximal finger phalanx in the area of the long fingers (e.g. a boxer fracture), ligament lesions in the long fingers in the proximiti of the metacarpophalangeal joints, injuries and diseases which make the immobilisation of long fingers necessary (for example infections)

SIZES

Hand length

$\overline{\mathcal{A}}$	Size 1	14 – 17 cm
	Size 2	17 – 22 cm

Hand



PRODUCT-BENEFIT

- Palm splint for immobilising the wrist while fingers are blocked
- Long, one-piece aluminium splint, individually adjustable
- Integrated slip-on aid
- Breathable material provides optimal comfort

COLOUR: grey

THE MED-BENEFIT

- Easy to put on, even if gripping ability is restricted
- Can also be worn at night



PRODUCT-BENEFIT

- Soft palm splint for immobilising the wrist while fingers and thumb are blocked
- Wide finger rest and thumb surround, individually adjustable
- Four Velcro fastenings to block the splint securely
- Extra-soft padding material
- Latex-free

COLOUR: grey

THE MED-BENEFIT

- Simple immobilisation of the thumb
- Safe and comfortable at day and night

INDICATIONS

Pre-operative, post-operative, posttraumatic, degenerative, e.g. carpal tunnel syndrome (CTS), tenosynovitis

INDICATIONS

Pre-operative, post-operative, posttraumatic, e. g. carpal tunnel syndrome (CTS), paralysis, for inflammatory and rheumatic diseases

SIZES

Circumference of wrist

н	x-small small medium large	 up to 15 cm 15 – 17 cm 17 – 19 cm 19 – 21 cm 21 – 23 cm
	x-large	 21 – 23 cm

please indicate left or right 1 cm = 0.3937 inch SIZES

Circumference of wrist

13	x-small	 13 – 15 cm
H	small medium	 15 – 17 cm 17 – 19 cm
	large	 19 – 21 cm

BORT Volar Resting Splint REF 103 100



PRODUCT-BENEFIT

- Palm rest splint for immobilising the wrist while fingers and thumb are blocked
- Adjustable stabilising element
- Soft foam padding





PRODUCT-BENEFIT

- Palm rest splint for immobilising the wrist while fingers and thumb are blocked
- Thermoplastic material
- Inside fully lined with textile flocking
- Washable

COLOUR: COLOUR:

THE MED-BENEFIT

- Individually adjustable without using tools
- Extra lightweight structure
- Lower arm hand rest with skinfriendly textile cover

COLOUR: beige

THE MED-BENEFIT

- For positioning the hand in the intrinsic-plus position
- Easy, individual adjustment using a hair-dryer at 60°C
- Can be applied by the patient

Pre-operative, post-operative, post-

INDICATIONS

Pre-operative, post-operative, posttraumatic, degenerative, e. g. carpal tunnel syndrome (CTS), for paralysis, for inflammatory and rheumatic diseases, sprain

SIZES

Circumference of wrist


traumatic, e. g. carpal tunnel syndrome (CTS), hand paresis, for inflammatory

and rheumatic diseases


INDICATIONS

SIZES

Circumference of wrist



please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Compression wrist support for easy stabilisation of the wrist
- Optimised fit in the palm area
- Anatomically contoured knit
- Hemmed thumb opening

COLOURS: beige blue

THE MED-BENEFIT

- Easy to use
- Can be worn on the left or right



PRODUCT-BENEFIT

- Compression wrist support for stabilisation of the wrist and metacarpus
- Accessory: Removable and freely positionable Stabilo[®] Strap

COLOURS: beige blue

THE MED-BENEFIT

- Easy to use
- Can be worn on the left or right

INDICATIONS

Minor sprains, strains, after removal of plaster cast, tenosynovitis

SIZES

Circumference of wrist

x-small ____ up to 13 cm small ____ 13 - 15 cm medium ____ 15 - 17 cm large ____ 17 - 19 cm x-large ____ 19 - 21 cm

right and left symmetrical 1 cm = 0.3937 inch
 small
 13 – 15 cm

 medium
 15 – 17 cm

 large
 17 – 19 cm

 x-large
 19 – 21 cm

Circumference of wrist

x-small ____ up to 13 cm

right and left symmetrical 1 cm = 0.3937 inch

SIZES



- Elastic wrist support
- Velcro fastening
- Width: 8 cm (3.2 inches)



THE MED-BENEFIT

- Easy to put on



PRODUCT-BENEFIT

- Elastic wrist support with Stabilo[®] Strap
- Velcro fastening
- Width: 8 cm (3.2 inches)

COLOURS: beige black

THE MED-BENEFIT

- Additional stabilisation thanks to the Stabilo[®] Strap
- Easy to put on

USE

Gentle Wrist Support

SIZES	Circumference of wrist	SIZES	Circumference of wrist
H	Size 1 up to 17 cm Size 2 17 - 19 cm Size 3 19 - 22 cm	8	Size 1 up to 17 cm Size 2 17 – 19 cm Size 3 19 – 22 cm

right and left symmetrical 1 cm = 0.3937 inch right and left symmetrical 1 cm = 0.3937 inch

For a strong back: Correction, stabilisation, relief.



104 700 180 500 180 100 104 730 104 760 P. 103 P. 113 P. 113 P. 114 P. 108 P. 0 P. 113 P. 117 P. 114 P. 108 P. 0 P. 113 P. 117 P. 114 P. 108 P. 0 P. 113 P. 116 P. 114 P. 108 P. 0 P. 116 P. 116 P. 114 P. 108 P. 0 P. 113 P. 116 P. 114 P. 108 P. 10 P. 116 P. 116 P. 116 P. 108 P. 10 P. 116 P. 116 P. 116 P. 108 P. 10 P. 116 P. 116 P. 116 P. 108 P. 10 P. 116 P. 116 P. 116 P. 108 P. 10 P. 116 P. 116 P. 116 P. 108 P. 10 P. 116 P. 116 P. 116 P. 117 P. 10 P. 116 P. 116 P. 116 P. 116 P. 10 P. 116 P. 116 P. 116	Overview of indications for selected products	Corset Sport Sport 104 660 SP 1	select * Stabilo ® 104 710 104 740	Stabilo- Basic Sport 104 690 SP 104 680 SP	Generation Back Support 215 850	Stabilo [®] Sacroiliac Brace 104 750	Vario-Plus special Width 112 690	Lumbo Xpress Back Support	select [*] Lumbar Spine Brace	Generation Stabilo [®] Back FA Brace Back 215 800 Brace		Stabilo [®] Lumbar Back Brace	Stabilo- Plus Back Support	Stabilo- Plus Bridging Brace	Generation Osteoporosis Brace 215 870
P.103 P.105 P.104 P.103 P.117 P.116 P.116 P.114 P.106 cic spine postue cic spine postue eic spine postue				Stabilo Basic special Width Padded Back Support 104 670				112 600	104 700			180 100			Generation Dynamic Osteo- porosis Brace 215 880
Thoracic spine postue •	Page	P. 103	P. 105	P. 106 / 108	P. 113	P. 112	P. 109	P. 111	P. 103	P. 113		P. 116	P. 114	P. 108	P. 118/119
ISG (sacrolitac joint) syndrome I <t< td=""><td>Thoracic spine posture correction</td><td>•</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>	Thoracic spine posture correction	•													
Immago Immago<	ISG (sacroiliac joint) syndrome					•									
Spondylose Imbar spine syndrome Imbar spine syndrom	Lumbago		•	•	•		•								
Immar spine syndrome Immar spi	Spondylose		•				•	•		•					
Low back pain Image: Schedung sched	Lumbar spine syndrome		•				•	•							
Ischieldia Image: Schieldia	Low back pain			•	•			•		•					
Eacet joint arthrosis	Ischialgia			•						•					
Prolaps (consentive/ post-operative) •	Facet joint arthrosis								•		•				
Osteoporosis 0 Thoracic spine fractures • Degenerative, lumbar change •	Prolaps (conservative/ post-operative)								•	•		•	•	•	
Thoracic spine fractures • • • Degenerative, lumbar change • • •	Osteoporosis														•
Degenerative, lumbar change	Thoracic spine fractures											•		•	
	Degenerative, lumbar change				•					•					



PRODUCT-BENEFIT

- Corset for correcting the posture of the thoracic spine
- Soft, non-elastic strap material
- For children too, see page 29



PRODUCT-BENEFIT

- Elastic corset for **active** straightening of the thoracic spine
- Breathable, skin-friendly material
- Slim design, not bulky

COLOURS: blue white

THE MED-BENEFIT

- Fastening in the chest area makes it possible to put on/remove the item independently
- Moulds itself perfectly to the body contours, especially in the critical underarm area
- Latex-free

INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine


COLOURS: □ white ■ grey

THE MED-BENEFIT

- Individual adjustment thanks to flexible elastic belts
- Easy to put on, thanks to the rucksack principle

INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine, back pain (pain in the region of the thoracic spine)



THORACIC SPINE/BACK 103

BORT Corset Sport REF 104 660 SP



PRODUCT-BENEFIT

- Corset for straightening and relief of the thoracic spine
- Straightening according to the rucksack principle
- Really easy to put on and off

COLOUR: I black/green

THE SPORT-BENEFIT

- Free shoulder blade movement due to the slim design
- Conforms perfectly to the body contours, especially in the critical armpit area
- Modern, sporty design (with a reflecting label)

INDICATIONS

Straightening and relief of the thoracic spine, e. g. misalignment of thoracic spine, back pain (pain in the region of the thoracic spine)

BORT select* Lumbar Spine Brace for Mobilisation REF 104 700

PRODUCT-BENEFIT

- Back brace for functional mobilisation of the lumbar spine
- Two anatomically pre-moulded splints in the dorsal area
- Easy to dose support and relief effect via inelastic belt tensioning
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

COLOUR: silver

THE SELECT-BENEFIT


- Removable restraint for adaptation to the course of therapy
- Easy to put on using the hand straps provided
- Easier to put on using the ergonomic latch plate

INDICATIONS

Moderate lumboischialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade l spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history), i. e. after intervertebral disc surgery with mild to moderate prolapse, facet joint arthrosis

SIZES	Underbust	circumference
	Size 1 Size 2 Size 3 Size 4	65 – 80 cm 80 – 95 cm 95 – 110 cm 110 – 125 cm

SIZES Circumference of waist





PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomicallycontoured knit
- Four supporting splints
- High-quality viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches). dorsal 23 cm (9.1 inches)





THE SELECT-BENEELT

- Hand straps allow support to be fitted easily and without effort
- Easier to put on using the ergonomic latch plate
- Comfortable to wear



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomicallycontoured knit
- Fitted at the waist: specific for women
- Four supporting splints
- High-quality viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches). dorsal 24 cm (9.4 inches)

COLOURS: silver black



THE SELECT-BENEELT

- Hand straps allow support to be fitted easily and without effort
- Easier to put on using the ergonomic latch plate
- Comfortable to wear

INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

SIZES

SIZES	Circumference of waist
(Waist)	Size 1 70 - 80 cm Size 2 80 - 90 cm Size 3 90 - 100 cm Size 4 100 - 110 cm Size 5 110 - 120 cm Size 6 120 - 130 cm Size 7 130 - 140 cm

1 cm = 0.3937 inch

1 cm = 0.3937 inch

Circumference of waist

(Waist)		 80 - 90 cm 90 - 100 cm 100 - 110 cm 110 - 120 cm
	Size 5 Size 6	



- Back brace to stabilise the lumbar spine
- Three anatomically pre-moulded splints in the dorsal area
- Viscoelastic lumbo-sacral pad with cross strap system
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

BORT select* Stabilo® Lady Back Support with Pad REF 104 740

PRODUCT-BENEFIT

- Back brace to stabilise the lumbar spine
- Fitted at the waist: specific for women
- Three anatomically pre-moulded splints in the dorsal area
- Viscoelastic lumbo-sacral pad with cross strap system
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)

COLOURS:	silver black



THE SELECT-BENEFIT

- Cross strap system for additional stabilisation and relief
- Easier to put on using the ergonomic latch plate
- COLOURS: silver ■ black




act,

THE SELECT-BENEFIT

- Cross strap system for additional stabilisation and relief
- Easier to put on using the ergonomic latch plate

INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, lumbosciatica, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery



1 cm = 0.3937 inch



PRODUCT-BENEEIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomicallycontoured knit
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)

COLOUR: black/green



THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body, for dryness and comfort
- Latex-free knitted fabric
- Practical hand straps allow support to be fitted easily and without effort

BORT StabiloBasic Lady Sport Back Support with Pad **RFF** 104 680 SP



PRODUCT-BENEELT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomicallycontoured knit
- Fitted at the waist: specific for women
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 24 cm (9.4 inches)

COLOUR: I black/green



THE SPORT-BENEFIT

- COOLMAX®: wicks moisture away from the body, for dryness and comfort - Latex-free knitted fabric
- Practical hand straps allow support to be fitted easily and without effort

INDICATIONS

01750

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), ischialgia, dorsalgia

SIZES	Circumterence of waist
(Waist)	Size 1 70 - 80 cm Size 2 80 - 90 cm Size 3 90 - 100 cm Size 4 100 - 110 cm Size 5 110 - 120 cm Size 6 120 - 130 cm Size 7 130 - 140 cm

$\begin{tabular}{ c c c c c c c c c c c c c c c c c c c$	SIZES	Circumference of waist
	(Waist)	$\begin{array}{cccccccccccccccccccccccccccccccccccc$

1 cm = 0.3937 inch

1 cm = 0.3937 inch

BORT Stabilo Basic Back Support with Pad REF 104 690



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomicallycontoured knit
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)

BORT Stabilo Basic Lady Back Support with Pad REF 104 680



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Active support with anatomically-contoured knit
- Cut fitted at the waist: especially for women
- Four supporting splints
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 24 cm (9.4 inches)



THE BASIC-BENEFIT

- Hand straps allow support to be fitted easily and without effort
- Easier to put on using the ergonomic latch plate
- Comfortable to wear

COLOURS: beige silver black/lilac saral pad THE BASIC-BENEFIT

- Hand straps allow support to be fitted easily and without effort
- Easier to put on using the ergonomic latch plate
- Comfortable to wear

INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

SIZES	Circumference of waist	SIZES	Circumference of waist
(Waist)	Size 1 70 - 80 cm Size 2 80 - 90 cm Size 3 90 - 100 cm Size 4 100 - 110 cm Size 5 110 - 120 cm Size 6 120 - 130 cm Size 7 130 - 140 cm	(Waist)	Size 0 60 - 70 cm Size 1 70 - 80 cm Size 2 80 - 90 cm Size 3 90 - 100 cm Size 4 100 - 110 cm Size 5 110 - 120 cm Size 6 120 - 130 cm
4 0.00	07: /	4 0.00	07: /

1 cm = 0.3937 inch

1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Anatomically contoured knit
- Viscoelastic sacral pad, which can be freely positioned
- Very flat construction, therefor hardly visible under clothing
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)



COLOUR: silver

THE BASIC-BENEFIT

- Additional elastic straps to lift the abdomen for relief of the lumbar spine for larger sizes
- Specially constructed fastening system simplifies putting on
- Breathable fabric



INDICATIONS

Pain in the lumbar spine, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis), dorsalgia

SIZES

Circumference of torso

Greener Size 1 Size 2	110 - 120 cm 120 - 130 cm 130 - 140 cm 140 - 150 cm 150 - 160 cm
---	--

BORT VarioPlus Back Support special Width REF 112 690



PRODUCT-BENEFIT

- Non-elastic back support for larger sizes
- Semi-elastic body belt with lateral Velcro fastening
- Fitted cut in the lumbar spinal area
- Lumbo-sacral friction pad
- Three supporting splints in the dorsal area
- Height: ventral 24 cm (9.4 inches) dorsal 30 cm (11.8 inches, sizes 1+2) ventral 28 cm (11.0 inches) dorsal 33 cm (13.0 inches, sizes 3+4)

COLOUR: Dlack



THE MED-BENEFIT

- Very comfortable to wear, even for larger torsos
- Additional elastic strap to lift abdomen for relief of the lumbar spine
- Robust construction

BORT Vario Back Support with Pad REF 112 610



PRODUCT-BENEFIT

- Non-elastic back support for stabilisation of lumbar spine
- Non-elastic zone in the lumbar area made of breathable textile
- Lumbo-sacral friction pad
- Two supporting splints in the dorsal area
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)



 Improved fit due to fitted cut
 Additional elastic pull straps from dorsal to ventral area

INDICATIONS


Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

SIZES	Circumference of torso	SIZES	Circumference of waist
(Tarso)	Size 1 105 - 120 cm Size 2 120 - 135 cm Size 3 135 - 150 cm Size 4 150 - 165 cm	(Waist)	small up to 89 cm medium 89 – 104 cm large 104 – 119 cm x-large 119 – 134 cm



PRODUCT-BENEFIT

- Active support for muscular stabilisation of lumbar spine
- Two splints in the dorsal area
- Viscoelastic sacral pad, which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)



- Practical hand straps allow support to be fitted easily and without effort
- Comfortable to wear

BORT Vario Basic Back Support with Pad REF 112 670



PRODUCT-BENEFIT

- Elastic support for stabilisation of lumbar spine
- Sacral pad, which can be freely positioned
- Four anatomical pre-molded splints in the dorsal area
- Elastic cross straps
- Height: ventral 17 cm (6.7 inches), dorsal 26 cm (10.2 inches)



- Grip-straps allow easy fitting and removal of support
- Comfortable to wear

BORT Vario Basic Back Support without Pad REF 112 650

INDICATIONS

Treatment of pain in the lumbar spine, low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, spondylolysis without spondylolisthesis, mild lumbar deformity (degenerative spondylolisthesis)

SIZES	Circumference of waist	SIZES	Circumference of waist
(Waist)	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	(Waist)	x-small 63 - 74 cm small 75 - 89 cm medium 90 - 105 cm large 106 - 120 cm x-large 121 - 140 cm xx-large 141 - 160 cm

BORT LumboXpress Back Support with Double Fastening **REF** 112 600



PRODUCT-BENEFIT

- Elastic support for stabilisation of lumbar spine
- Four anatomical pre-molded splints in the dorsal area
- Two Velcro belts for safe stabilisation
- Two additional elastic pull straps from dorsal to ventral area enable individual dosage of the pressure
- Height: ventral 16 cm (6.3 inches), dorsal 26 cm (10.2 inches)

COLOUR: silver

THE MED-BENEFIT

- Specially constructed fastening system simplifies putting on
- Very flat construction, therefor hardly visible under clothing

INDICATIONS

Treatment of increased pain in the lumbar spine area, severe low back pain, recurrent low back pain, osteochondrosis, spondylarthrosis, mild radicular syndrome, spondylolysis without spondylolisthesis, mild to moderate lumbar deformity (degenerative spondylolisthesis), after intervertebral disc surgery

SIZES	Circumference of waist
(Waist)	Size 1 70 - 80 cm Size 2 80 - 90 cm Size 3 90 - 100 cm Size 4 100 - 110 cm Size 5 110 - 120 cm Size 6 120 - 130 cm Size 7 130 - 140 cm Size 8 140 - 150 cm
1 cm = 0.3937	' inch

BORT Back Support for Pregnant Woman **REF** 104 580



PRODUCT-BENEFIT

- Elastic back support for relief of lumbar spine and to support the abdomen
- Individual adjustment of width using flexible adjustable straps
- Thin spring rods for stabilisation of the lumbar region
- Dorsal height: 21 cm (8.3 inches)

COLOUR: white

THE MED-BENEFIT

- Adaptable to varying abdominal circumferences

INDICATIONS

For stabilisation of lumbar spine with metered pressure during pregnancy, for the treatment of pain in the lumbar area

SIZES

Circumference of hip

\ /	Size 0	75 – 91 cm
(Hip)	Size 1	91 - 105 cm
12QA	Size 2 1	105 – 119 cm
1197931	Size 3 1	119–140 cm



PRODUCT-BENEELT

- Non-elastic support for stabilisation of pelvis
- Pressure on the symphysis can be individually adjusted
- Anatomically designed cut
- Padded on the inside
- Height: ventral 8,5 cm (3.3 inches), dorsal 15 cm (5.9 inches)

COLOUR: silver

THE MED-BENEFIT

- Circular stabilisation and relief of pelvis and symphysis
- Comfortable to wear thanks to narrow closure section (for example when sitting)

INDICATIONS

Symphysis insufficiency, symphysis rupture, structural loosening in the area of the iliosacral joints/pubic symphysis. after birth

SIZ	ZES
-----	-----

Circumference of hip



small	85 – 95 cm
medium	95 – 105 cm
large	105 – 115 cm
x-large	115 – 125 cm
xx-large	125 – 135 cm

BORT Stabilo® Sacroiliac Brace **REF** 104 750



PRODUCT-BENEFIT

- Brace for the stabilisation of sacroiliac ioints
- Special combination of elastic and non-elastic knit fabric
- Two-part sacroiliac pad to stimulate the circulation and encourage the healing process
- Circular tension adjustable with additional tension straps

COLOUR: grey

THE MED-BENEFIT

- Circular stabilisation and relief of pelvis and symphysis
- Very comfortable to wear thanks to narrow closure section (for example when sitting)

INDICATIONS

Sacroiliac joint arthrosis, sacroiliac ioint instability, sacroiliac joint blockage, myalgias and tendopathies in the pelvic region, structural abnormalities after lumbar spinal spondylodeses

SI7FS Circumference of iliac crest



	-				-	
Sizes	1	_	7	0 – 8	35	cm
Sizes	2		85	i — 10)0	ст
Sizes	3		100	1 – 11	15	ст
Sizes	4		115	i – 13	30	ст
Sizes	5		130	1 – 14	15	ст



- Back brace for stabilisation of lumbar spine
- Integrated stabilisation elements
- Non-elastic, breathable carrier material
- Viscoelastic pads which can be freely positioned
- Height: ventral 17 cm (6.7 inches), dorsal 21 cm (8.3 inches)



THE GENERATION-BENEELT

COLOUR: blue

- Wide grip straps on the latch plate
- Innovative light brace construction
- Energy-saving pulley system
- Simple and easy to put on
- Easier to put on using the ergonomic latch plate

INDICATIONS

Degenerative lumbar change, mild to moderate lumboischialgia, slipped disc with mild deficits, spondylosis, spondylolisthesis, spinal canal stenoses, radicular irritation syndrome, condition after inter-vertebral disc surgery

SIZES	Circumference of waist
(West)	$\begin{array}{cccccccccccccccccccccccccccccccccccc$

1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Provides relief to the lumbar spine with lordosis reduction
- Non-elastic, breathable base material
- Back plate heights: 27 cm or 31 cm

COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Pelvic girdle, adjustable tapered
- Innovative, light brace construction
- Wide handle flaps for easy handling
- Energy-efficient pulley system
- Easier to put on using the ergonomic latch plate

INDICATIONS

Facet syndrome, lumbosciatica, invertebral disc protrusions/prolapse, radicular irritation syndrome, degenerative lumbar change, condition after inter-vertebral disc surgery, stable fractures, spondylosis, spondylolisthesis, spinal canal stenoses

SIZES	Circumference of waist
(Waist)	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$

Please indicate height 27 cm (10.6 inches) or 31 cm (12.2 inches) 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Lumbar brace for bridging and for relief of the lumbar spine
- Stimulation of abdominal muscles and straightening of pelvis
- Semi-elastic stabilisation strap, individually adjustable
- Height: ventral 17 cm (6.7 inches), dorsal 32 cm (12.6 inches)

THE MED-BENEFIT

- Anatomical cut for precise fit, flexible section in the hip-bending area
- The degree of stabilisation can be adjusted to suit the particular course of therapy
- Easy to put on using the hand straps provided

COLOUR: silver



INDICATIONS

Moderate lumboischialgia with mild deficits with herniated discs/prolapse, moderate radicular, pseudo-radicular lumbar syndrome, grade I spondylolisthesis with low back pain, moderate lumbar deformity (recurrent treatment, complicated history) i. e. after inter-vertebral disc surgery with mild to moderate prolapse, spinal canal stenoses, lumbar instability syndrome

SIZES

Circumference of waist

Size 1	70 – 80 cm 80 – 90 cm
Size 3	90 – 100 cm 100 – 110 cm
Size 5	110 – 120 cm
Size 6 Size 7	120 – 130 cm 130 – 140 cm



- Lumbar brace for bridging and for relief of the lumbar spine
- Classical brace-shaped stable plastic bridging frame, anatomically pre-moulded, thermoplastically adjustable
- Non-elastic stabilisation strap, individually adjustable
- Height: ventral 17 cm (6.7 inches), dorsal 32 cm (12.6 inches)

THE MED-BENEFIT

- High degree of stabilisation, can be gradually adjusted during the course of therapy (mobility function)
- Additional ventral abdominal reinforcing plate
- Basic support with an anatomic cut for optimum fitting
- Easy to put on using the hand straps provided

COLOUR: 🔳 grey

INDICATIONS

Severe lower back pain with muscular weakness and significant hernia protrusion/ prolapse, without absolute indication for surgery, severe lumbar radicular or pseudoradicular syndrome (refractory to conservative treatment), grade 2 spondylolisthesis with recurring lower back pain, severe lumbar deformity with facet syndrome/arthrosis, postsurgery for herniated disks without resulting paralysis with removal of the hernia, decompression of the spinal canal, single level, without deformity, single level herniated disk surgery, short stenosis of the spinal canal without significant paralysis, vertebral fractures resulting from minor trauma, superior end plate deformities in the lumbar spine

SIZES	Circumference of waist
Waist Size 2 Size 3 Size 4 Size 5 Size 6	70 - 80 cm 80 - 90 cm 90 - 100 cm 100 - 110 cm 110 - 120 cm 120 - 130 cm 130 - 140 cm



PRODUCT-BENEFIT

- Lumbar brace with integrated plastic pelvic girdle to secure the LV 1 LV 5 segments
- Continuously adjustable stabilisation of the lumbar spine with BOA® Fit System
- Integrated lumbar support with abdominal closing plate

THE MED-BENEFIT

- Rotation of the lumbar spine is effectively restricted
- Modular design: can be downgraded during the course of therapy
- Easily adjustable: Setting lumbar lordosis possible via cold deformation

COLOUR: grey

INDICATIONS

Major lumboischialgia without absolute indication for surgery, when surgery is contraindicated, when surgery is rejected, major radicular and pseudo radicular lumbar syndrome (cautionary, therapy-resistant), spondylolisthesis with severe lumboischialgia, severe lumbar deformity with hypermobility associated with Facet Syndrome / Osteoarthritis, condition after spine surgery, i.e. dorso-ventral/dorso-lateral fusion, spinal column compression (single level with deformity and multi-level), thoracic-surgery with severe prolapse/lumbar instability, spinal canal stenosis with paresis, post-traumatic vertebral fractures with significant spinal damage

SIZES	Circumference of waist
/Waati	70 – 80 cm 80 – 90 cm 90 – 100 cm 100 – 110 cm 110 – 120 cm 120 – 130 cm 130 – 140 cm



- Lumbar brace for partial relief of the lumbar spine and **facet joints**
- Two lateral supporting splints
- Height: ventral 15,5 cm (6.1 inches), dorsal 30 cm (11.8 inches)

THE MED-BENEFIT

- Elastic belt system to reduce lordosis of the lumbar spine
- Hand straps allow easy, effortless fitting
- Adjustable width using ergonomic latch plates

COLOUR: grey

INDICATIONS

Facet joint arthrosis, invertebral disc protrusions/prolapse, severe low back pain, lumboischialgia, lumbago, lumbago with sacroiliac joint involvement, chronic pain in the lumbosacral area, spondylolisthesis, muscular imbalance, segment instability, lumbar spinal canal stenoses

SIZES	Circumference of waist
Near A Size 1 Size 2 Size 3 Size 4 Size 5 Size 6 Size 7	70 – 80 cm 80 – 90 cm 90 – 100 cm 100 – 110 cm 110 – 120 cm 120 – 130 cm 130 – 140 cm



PRODUCT-BENEFIT

- Activating and spine-straightening back brace
- Stable aluminium back splint, anatomically pre-shaped and individual, without tools, adjustable
- Relief and stabilisation of the lumbar spine
- Active straightening by activating the back muscles
- Increasing mobility

THE GENERATION-BENEFIT

- Strap system according to the rucksack principle
- Guide loops in the chest area for easy, stepless straightening of the upper body
- Lightweight and slim structure for good wearing comfort
- Easy to handle
- Fastener with a large finger hold and hand straps for effortless fitting



BORT Generation:

For better mobility for the aged Product range with special handling advantages

COLOUR: blue

INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, damage in the thoracic and lumbar spine region, stable, osteoporotic vertebral fractures, osteoporosis

SIZES			Back length
torso	x-small . small . medium . large . x-large .		up to 43 cm 44 – 49 cm 50 – 55 cm 56 – 61 cm 62 – 67 cm
guilea		hs are individually adjustable for a girth 5 cm	

1 cm = 0.3937 inch

BORT Generation Dynamic Osteoporosis Brace REF 215 880



PRODUCT-BENEFIT

- Activating and spine-straightening back brace
- Back splint and limitation restraint, stable and dynamically adjustable
- Two-part, stable aluminium back splint, anatomically pre-shaped and individually adjustable without tools
- Relief and stabilisation of the lumbar spine
- Active straightening by activating the back muscles
- Increasing mobility
- COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Strap system according to the rucksack principle
- Guide loops in the chest area for easy, stepless straightening of the upper body
- Lightweight and slim structure for good wearing comfort
- Fastener with a large finger hold and hand straps for effortless fitting



BORT Generation: For better mobility for the aged Product range with special handling advantages

INDICATIONS

Active straightening for relief and correction of the lumbar/thoracic spine associated with restricted movement in the sagittal plane, damage in the thoracic and lumbar spine region, stable, osteoporotic vertebral fractures, osteoporosis

SIZES			Back length
torso gluteau	Belt lengt	hs are individually adjustable for a girth	up to 43 cm 44 – 49 cm 50 – 55 cm 56 – 61 cm 62 – 67 cm
	of 70 – 16	65 cm	

ää

1 cm = 0.3937 inch

Active and secure for everyday life: Straighten, compress, stabilise.



Overview of Indications for selected products	Rib Belt for Men 102 800	Scar Rupture Support 104 050	Umbilical Hernia Support 104 070	Stoma Support 104 090	Vario Stoma Support 104 010	PostOban® 104 150 PostOban®	PostOban [®] SOFT Thorax-Abdominal- Support	Abdominal Support for Pregnant Women 104 620	CoxaPro Hip Brace 201 400
	Rib Belt for Women 102 900					Special Width 104 180	104 100		
Page	P.123/122	P. 124	P.124	P.125	P. 125	P.126/127	P. 126	P. 127	P. 130
Rib fracture/contusion	•								
Incisional hernia		•							
Umbilical hernia			•						
Weakness of the connective tissues						•			
Paralysis of abdominal wall				•	•				
Condition after ileostomy and colostomy				•	•				
Relief of abdomen								•	
Pendulous abdomen							•		
Post-operative		•				•			
Tractus iliotibialis syndrome									•
Coxarthrosis									•
Piriformis syndrome									•
Condition after hip joint luxation									•

_____ OVERVIEW OF INDICATIONS 121

Torso Groin Hip



PRODUCT-BENEFIT

- Non-elastic circular support for stabilisation of the thorax
- Non-elastic side and front section with Velcro fastening, elastic back section
- Breathable fabric with Velcro fastening
- Feminine cut fitted at the waist
- Height: dorsal 16 cm (6.3 inches), ventral 12 cm (4.7 inches)

COLOUR: White

THE MED-BENEFIT

- Hardly noticeable, invisible under clothing
- Very comfortable to wear due to the use of cotton satin and fine terry cloth material



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of the thorax
- Elastic, breathable fabric with Velcro fastening
- Feminine cut fitted at the waist
- Height: dorsal 16 cm (6.3 inches), ventral 12 cm (4.7 inches)

COLOUR: White

THE MED-BENEFIT

- Hardly noticeable, not visible under clothing
- Good wearing comfort

INDICATIONS

Rib fracture, trauma consequences in the rib area (e. g. contusion, sprain)

SIZES

Underbust circumference

SIZES

Underbust circumference

x-large ____ 106 - 125 cm

65 – 78 cm

____ 78 – 91 cm ____ 91 – 106 cm

Small up to 66 cm medium 66 - 76 cm large 76 - 89 cm x-large 89 - 102 cm xx-large 102 - 125 cm
--

1 cm = 0.3937 inch

1 cm = 0.3937 inch

small _ medium _

large



- Non-elastic circular support for stabilisation of the thorax
- Non-elastic side and front section, elastic back section
- Breathable fabric with Velcro fastening
- Straight cut
- Height: 16 cm (6.3 inches)



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of the thorax
- Elastic, breathable fabric with Velcro fastening
- Straight cut
- Height: 16 cm (6.3 inches)

$COLOUR: \Box$ white

THE MED-BENEFIT

- Hardly noticeable, not visible under clothing
- Very comfortable to wear due to the use of cotton satin and fine terry cloth material

COLOUR: White

THE MED-BENEFIT

- Hardly noticeable, not visible under clothing
- Good wearing comfort

INDICATIONS

Rib fracture, trauma consequences in the rib area (e. g. contusion, sprain)

SIZES	Circumference of chest	SIZES	Circumference of chest
je t	small up to 82 cm medium 82 – 96 cm large 96 – 112 cm x-large 112 – 128 cm xx-large 128 – 144 cm		small 65 – 78 cm medium 78 – 91 cm large 91 – 106 cm x-large 106 – 125 cm xx-large 125 – 145 cm
1 cm = 0.39	37 inch	1 cm = 0.39	37 inch



PRODUCT-BENEFIT

- Elastic compression support for incisional hernia
- Round silicone pad reinforced with plastic, for individual positioning
- Two supporting splints
- Elastic, breathable fabric with Velcro fastening
- Heights: 13 cm (5.1 inches), 21 cm (8.3 inches), 26 cm (10.2 inches)

$COLOUR: \square$ white

Schematic representation of pad

THE MED-BENEFIT


- Removable pad cover made of soft terrycloth
- Height 13 cm (5.1 inches): round pad with 10 cm (3.94 inches) diameter
- Height 21 cm (8.3 inches): pad length × width: 18 (7.09) × 10 cm (3.94 inches), oval
- Height 26 cm (10.2 inches): pad length × width: 22 (8.7) × 10 cm (3.94 inches), oval

INDICATIONS

Scar rupture, scar insufficiency, abdominal muscle infirmity/paralysis

SIZES

Circumference of waist



please indicate height 13 cm (5.1 inches), 21 cm (8.3 inches) or 26 cm (10.2 inches) 1 cm = 0.3937 inch

BORT Umbilical Hernia Support REF 104 070



PRODUCT-BENEFIT

- Elastic compression support for umbilical hernia
- Anatomically concaved umbilical hernia pad, with adhesive layer for individual positioning
- Two supporting splints
- Elastic, breathable fabric with Velcro fastening
- Height: 13 cm (5.1 inches)

COLOUR: White



Schematic representation of pad

THE MED-BENEFIT

- Removable pad cover made of soft terry cloth
- Round pad with 10 cm (3.94 inches) diameter

INDICATIONS

Umbilical hernia, reducible umbilical hernias

SIZES	Circumference	of waist

Size 1	 60 – 95 cm
Size 2	 95 – 125 cm
Size 3	 125 – 150 cm

1 cm = 0.3937 inch



- Circular abdominal support with a variable stoma recess
- Non-elastic composite velour material, elastic insert in the back section
- Two supporting splints
- Plastic reinforcing plate to be inserted
- Heights: 15 cm (5.9 inches), 22 cm (8.7 inches)

COLOUR: White

BORT Vario Stoma Support REF 104 010



PRODUCT-BENEFIT

- Circular abdominal support with a variable stoma recess
- Non-elastic composite velour material, elastic insert in the back section
- Two supporting splints
- The size, shape and positioning of the stoma recess can be individually determined
- Heights: 15 cm (5.9 inches), 22 cm (8.7 inches)

COLOUR: White

THE MED-BENEFIT

recess

area

Washable at 90°C

- Stencil provided for easy, individual

size determination of the stoma

- Large, skin-friendly lining in stoma


THE MED-BENEFIT

- The variable stoma recess width can be adjusted using a Velcro fastening
- Large, skin-friendly lining in stoma area
- Easy to put on and position the stoma recess using the Velcro fastener on the right and left in the abdominal region
- Washable at 90°C

ACCESSORIES: REF 104 099 | BORT Stoma Protector

INDICATIONS

Paralysis of abdominal wall, hernia of abdominal wall, condition after ileostomy and colostomy with accompanying weakness of connective tissue





PRODUCT-BENEFIT

- Elastic circular support for stabilisation of thoracic and abdominal regions
- Elastic, soft fabric
- Wide, two-piece Velcro fastening
- Heights: 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: beige

THE MED-BENEFIT

- Individually adjustable
- Comfortable to wear

BORT PostOban® Thorax-Abdominal-Support REF 104 150



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of thoracic and abdominal regions
- Elastic, breathable fabric
- Wide, two-part Velcro fastening
- Heights: 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: White

- THE MED-BENEFIT
- Individually adjustable
- Two-part Velcro fastening makes the support easy to put on and ensures a better fit

INDICATIONS

After abdominal surgery, weakness of the connective tissue in the abdominal area, pendulous abdomen

INDICATIONS

After abdominal/thoracic surgery, weakness of the connective tissue in the abdominal area, pendulous abdomen

SIZES	Circumference of waist

1 /	small	 up to 87 cm
(Waist)	medium	 87 – 102 cm
	large	 102 – 117 cm
KMAY	x-large	 $117-132\ \text{cm}$

please indicate height 21 cm (8.3 inches) or 26 cm (10.2 inches) 1 cm = 0.3937 inch SIZES

Circumference of waist

\ /	Size 1	 75 – 85 cm
(Waist)	Size 2	 85 – 95 cm
2. La	Size 3	 95 – 105 cm
119731	Size 4	 105 – 115 cm
	Size 5	 115 – 125 cm
	Size 6	 125 – 135 cm
	Size 7	 135 – 145 cm

please indicate height 21 cm (8.3 inches) or 26 cm (10.2 inches) 1 cm = 0.3937 inch

BORT PostOban® Special Width Thorax-Abdominal-Support [REF] 104 180



PRODUCT-BENEFIT

- Elastic circular support for stabilisation of thoracic and abdominal regions for larger sizes
- Two-piece elastic fastener
- Elastic, breathable fabric
- Heights: 21 cm (8.3 inches), 26 cm (10.2 inches)

COLOUR: White

THE MED-BENEFIT

- Very comfortable to wear, even for larger torsos
- Individually adjustable
- Integrated dressing aids
- Easy to put on

BORT Abdominal Support for Pregnant Women REF 104 620



PRODUCT-BENEFIT

- Elastic bandage to support the abdomen
- Elastic, breathable fabric
- Four supporting rods which can be removed, if required
- Dorsal height: 24 cm (9.5 inches)

COLOUR: White

THE MED-BENEFIT

- Back relief due to additional elastic strap
- Adjustable to fit varying abdominal circumferences

INDICATIONS

After abdominal surgery, weakness of the connective tissue in the abdominal area, pendulous abdomen

INDICATIONS

Lumbar back pain in pregnancy, obesity, multiple pregnancy

SIZES	Circumference of torso
(Terce)	Size 0 90 - 105 cm Size 1 105 - 120 cm Size 2 120 - 135 cm Size 3 135 - 150 cm Size 4 150 - 165 cm Size 5 165 - 180 cm

please indicate height 21 cm (8.3 inches) or 26 cm (10.2 inches) 1 cm = 0.3937 inch 1 cm = 0.3937 inch

SIZES Circumference of iliac crest





PRODUCT-BENEFIT

- Single-sided, rounded pressure pad
- Non-elastic hip belt with approx.
 50 cm (19.7 inches) long spring made of spring steel

COLOUR: beige

BORT HerniaFix® Elastic Spring Truss with Rounded Pad, Short Spring REF 109 400

PRODUCT-BENEFIT

- Single-sided, rounded pressure pad
- Elastic hip belt with short spring 20 cm (7.9 inches) made of spring steel

with Anatomically-Shaped Pad REF 109 200

BORT HerniaFix® Spring Truss

PRODUCT-BENEFIT

- Single-sided, anatomically-shaped pad
- Non-elastic hip belt with approx.
 50 cm (19.7 inches) long spring made of spring steel
- Thigh strap prevents pad from slipping

COLOUR: beige

BORT HerniaFix* Elastic Spring Truss with Anatomically-Shaped Pad, Short Spring **REF** 109 700

PRODUCT-BENEFIT


- Single-sided, anatomically-shaped pad
- Elastic hip belt with short spring 20 cm (7.9 inches) made of spring steel
- Thigh strap prevents pad from slipping

COLOUR: beige

COLOUR: beige

INDICATIONS

Reducible inguinal rupture



please indicate left or right 1 cm = 0.3937 inch SIZES Circumference of iliac crest

\ /	85	83 – 88 cm
(lliac+	90	88 – 93 cm
West ALS	95	93 – 98 cm
KIGAN	100	98 – 103 cm
111 11 111	105	103 – 108 cm

please indicate left or right 1 cm = 0.3937 inch

Groin



PRODUCT-BENEFIT

- Double-sided, rounded pressure pad
- Elastic hip belt with approx. 20 cm (7.9 inches) long spring made of spring steel

BORT HerniaFix® Special Spring Truss with Adjustable Anatomically-Shaped Pad [REF] 110 300



PRODUCT-BENEFIT


- Special spring truss, no pressure or spring, single sided
- Anatomically adjustable pad
- Elastic hip belt
- Thigh strap

COLOUR: beige

COLOUR: beige

INDICATIONS

Reducible inguinal rupture





PRODUCT-BENEFIT

- Hip joint brace for securing the function and/or mobilisation of the hip joint
- Simple restriction of movement for bending and stretching
- Tool-free extension and flexion setting 0°, 15°, 30°, 45°, 60°, 75°, 90°
- Trochanter pad for intermittent massage
- Joint splint in the hip region and thigh can be adapted using cold-shaping (aluminium splint)

COLOUR: grey

THE MED-BENEFIT

- Good functional security for the hip joint via two wide, opposing restraints on the thigh
- Light, slim design for very comfortable wearing
- High-quality functional cotton knitted fabric, particularly Velcro-resistant



ACCESSORIES: REF 991 126 Piriformis pad – For trigger point massage of the piriformis muscle

INDICATIONS

All indications for which functional security and/or mobilisation of the hip joint is necessary, such as state after hip joint luxation, loosening hip endoprostheses, medium to severe coxarthrosis, muscular imbalance of the lumbar pelvic hip region, e.g. tractus iliotibialis syndrome, piriformis syndrome

SIZES	Circumference of iliac crest	f-measure
(lliac)	small 85 - 100 cm medium 100 - 115 cm large 115 - 130 cm x-large 130 - 145 cm	. 39 – 52 cm . 45 – 65 cm . 54 – 66 cm . 60 – 76 cm

please indicate left or right 1 cm = 0.3937 inch

BORT Generation Hip Protector REF 215 900



PRODUCT-BENEFIT

- Hip protection beltwith a shock absorbing effect
- Protectors made of soft visco-elastic foam can be removed
- Strap with an anatomically designed cut can be fully opened

COLOUR: 🗖 blue

THE GENERATION-BENEFIT

- Easy to put on/remove thanks to Velcro fastening
- Wide frontal grip straps
- Safe handling and comfortable to wear
 Ideal for temporary use during phases of activity
- Attractive, functional design

ACCESSORIES:

- Hip Protector set (REF) 201 120)

INDICATIONS

Tendency to fall for various reasons, danger of falling



PRODUCT-BENEFIT

- Hip protection pants with shock absorbing effect
- Protectors made of soft visco-elastic foam can be removed

$COLOUR: \square$ white

THE MED-BENEFIT

- Highly elastic waistband and leg edge
- Protectors can be easily removed
- High-quality functional pants made of cotton and Lycra

ACCESSORIES:

- Hip Protector set (REF 201 120)
- Pants (REF 201 200)

INDICATIONS

Tendency to fall for various reasons

SIZES	Circumference of iliac crest
(lliac crest)	small 80 - 90 cm medium 90 - 100 cm large 100 - 110 cm x-large 110 - 120 cm xx-large 120 - 130 cm

	Dress	Women	Men
SIZES	Size	GB	GB
small		38	46
medium		_ 40 _	48
large		_ 42	50
x-large		44	52
xx-large		46/48	54/56
xxx-large		_ 50/52	58/60

Walking into the future: functional stabilisation with maximum mobility.



300 e a ° e -	8												
Stab Prr Stab Prr Stab 182 Sty Sty	P. 148							•	•				
OA- Genu- Stabilo- Xpress Xpress Pro® 101 600 100 700 182 150 Genu- Xtabilo- Xtabilo- Xpress Pro® 100 800 Xpress Pro® 182 300	P. 152							•					•
0A- Xpress 101 600	P. 151									•			
Genera- tion Soft OA 1 Knee Brace 215 440	P. 150									•			
Genera- tion Knee Brace 215 400	P: 149						•	•	•				
Immob Splint 145 000 Immo- Xpress 145 300	P. 154/155										•		
Stabilo® Knee Support with Hinge 150 140	P. 147							•					
Osgood- Schlatter Knee Support 114 510	P. 141	•										•	
Centeration Asymmetrice Osgood-Stabilo® Knee 114 700 Schlatter Knee Nupport 215 420 Plus Nupport 215 420 Plus 114 510 Hinge Reneration 114 900 114 510 Hinge Knee Asymmetrice Support Vith Plus 114 900 114 510 Hinge Knee Asymmetrice Support Vith Plus 6rip Top 114 400 150 140	P. 146/145		•	•									
Generation Knee Support 215 420 Generation Knee Support Plus 114 490 G	P. 137	•					•		•				
Stabilo Gen [®] Eco 114 520	P. 138	•					•	•	•				
Patella Tendon Strap Sport 114 010 SP	P. 143				•	•							
StabiloGen [®] Sport 114 520 SP StabiloGen [®] Sport with a Grip Top 114 520 SP-H	P. 136	•					•	•	•				
select* Stabilo- Gen® 114 420 select* GenuZip® 114 390	P. 134/135	•					•	•	•				
Overview of indications for selected products	Page	Irritations	Patellar lateralisation	Patellar (sub) luxation	Patellar chondropathy	Patellofemoral arthrosis	Meniscopathy	Ligament instability	Gonarthrosis	Varus/Valgus Gonarthrosis	Complete immobilisation	Osgood-Schlatter disease	ACL/PCL ruptur

_____ OVERVIEW OF INDICATIONS 133



PRODUCT-BENEFIT

- High-quality knee support for compression of soft tissue with a viscoelastic pad, for stabilisation and relief of the knee joint
- Lengthwise-elastic stretch zone above the viscoelastic pad for an ideal fit at all bending angles
- With an highliy elastic pad for the popliteal fossa, providing the best possible comfort for wearing
- One spiral stay on each side for lateral support
- Compression class II

COLOURS: beige silver

THE SELECT-BENEFIT

- SOFTflex Plus knit fabric:
 - \cdot SENSIL® Breeze: Luxury fibre with a cooling effect
 - Soft-velvety surface structure for top wearing comfort
- Extra-soft, snag proof material edge
- Anatomically and functionally optimised pad form for safe guidance of the patella
- Plus-size system results in a good fit

BORT INNOVATION: - Woven silicone bonding zone in the thigh area

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

SIZES	d-measure	f-measure	d-measure	f-measure
Jis	Size 1 26 - 29 Size 2 29 - 33 Size 2+ 29 - 33 Size 3 33 - 37 Size 3+ 33 - 37 Size 4 37 - 41	44 – 48 cm		

right and left symmetrical 1 cm = 0.3937 inch



- High-quality knee support for compression of soft tissue with a viscoelastic pad, for stabilisation and relief of the knee joint
- Side zip fastener
- Lengthwise-elastic stretch zone above the viscoelastic pad ensures an ideal fit at all bending angles
- With an highliy elastic pad for the popliteal fossa, providing the best possible comfort for wearing
- One spiral stay on each side for lateral support
- Compression class II

THE SELECT-BENEFIT

- SOFTflex Plus knit fabric:
 - · SENSIL® Breeze: Luxury fibre with a cooling effect
 - Soft-velvety surface structure for top wearing comfort
- Side zip fastener makes it easy to put on/remove
- Extra-soft, snag proof material edge
- Very comfortable to wear
- Anatomically and functionally optimised pad form for safe guidance of the patella
- Plus-size system results in a very good fit

COLOUR: silver

BORT INNOVATION:

Woven silicone bonding zone in the thigh area

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrose, arthritis, e. g. ligament instability, meniscopathy

SIZES	d-measu	ure f-measur	е	d-measure	f-measure
	Size 1 26 Size 2 29 Size 2+ 29 Size 3 33 Size 3+ 33 Size 4 37	37 40 – 44 ci 37 44 – 48 ci	n Size 6 n Size 6+	45 - 49	60 – 64 cm

please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Knee support for compression of soft tissue with a viscoelastic pad for stabilisation and relief
- SOFTflex 3D knitted fabric with a velvety soft surface structure
- One lateral spiral stay on each side
- Plus-size system results in a good fit
- Highly elastic, soft flexion region in the hollow of the knee

COLOUR: black/green

THE SPORT-BENEFIT

- Woven silicone bonding zone in the thigh area
- COOLMAX[®]: wicks moisture away from the body
- Latex-free knitted fabric
- Anatomically and functionally optimised pad form for safe guidance of the patella during intense activity

INDICATIONS



PRODUCT-BENEELT

- Knee support for compression of soft tissue with a viscoelastic pad for the stabilisation and relief
- With additional silicone grip top in all sizes
- SOFTflex 3D knitted fabric with a velvety soft surface structure
- One lateral spiral stay on each side
- Double plus size system resulting in very good fitting

COLOUR: black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body, for dryness and comfort
- Latex-free fabric
- Anatomically and functionally optimised pad form for safe guidance of the patella during intense activity
- Support length in size 3: 35 cm

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. gonarthrose, arthritis, ligament instability, meniscopathy

SIZES	d-measure	f-measure		d-measure	f-measure
Size 2 Size 2 Size 2 Size 3 Size 3 Size 3 Size 4 Size 4 Size 4	29 - 33 29 - 33 33 - 37 33 - 37 33 - 37 37 - 41	$\begin{array}{c} 36-40\ \text{cm} \\ 40-44\ \text{cm} \\ 40-44\ \text{cm} \\ 44-48\ \text{cm} \\ 44-48\ \text{cm} \\ 48-52\ \text{cm} \end{array}$	Size 5+ *Size 5++ Size 6 *Size 6+ *Size 6++ *Size 7 *Size 7+	41 - 45 45 - 49 45 - 49 45 - 49 45 - 49 49 - 53	$\begin{array}{c} 60-64\ cm\\ 60-64\ cm\\ 64-68\ cm\\ 68-72\ cm\\ 68-72\ cm\\ \end{array}$

right and left symmetrical 1 cm = 0.3937 inch



- Knee support with adjustable thigh width for stabilisation of the knee joint
- One lateral spiral stay each side
- Functional knitted fabric, very comfortable to wear
- Anatomical, viscoelastic pad
- Integrated plus size system

COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Easy to put on/remove thanks to closed design on the shank and open design on the thigh
- Grip straps allow safe and energysaving handling, also in case of restricted gripping ability
- Simple and easy to put on
- Attractive, functional design

BORT Generation Knee Support Plus REF 114 490 G



PRODUCT-BENEFIT

- Knee support with adjustable thigh width for stabilisation of the knee joint
- Circular closure strap for improved stabilisation on the thigh
- One lateral spiral stay each side
- Functional knitted fabric, very comfortable to wear
- Anatomical, viscoelastic pad
- Integrated plus size system

COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Easy to put on/remove thanks to closed design on the shank and open design on the thigh
- Grip straps allow safe and energysaving handling, also in case of restricted gripping ability
- Simple and easy to put on
- Attractive, functional design

INDICATIONS

Gonarthrosis, recurrent articular effusion, arthritis, slight ligament instability

SIZES	d-measure	f-measure
Size 2 Size 3 Size 4 Size 5 Size 6	29 – 33 cm 33 – 37 cm 37 – 41 cm	36 - 40 cm 40 - 48 cm 44 - 52 cm 52 - 60 cm 60 - 68 cm



PRODUCT-BENEFIT

- Classical knee support for compression of soft tissue with a viscoelastic pad for stabilisation and relief of the knee joint
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- One lateral spiral stay each side
- Plus-size system resulting in very good fitting
- Compression class II
- For children too, see page 32

COLOURS: beige silver black denim

THE ECO-BENEFIT

- Breathable functional comfort zone in the popliteal fossa
- Velvety soft surface structure
- Optimal comfort
- Siliconeadhesive border from 5 size
- Length of the support in size 3: 32 cm

BORT StabiloGen® Latex-Free REF 114 250

PRODUCT-BENEFIT

- Latex-free knee support for compression of soft tissues
- Stabilisation and relief of the knee joint
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Highly elastic, soft flexion region in the hollow of the knee
- Viscoelastic pad
- One spiral stay on each side
- Compression class II

COLOUR: black

THE MED-BENEFIT

- Latex-free knitted fabric
- Velvety soft surface structure
- Plus-size system resulting in very good fitting
- Silicone adhesive border from 5 size

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e.g. gonarthrose, arthritis, ligament instability, meniscopathy

SIZES		d-measure	f-measure	d-measure f-measure
is m	Size 1	23 - 26 26 - 29 29 - 33 29 - 33 33 - 37 37 - 41 37 - 41	32 - 36 cm 36 - 40 cm 40 - 44 cm 40 - 44 cm 44 - 48 cm 44 - 48 cm	$\begin{array}{cccccccccccccccccccccccccccccccccccc$

right and left symmetrical 1 cm = 0.3937 inch *Children size only available for REF 114 520



- Classical knee support for compression of soft tissue with a viscoelastic pad for stabilisation and relief of the knee joint
- With additional silicone grip top in all sizes
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- One lateral spiral stay each side
- Double plus size system resulting in very good fitting

COLOUR: silver

THE ECO-BENEFIT

- Breathable functional comfort zone in the popliteal fossa
- Optimal adhesion and pressure distribution on the edge of the thigh
- Velvety soft surface structure
- Optimal comfort
- Length of the support in size 3: 37 cm

BORT StabiloGen® Eco with a Grip Top Short IREF 114 520 HK

PRODUCT-BENEFIT

- Classical knee support for compression of soft tissue with a viscoelastic pad for stabilisation and relief of the knee joint
- With additional silicone grip top in all sizes
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- One lateral spiral stay each side
- Double plus size system resulting in very good fitting

COLOURS: silver black

THE ECO-BENEFIT

- Breathable functional comfort zone in the popliteal fossa
- Short version in the thigh area, no extension thanks to the silicone grip top
- Optimal comfort
- Length of the support in size 3: 32 cm

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, gonarthrosis, arthritis, ligament instability, meniscopathy

SIZES	d-measure	f-measure	d-measure	f-measure
	Size 1 $26 - 29$ Size 2 $29 - 33$ Size 3 $33 - 37$ Size 3+ $33 - 37$ Size 4 $37 - 41$ Size 5 $41 - 45$	36 - 40 cm 40 - 44 cm 40 - 44 cm 44 - 48 cm 44 - 48 cm 48 - 52 cm	Size 6 45 - 49 Size 6+ 45 - 49 Size 6++ 45 - 49 Size 7 49 - 53	60 - 64 cm 60 - 64 cm 64 - 68 cm 68 - 72 cm 68 - 72 cm



PRODUCT-BENEFIT

- Knee support for compression of soft tissue with a viscoelastic pad for stabilisation and relief of the knee joint
- Anatomically contoured knit
- One lateral spiral stay each side
- Compression class II

COLOURS: beige mineral grey

THE MED-BENEFIT

- Plus-size system resulting in very good fitting
- Soft surface structure



PRODUCT-BENEFIT

- Knee support for compression of soft tissues with V-shaped stabilisation bands for stabilisation and relief of the knee joint
- Viscoelastic pad
- Silicone non-slip border in the proximal and distal area
- Long femoral part (22 cm, 8.7 inches)
- Length: 37 cm (14.6 inches)

COLOUR: silver

THE MED-BENEFIT

- No lateral splints, there is no pressure in this area
- Relax zone in the knee bend
- Plus-size system resulting in very good fitting

INDICATIONS

Chronic, post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, e. g. arthritis, gonarthrose, ligament instability

INDICATIONS

Chronic post-traumatic or post-operative irritations of knee joint soft tissues, recurrent articular effusion, arthritis, arthrosis


SI	ZE	S
----	----	---

d-measure f-measure

d-measure f-measure

	Size 3+		33 – 37		Size 4+ Size 5 Size 5+ Size 6 Size 6+		
--	---------	--	---------	--	---	--	--

right and left symmetrical 1 cm = 0.3937 inch



- Knee support with patella recess for stabilisation and relief of the knee joint
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Viscoelastic pad
- One lateral spiral stay each side
- Silicone adhesive border from 5 size

COLOUR: silver

THE MED-BENEFIT

- Breathable functional comfort zone in the popliteal fossa
- Velvety soft surface structure
- Complete pressure relief of the patella
- Plus-size system resulting in very good fitting

INDICATIONS

Patellar chondropathy, meniscopathy, arthritis, bursitis

BORT Osgood-Schlatter Knee Support International Internationea International International Internati

PRODUCT-BENEFIT

- Knee support for local decompression of tibial tuberosity
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Anatomically contoured knit
- Multi-zone silicone pad
- One lateral spiral stay on each side
- For children too, see page 31

COLOUR: silver

THE MED-BENEFIT

- Velvety soft surface structure
- Plus-size system resulting in very good fitting
- Adjustable setting using V-shaped strap correction

INDICATIONS

Facilitation of the movement of the patella, Osgood-Schlatter disease, tendinosis of the tibial tuberosity ligament, gonarthrosis, arthritis, articular effusion, irritations

SIZES		d-measure	f-measure	d-measure	f-measure
Astronom d	*Children (blue) Size 1 Size 2 Size 2+ Size 3 Size 3+	26 - 29 29 - 33 29 - 33 33 - 37	32 – 36 cm 36 – 40 cm 40 – 44 cm	$\begin{array}{rrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrr$	48 – 52 cm 52 – 56 cm 56 – 60 cm 60 – 64 cm

right and left symmetrical 1 cm = 0.3937 inch *Children size only available for REF 114 510



PRODUCT-BENEFIT

- Active lower leg support with circular straps and cross restraints
- Anatomically contoured knit
- Relief zones in the injury region
- Length: 24 cm (9.4 inches) (small)
 25 cm (9.8 inches) (medium)
 26 cm (10.2 inches) (large; x-large)

COLOUR: I black/green

THE SPORT-BENEFIT

- Its operating principle consists of compression, conventional tape and Kinesio tape
- Proprioceptive effect
- Modern, sporty design
- Latex-free knitted fabric

BORT MyoActiv Sport Thigh Support REF 116 600 SP



PRODUCT-BENEFIT

- Active thigh support with circular straps and cross restraints
- Anterior/posterior application is possible
- Anatomically contoured knit
- Relief zones in the injury region
- Length: 23 cm (9.1 inches)

COLOUR: I black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body
- Its operating principle consists of compression, conventional tape and Kinesio tape
- Proprioceptive effect
- Modern, sporty design
- Latex-free knitted fabric

INDICATIONS

For soft tissue compression on the lower leg in case of acute muscle fiber tears and/or acute muscle strains

INDICATIONS

Preventive: For highly strained muscle groups, Acute: Swollen contusions with haematoma, muscle injury, Chronic: following muscle injury, protection against new muscle lesions, reduction of the burden on damaged muscle structures (e.g. scarring after old injuries)

SIZES

Circumference of calf

	small medium large x-large	up to 33 cm 33 – 38 cm 38 – 43 cm 43 – 47 cm
لنسكه	x-large	43 – 47 cm

right and left symmetrical 1 cm = 0.3937 inch SIZES

Circumference of thigh

small up to 43 medium 43 – 50 large 50 – 56 x-large 56 – 61 xx-large 61 – 66	cm cm cm

right and left symmetrical 1 cm = 0.3937 inch

THIGH/CALF/KNEE 143



PRODUCT-BENEFIT

- So-called Kassel Patellar Tendon Strap for target relief of the patella tendon
- Sports version without a guide loop
- Optimised construction for high activity levels
- Soft, elastic strap material

COLOUR: I black/green

THE SPORT-BENEFIT

- Also suitable for contact sport types
- Comfortable to wear during physical activity
- Individually controllable pad pressure
- Anti-slip positioning of the silicone pad via direct skin contact
- Modern, sporty design

BORT Patella Tendon Strap (REF 114 010

PRODUCT-BENEFIT

- So-called Kassel Patellar Tendon Strap for target relief of the patella tendon
- Short stretch elastic material
- Anatomical silicone proprioceptive pad
- Soft bordering
- Reversible strap with Velcro fastening

COLOUR: grey

THE MED-BENEFIT

- Pad pressure adjustable without graduation
- Improved joint coordination and pain relief
- No pressure in the hollow of the knee

INDICATIONS

Patellar chondropathy, femoropatellar arthrosis, patellar tendonitis

SIZES	Ci	rcumference of knee joint below patella
-	Size 1 Size 2 Size 3	up to 34 cm 34 – 41 cm 41 – 48 cm

right and left symmetrical 1 cm = 0.3937 inch
 SIZES
 below patella

 JUL/ Size 1
 up to 34 cm

 Size 2
 34 – 41 cm

 Size 3
 41 – 48 cm

Circumference of knee joint

right and left symmetrical 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Knee support with side joint splints, which can be fully opened
- Four Velcro fastenings
- Ideal in case of restrictions of movement and injuries that do no allow it to be slipped on
- Non-slip fit

COLOUR: silver

THE MED-BENEFIT

- The support can be opened completely thanks to its Velcro fastenings
- Support can be put on/removed easily (ideal for use in the elderly)

BORT Stabilo® Knee Support Special Width REF 114 490



PRODUCT-BENEFIT

- Knee support with adjustable thighwidth for guidance and stabilisation of the knee joint
- With two lateral spiral stays on each side and lateral silicone pads
- Patella recess
- Two Velcro fastening straps for optional adjustment of upper thigh width

COLOUR: beige

THE MED-BENEFIT

- Particularly suitable for patients with broader thighs
- Support can be put on/removed easily
- Non-slip fit

INDICATIONS

Mild to moderate lateral ligament instability of the knee joint, arthrosis, arthritis

INDICATIONS

Recurrent articular effusion, arthrosis, arthritis, ligament instability

Circumference

SIZES	Circumference of knee	d-measure	f-measure
	Size 1 _ up to 34 Size 2 _ 34-37 Size 3 _ 37-40 Size 4 _ 40-43 Size 5 _ 43-46 Size 6 _ 46-49	30-35 32-38 36-43 40-46 44-50	36-40 cm 40-44 cm 44-48 cm 48-52 cm 52-56 cm 56-60 cm
	Size 7 _ 49-52	48-54	60-64 cm

right and left symmetrical 1 cm = 0.3937 inch

SIZES	of knee	d-measure f-measure		
 	Size 1 _ 34-37 Size 2 _ 37-40 Size 3 _ 40-43 Size 4 _ 43-46 Size 5 _ 46-49 Size 6 _ 49-52 Size 7 _ 52-55 Size 8 _ 55-58	29–33 36–45 cm 33–37 40–50 cm 37–41 45–55 cm 41–45 50–60 cm 45–49 55–65 cm 49–53 58–66 cm 53–57 63–78 cm 57–61 67–82 cm		
vielet and laft assure thing!				

right and left symmetrical 1 cm = 0.3937 inch

BORT Asymmetric[®] Sport with Grip Top REF 114 700 SP H



PRODUCT-BENEFIT

- Knee support for influencing the patella glide path
- SOFTflex 3D knitted fabric for more elasticity, stability and very comfortable to wear
- With additional silicone grip top in all sizes
- Double Plus size system, therefore an excellent fit
- One lateral spiral spring stay on each side
- Compression class II

COLOUR: I black/green

THE SPORT-BENEFIT

- Pad pressure individually and directly adjustable using two Velcro straps
- COOLMAX[®]: Wicks moisture away from the body
- Ideal fit for sports activities
- Latex-free fabric
- Modern sporty design (with a reflecting label)
- Support length in size 3: 35 cm

INDICATIONS

Patellofemoral pain syndrome, high riding patella, patellar lateralisation, condition after lateral release surgery, recurrent or chronic patellar instability, meniscopathy

SIZES	d-measure	f-measure	d-measure	f-measure
lisem Hisem	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	$\begin{array}{c} 36-40 \text{ cm} \\ 40-44 \text{ cm} \\ 40-44 \text{ cm} \\ 44-48 \text{ cm} \\ 44-48 \text{ cm} \\ 48-52 \text{ cm} \end{array}$	$\begin{array}{rrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrr$	60 - 64 cm 60 - 64 cm 64 - 68 cm 68 - 72 cm 68 - 72 cm



PRODUCT-BENEFIT

- Knee support to facilitate movement of the patella
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- One lateral spiral stay on each side
- Compression class II

BORT Asymmetric® Plus REF 114 900



PRODUCT-BENEFIT

- Soft brace made of textile with joints to facilitate movement of the patella
- Lateral silicone half ring
- Side joint splints
- Circular fastening straps on the thigh and lower leg
- Compression class II

COLOUR: silver

THE MED-BENEFIT

- Pad pressure individually and directly adjustable using two Velcro straps
- Woven silicone bonding zone in the thigh area
- Plus size system, resulting in very good fitting

INDICATIONS

Patellofemoral pain syndrome, high riding patella, patellar lateralisation, condition after lateral release surgery, recurrent or chronic patellar instability, meniscopathy

COLOUR: silver

THE MED-BENEFIT

- Vector correction filaments: Innovative strap system which enables individual gradual guidance of the patella
- Individually adjustable joint splints
- Plus size system, resulting in very good fitting

INDICATIONS

High riding patella, patellar lateralisation, patellofemoral pain syndrome, patellar luxation or subluxation combined with lateral ligament instability, post-operative

SIZES		d-measure	f-measure		d-measure	f-measure
135 - 4 15 cm	Size 2 Size 2+ Size 3 Size 3+	26 - 29 29 - 33 29 - 33 33 - 37 33 - 37 37 - 41	36 – 40 cm 40 – 44 cm 40 – 44 cm 44 – 48 cm	Size 5 Size 5+ Size 6	37 - 41 41 - 45 41 - 45 41 - 45 45 - 49 45 - 49	52 - 56 cm 56 - 60 cm 60 - 64 cm

please indicate left or right

1 cm = 0.3937 inch

BORT Knee Support with Adjustable Hinge REF 114 460



PRODUCT-BENEFIT

- Knee support with side joint splints for guidance of the knee joint
- Extension limit at 0° and 10°
- Flexion limitation at 30°, 60° and 90°
- Circular fasteningstraps on the thigh and lower leg
- Plus-size system resulting in very good fitting

COLOUR: silver

THE MED-BENEFIT

- Adjustable joint for limiting extension and flexion
- Easily anatomically adjustable joint splints with effective side stabilisation
- Slim joint splint design

BORT Stabilo® Knee Support with Hinge REF 150 140



PRODUCT-BENEFIT

- Knee support with side joint splints for guidance of the knee joint
- Extension limit at 0°
- Circular fastening straps on the thigh and lower leg
- Pressure-reduced zones in the edge and popliteal fossa areas
- Plus-size system resulting in very good fitting
- For children too, see page 32

COLOUR: silver

THE MED-BENEFIT

- Easily anatomically adjustable joint splints with effective side stabilisation
- Slim joint splint design

INDICATIONS

Meniscus injuries, mild to moderate lateral ligament instability of the knee joint, gonarthrosis, arthritis, bursitis

INDICATIONS

Mild to moderate lateral ligament instability of the knee joint, gonarthrosis, arthritis

SIZES		d-measure	e f-measure	d-me	asure	f-measure
	Size 2+ Size 3	d-measure 23 - 24 26 - 29 29 - 33 29 - 33 33 - 33 33 - 33	40 - 44 cm 40 - 44 cm	Size 5+ 4 Size 6 4	7 — 41 1 — 45 1 — 45 5 — 49	48 - 52 cm 52 - 56 cm 56 - 60 cm 60 - 64 cm

please indicate left or right 1 cm = 0.3937 inch *Children size only available for REF 150 140



PRODUCT-BENEFIT

- Soft brace with movement limitation for guidance and relief of the knee joint
 Medial and lateral stabilisation
- Adjustable extension: 0°, 5°, 10°, 15°, 20°, 25°, 30°
- Adjustable flexion: 0° (180°), 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Joint splints cold-mouldable
- Breathable, elastic material for optimal comfort

COLOUR: grey

THE MED-BENEFIT

- With four circular Velcro straps for particularly secure support
- It can be fully opened so is very easy to put on
- Plus-size system resulting in very good fitting
- Tool-free setting of extension and flexion

BORT StabiloPro® Knee Brace REF 182 150

PRODUCT-BENEFIT

- Soft brace with limitation of mobility for guidance and relief of the knee joint
 With open patella
- Adjustable extension: 0°, 5°, 10°, 15°, 20°, 25°, 30°
- Adjustable flexion: 0° (180°), 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Joint splints cold-mouldable
- Breathable, elastic material for optimal comfort

COLOUR: 🔳 grey

THE MED-BENEFIT

- With two circular Velcro straps for additional support
- Plus-size system
- Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint
- Tool-free setting of extension and flexion

BORT INNOVATION:

 Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint

INDICATIONS

Moderate ligament instability of the knee joint, meniscus injuries, after-care in arthroscopic procedures involving the knee, chronic, postoperative and post-traumatic inflammations

INDICATIONS

Mild to moderate ligament instability of the knee joint, meniscus injuries, gonar-throsis

SIZES	d-measure	f-measure		d-measure	f-measure
Size 2 Size 2 Size 2 + Size 3 Size 3 + Size 4	26 - 29 29 - 33 29 - 33 33 - 37 33 - 37 37 - 41 37 - 41	36 - 40 cm 40 - 44 cm 40 - 44 cm 44 - 48 cm 44 - 48 cm	Size 5+ Size 6	41 – 45 41 – 45 45 – 49 45 – 49	56 - 60 cm 60 - 64 cm

right and left symmetrical 1 cm = 0.3937 inch



- Soft brace with movement limitation for guidance and relief of the knee joint
 Adjustable extension: 0°, 5°, 10°, 15°,
- Adjustable extension: 0°, 5°, 10°, 15°, 20°, 25°, 30°
- Adjustable flexion: 0° (180°), 10°, 20°, 30°, 45°, 60°, 75°, 90
- Joint splints cold-mouldable
- Breathable, elastic material for optimal comfort

THE GENERATION-BENEFIT

- Can be fully opened
- Extremely easy handling and putting on the brace
- Wide grip straps allow safe and energy-sparing handling also if gripping ability is restricted
- Popliteal opening facilitates flexion when sitting
- Extra-wide guide loops
- Attractive, functional design
- Tool-free setting of extension and flexion

COLOUR: 🗖 blue

BORT INNOVATION:

 Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint

INDICATIONS

Degenerative diseases, e. g. mild to moderate ligament instability of the knee joint, gonarthrosis, meniscus injuries

SIZES d	l-measure	f-measure		d-measure	f-measure
Size 1 Size 2 Size 2+ Size 3 Size 3+ Size 4 Size 4+	29-33 29-33 33-37 33-37 33-37 37-41	36 - 40 cm 40 - 44 cm 40 - 44 cm 44 - 48 cm 44 - 48 cm	Size 5+ Size 6 Size 6+ Size 7	41 - 45 41 - 45 45 - 49 45 - 49 45 - 49 49 - 53 49 - 53	56 - 60 cm 60 - 64 cm 64 - 68 cm 68 - 72 cm

right and left symmetrical 1 cm = 0.3937 inch

BORT Generation Soft OA Knee Brace REF 215 440



PRODUCT-BENEFIT

- Soft relief brace for isolated medial/ lateral gonarthrosis axial differenceup to 5° axial difference
- Adjustable, polycentric joint
- Limitation of extension: 0°, 5°, 10°, 15°, 20° and 25°
- Limitation of flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°
- Brace length: 39 cm

COLOUR: 🔳 blue

- Application of model A – Right leg, hinge is positioned lateral: For varus gonarthrosis (for valgisation in case of bowleg position)
- Left leg, hinge is positioned medial: For valgus gonarthrosis (for varisation in case of knock knee position)

THE GENERATION-BENEFIT

- Breathable material provides good wearing comfort
- Targeted relief of the medial/lateral compartment
- Varus / valgus correction freely adjustable
- Lightweight and slim structure for ideal wearing comfort



BORT Generation:

For better mobility for the aged Product range with special handling advantages

Application of model B

- Right leg, hinge is positioned medial: For valgus gonarthrosis (for varisation in case of knock knee position)
- Left leg, hinge is positioned lateral: For varus gonarthrosis (for valgisation in case of bowleg position)

INDICATIONS

Partial relief of the medial or lateral compartment in case of level II-III gonarthrosis, Bone bruise (post-traumatic bone marrow edema), medial/lateral gonarthrosis as a result of a tibial plateau fracture

SIZES		d-me	asure	f-measure
15 cm	Size 1 Size 2 Size 3 Size 4 Size 5		36 cm 42 cm 48 cm 54 cm 50 cm	56 – 62 cm
	Modell A Modell B			/ valgus left t / varus left

BORT OA-Xpress REF 101 600



PRODUCT-BENEFIT

- Knee relief brace for isolated medial / lateral gonarthrosis up to 5° axial deviation
- Varus / valgus correction freely adjustable
- Adjustable, polycentric joint
- Limitation of extension: 0°, 5°, 10°, 15°, 20° and 25°
- Limitation of flexion: 0° (180°), 30° 45°, 60°, 75° and 90°
- Brace length: 40 cm

COLOUR: grey

Application of model A

- Right leg, hinge is positioned lateral: For varus gonarthrosis (for valgisation in case of bowleg position)
- Left leg, hinge is positioned medial:
 For valgus gonarthrosis (for varisation in case of knock knee position)

Application of model B

THE MED-BENEFIT

ideal wearing comfort

compartment

- Partially elastic shells for good fitting

- Targeted relief of the medial / lateral

to different leg circumferences

- Lightweight and slim structure for

- Right leg, hinge is positioned medial:
 For valgus gonarthrosis (for varisation in case of knock knee position)
- in case of knock knee position) - Left leg, hinge is positioned lateral: For varus gonarthrosis (for valgisation in case of bowleg position)

INDICATIONS

Partial relief of the medial or lateral compartment in case of level II-III gonarthrosis, Bone bruise (post-traumatic bone marrow edema), medial/lateral gonarthrosis as a result of a tibial plateau fracture

IMPLEMENTATION

Model A	varus right/valgus left
Model B	valgus right/varus left



PRODUCT-BENEFIT

- Frame knee orthosis by the 4-point principle for guidance and stabilisation of the knee joint
- Limitation of extension: 0°, 5°, 10°, 15°. 20° and 25°
- Limitation of flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

COLOUR: qrey

THE MED-BENEFIT

- High frame strength combined with low self-weight thanks to ultramodern material
- Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint
- Tool-free setting of extension and flexion

INDICATIONS

Serious and/or complex instability of the knee joint, meniscus injuries, genu recurvatum, functional pre-and/or postoperative treatment of ruptured ligaments (ACL, PCL, collateral ligaments)



PRODUCT-BENEFIT

- Frame knee orthosis by the 4-point principle for guidance and stabilisation of the knee joint
- Limitation of extension: 0°, 5°, 10°, 15°, 20° and 25°
- Limitation of flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

COLOUR: grey

THE MED-BENEFIT

- High frame strength combined with low self-weight
- Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint
- Self-adjusting tibia fitting
- Tool-free setting of extension and flexion

INDICATIONS

Serious and/or complex instability of the knee joint, meniscus injuries, functional pre-and/or post-operative treatment of ruptured ligaments (ACL, PCL, collateral ligaments)

SIZES		f-measure	Width of knee
With the second	medium plus large large plus x-large x-large plus	47-54 47-54 54-61	9,5-10,5 cm 9,5-10,5 cm 10,5-11,5 cm 10,5-11,5 cm 11,5-12,5 cm 11,5-12,5 cm

please indicate left or right 1 cm = 0.3937 inch



- Simple knee support for easy stabilisation
- Anatomically contoured knit
- Compression class II



PRODUCT-BENEFIT

- Simple round knitted knee support for easy stabilisation
- Compression class II

COLOURS: beige black

THE MED-BENEFIT

- Easy to use
- Not bulky
- Short design

COLOURS: beige blue

THE MED-BENEFIT

- Easy to handle

- Not bulky

INDICATIONS

Sports injuries, articular effusion, sprains, contusions

INDICATIONS

Sports injuries, articular effusion, sprains, contusions

SIZES

Circumference of knee

4	small medium large	 up to 32 cm 32 – 37 cm 37 – 42 cm
1	x-large	 42 – 47 c

SIZES	C

Circumference of knee

m la x-	nall edium rge large :-large	 $\begin{array}{c} \text{up to 31 cm} \\ 31-37 \text{ cm} \\ 37-42 \text{ cm} \\ 42-48 \text{ cm} \\ 48-54 \text{ cm} \end{array}$



PRODUCT-BENEFIT

- Knee brace for immobilisation of the knee joint with a 0° - 20° flexion
- 3-piece with open patella
- Medial and lateral aluminium splints, which can be freely positioned
- For children too, see page 33

COLOUR: 🗖 black

THE MED-BENEFIT

- The support adapts to body shape thanks to the 3-piece design
- Lined with terry cloth on the inside for excellent comfort
- The aluminium splints can be removed for washing

PRODUCT-BENEFIT

- Knee brace for mobilisation
- Brace length can be set using lateral telescopic splints
- Limitation of flexion: 0°, 10°, 20°, 30°, 40°, 50°, 60°, 70°, 80°, 90°, 100°, 110° and 120°
- Limitation of extension:
 -10°, 0°, 10°, 20°, 30°, 40°, 50°, 60°, 70°, 80° and 90°
- Immobilisation: -10°, 0°, 10°, 20°, 30°

COLOUR: 🔳 grey

THE MED-BENEFIT

- Setting joint limit stops and brace length without tools
- Lightweight and slim structure for ideal wearing comfort
- easyClick fasteners make it easy to put on/take off
- Easy to handle

INDICATIONS

Pre-operative, post-operative, posttraumatic, joint, muscle or ligament injuries

SPLINT LENGTH

Fits height

Children (38 cm)	80 – 145 cm
45 cm	145 – 155 cm
50 cm	155 – 165 cm
55 cm	165 – 180 cm
60 cm	
Please indicate 0° or 20	° flexion

right and left symmetrical 1 cm = 0.3937 inch

INDICATIONS

Postoperative, posttraumatic early functional mobilisation with restricted movement of the knee joint, Condition after damage to bone/capsule/ligament structures of the knee joint

SIZE

One size only

right and left symmetrical 1 cm = 0.3937 inch



- Knee brace for immobilisation of the knee joint with a 0° - 20° flexion
- 3-piece with open patella
- Medial and lateral aluminium splints, which can be freely positioned

BORT adjustable Immobilisation Splint REF 145 080



PRODUCT-BENEFIT

- Knee brace for immobilisation of the knee joint
- Adjustable, 0° or 15° flexion
- With open patella, 1-piece
- Medial and lateral aluminium splints
- Length: 55 cm (21.7 inches)

COLOUR: Dlack

THE MED-BENEFIT

- Variable adaption thanks to 3-piece design
- On the inner side of the lower leg with an extensive anti-slip zone
- The aluminium splints can be removed for washing

COLOUR: Dlue/black

THE MED-BENEFIT

- Easy to handle
- Lined with terry cloth on the inside for excellent comfort
- The aluminium splints can be removed for washing

BORT Immob Splint Single-Piece REF 145 070

 Knee brace for immobilisation of the knee joint 0° or 15° flexion

INDICATIONS

Pre-operative, post-operative, post-traumatic, joint, muscle or ligament injuries

Fits height

40 cm _____ 155 – 165 cm 50 cm _____ 165 – 180 cm 60 cm _____ over 180 cm SIZE

One length only, 55 cm

For REF 145 070 please indicate 0° or 15° flexion

right and left symmetrical 1 cm = 0.3937 inch

Please indicate 0° or 20° flexion

right and left symmetrical 1 cm = 0.3937 inch

Sure footed: with all-round protection and stability at every step.



Application area				Upper	Upper ankle joint					Lower ankle joint		Foot	
Overview of indications for selected products	select* TaloStabil® 054 700 select* TaloStabil® Plus 054 100	select* Achillo Stabil® Plus 053 100	TaloStabil® Sport 054 600 SP TaloStabil® Plus Sport 054 100 SP	Achillo Stabil® Plus Sport 053 100 SP	Generation Ankle Support 215 600	Generation TaloStabil* Ankle Eco Support 054 600 215 600 TaloStabil* Eco Plus 054 650	Achillo Stabil® Eco 054 900	Stabilo® Ankle Support 114 100	TaloFX® 054 800 Talo Talo Xpress 053 900	Spiratdynamik. 053 200	Malleo Stabil≋ Brace 100 500	Malleo Xpress Active 100 580 Malleo Xpress 100 560	Foot Levator Orthotic Device 054 300
Page	P. 158/159	P. 158	P. 160 / 159	P. 160	P. 161	P. 162/163	P. 163	P. 167	P. 167/170	P. 166	P. 169	P. 168	P. 170
Ligament lesions	•		•		•	•		•	•		•	•	
Distorsion of upper ankle joint	•		•		•	•		•	•		•	•	
Achillodynia		•		•			•						
Arthrosis	•		•		•	•							
Supination/Pronation Trauma									•		•	•	
Tendinitis								•					
Flaccid peroneal paralysis													•
Paralysis of the dorsiflexors and the big-toe-lift													•
Foot problems due to static ano- malies and defects in the shape of the foot due to dysfunction of the lower ankle joint, e.g. skew and splayfoot										•			

OVERVIEW OF INDICATIONS 157

Ankle Foot



PRODUCT-BENEFIT

- High quality active ankle support
- Highly elastic stretch zone in the sensitive instep area simplifies putting on and removal
- Two structured silicone pads for intermittent massage



THE SELECT-BENEFIT

- SOFTflex Plus knit fabric

- SENSIL[®] Breeze: Luxury fibre with a cooling effect
- Soft-velvety surface structure for top wearing comfort
- Extra soft, snag proof material edge at the calf
- Separate slip-on aid included

INDICATIONS

Chronic, post-traumatic or postoperative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, tendomyopathies

SIZES	
-------	--

Circumference of ankle

breasure	up to 21 cm
breasure	21 – 23 cm
large	23 – 25 cm
x-large	25 – 27 cm
xx-large	27 – 29 cm

please indicate left or right 1 cm = 0.3937 inch





PRODUCT-BENEFIT

- High-quality active Achilles tendon support for compression of soft tissue
- Highly elastic stretch zone in the sinstep area
- Extended, structured silicone pad for intermittent massage of the achilles tendon and the muscletendon transitional area towards the calf muscle

COLOURS: silver black

THE SELECT-BENEFIT

- Its half-open structure makes it easier to put on and take off
- Easy and simple fixing of the support using three Velcro fasteners
- Heel wedge to relieve the Achilles tendon during the acute phase

INDICATIONS

SI

Chronic or post-traumatic irritation or tenderness in the achilles tendon area, haglund syndrome heel irritation, achillodynia, irritation and pain of the achilles tendon and surrounding areas, after surgery, e. g. achilles tendon rupture, bursitis subachillea

ZES	Circumference	of ankle

b-measure	small medium large x-large xx-large	 up to 21 cm 21 – 23 cm 23 – 25 cm 25 – 27 cm 27 – 29 cm
	xx-large	 27 – 29 cm

right and left symmetrical 1 cm = 0.3937 inch



- High quality active ankle support with Stabilo[®] Strap, partially elastic, it can be freely positioned
- Highly elastic stretch zone in the sensitive instep area simplifies putting on and removal
- Treatment-supporting positioning of the foot in the pronation / supination / neutral position
- Two structured silicone pads for intermittent massage

COLOURS: silver black

THE SELECT-BENEFIT

- SOFTflex Plus knit fabric:

- \cdot SENSIL® Breeze: Luxury fibre with a cooling effect
- · Soft-velvety surface structure for top wearing comfort
- Extra soft, snag proof material edge at the cal
- Separate slip-on aid included

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

SIZES Circumfe	rence of ankle
----------------	----------------

small up to 21 cm medium 21 – 23 cm large 23 – 25 cm x-large 25 – 27 cm xx-large 27 – 29 cm

please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Active ankle support for compression of soft tissues, with Stabilo[®] Strap, partially elastic can be freely positioned
- SOFTflex 3D knitted fabric with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II

COLOUR: black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body, for dryness and comfort
- Treatment-supporting positioning of the foot in the pronation/supination/ neutral position
- Ideal fit for sports activities
- Latex-free knitted fabric

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

SIZES

Circumference of ankle

b-measure	small _	up to 21 cm
	medium _	21 – 23 cm
لى مەلىرىكى ئەرىپىلى ئەرىپىلى ئەرىپىلى ئەرىپىلى ئەرى	large _	23 – 25 cm
with the second s	x-large _	25 – 27 cm

right and left symmetrical 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Ankle support for compression of soft tissues
- SOFTflex 3D knitted fabric with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II

COLOUR: I black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body, for dryness and comfort
- Ideal fit for sports activities
- Latex-free knitted fabric
- Modern, sporty design (with a reflecting label)

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, arthrosis

SIZES

Circumference of ankle

b-measure	small medium large x-large xx-large	 up to 21 21 – 23 23 – 25 25 – 27 27 – 29	cm cm cm

right and left symmetrical 1 cm = 0.3937 inch

BORT AchilloStabil® Plus Sport REF 053 100 SP



COOLMAX

PRODUCT-BENEFIT

- Active Achilles tendon support
- SOFTflex 3D knitted fabric with a velvety soft surface structure very comfortable to wear
- Extended structured silicone pads for intermittent massage of the Achilles tendon
- Extra soft, snag proof material edge at the calf

COLOUR: I black/green

THE SPORT-BENEFIT

- COOLMAX[®]: wicks moisture away from the body
- Heel wedge to relieve the Achilles tendon during the acute phase
- Integrated comfort zone thanks to Dynamic-Flex-fabric
- Separate slip-on aid included
- Latex-free knitted fabric

INDICATIONS

SIZES

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, in the case of Haglund syndrome heel irritation, achillodynia, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

Circumference of ankle

right and left symmetrical 1 cm = 0.3937 inch

BORT Generation Ankle Support REF 215 600



PRODUCT-BENEFIT

- High-quality active ankle support for compression of soft tissue
- Highly elastic stretch zone in the sensitive instep area simplifies putting on and removal
- Functional knitted fabric, very comfortable to wear
- Two structured silicone pads for intermittent massage

COLOUR: 🔳 blue

THE GENERATION-BENEFIT

- Easy to put on/remove thanks to open design, also with restricted movement
- Easy and simple fixing of the support using two Velcro fasteners
- Attractive, functional design



BORT Generation: For better mobility for the aged Product range with special handling advantages

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

SIZES	
-------	--

Circumference of ankle

	nall	up to 21 cm 21 – 23 cm
ے ا	rge	23 – 25 cm
	large -large	25 – 27 cm 27 – 29 cm
1.1.		



PRODUCT-BENEFIT

- Active ankle support for compression of soft tissues
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Two silicone pads for intermittent massage
- Compression class II
- Also for children, see page 34

COLOURS: beige silver black

THE ECO-BENEFIT

- Velvety soft surface structure
- Pressure relief at the base of the ankle
- Plus size system resulting in very good fitting

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

SIZES		b-measure	b1-measure**
b-measure	Children (blue) small medium large x-large xx-large xx-large plus* xxx-large plus*	$\begin{array}{c} 17 - 19 \text{ cm} \\ \text{up to } 21 \text{cm} \\ 21 - 23 \text{ cm} \\ 23 - 25 \text{ cm} \\ 25 - 27 \text{ cm} \\ 27 - 29 \text{ cm} \\ 27 - 29 \text{ cm} \\ 29 - 31 \text{ cm} \\ 29 - 31 \text{ cm} \end{array}$	32 – 34 cm 34 – 36 cm 36 – 38 cm 38 – 40 cm

right and left symmetrical

*only in colour silver

1 cm = 0.3937 inch ** the b1 measurement is required for differentiation of the Plus size



- Active ankle joint support with Stabilo[®] strap, partially elastic, freely positionable
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Treatment-supporting positioning of the foot in pronation, supination or neutral position
- Two silicone pads for intermittent massage
- Compression class II

COLOUR: silber

THE ECO-BENEELT

- Soft-velvety surface structure
- Pressure-reduced edge in the ankle area

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

SIZES	Circumference of ankle		
b-measure	small medium large x-large xx-large xxx-large		bis 21 cm 21 – 23 cm 23 – 25 cm 25 – 27 cm 27 – 29 cm 29 – 31 cm



PRODUCT-BENEFIT

- Achilles tendon support
- SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- Viscoelastic pad alongside the Achilles tendon for intermittent massage
- Heel wedge to relieve the Achilles tendon during the acute phase
- Compression class II

COLOURS: beige silver black

THE ECO-BENEELT

- Velvety soft surface structure
- Reduced pressure in the ankle area

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area. Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

SIZES	Circumference of ankle	
b-measure	small medium large x-large xx-large	_ up to 21 cm _ 21 – 23 cm _ 23 – 25 cm _ 25 – 27 cm _ 27 – 29 cm



PRODUCT-BENEFIT

- Anatomically knitted ankle joint support for soft tissue compression
- Two silicone pads for intermittent massage
- Microfiber knitted fabric, soft to the touch
- Compression class II
- COLOURS: beige mineral grey

THE MED-BENEFIT

- Reduced pressure in the ankle area
- Anatomical knitted heel for optimal fit

INDICATIONS

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, sprain, ligament weakness, after removal of plaster cast, ligament injuries, arthrosis

C	17	EC
J	12	LJ

Circumference of ankle

large 23 - 25 c	m la	nedium 2 arge 2	p to 21 cm 21 – 23 cm 23 – 25 cm 25 – 27 cm
-----------------	---------	--------------------	--

right and left symmetrical 1 cm = 0.3937 inch



- Anatomically contoured ankle support with Stabilo[®] Strap, partially elastic, freely positionable
- Treatment-supporting positioning of the foot in pronation, supination or neutral position
- Two silicone pads for intermittent massage
- Microfibre knitted fabric, soft grip
- Compression class II

COLOUR: graphit

THE MED-BENEFI

- Pressure-reduced edge in the ankle area
- Anatomically contoured heel with a good fit

INDICATIONS

SIZES

Chronic, post-traumatic or post-operative soft tissue irritations in the ankle joint area, arthrosis, ligament weakness, after removal of plaster cast, ligament injuries, distortion of upper ankle joint, supination prophylaxis

BORT activemed Achilles Tendon Support REF 220 900



PRODUCT-BENEFIT

- Achilles tendon support
- Viscoelastic pad alongside the Achilles tendon for intermittent massage
- Heel wedge to relieve the Achilles tendon during the acute phase
- Compression class II

COLOURS: beige mineral grey graphite

THE MED-BENEFIT

- Reduced pressure in the ankle area

- Anatomical knitted heel for optimal fit

INDICATIONS

Chronic or post-traumatic irritation or tenderness in the Achilles tendon area, Haglund syndrome heel irritation, achillodynia, irritation and pain of the Achilles tendon and surrounding areas, after surgery, e. g. Achilles tendon rupture, bursitis subachillea

Circumference of ankle

b-measure small	up to 21 cm
medium	21 – 23 cm
large	 23 – 25 cm
x-large	 25 – 27 cm





PRODUCT-BENEFIT

- Ankle bandage to support torsion of the back of the foot and to raise the instep
- Ligament separation thanks to longitudinal cuts for an anatomical fit and optimal function
- Thin, elastic straps
- 3-D contoured knit
- Two silicone pads
- Extra soft, snag proof material edge at the calf
- Also for children, see page 33

COLOUR: black

THE MED-BENEFIT

- Stabilisation of the ankle joint and activation of foot muscles through compression
- Functional supportive restraint in the case of classical foot deformities (i.e. skew and splayfoot)
- Optimal comfort even during intense activity
- Intermittent massage
- Simplified handling
- Separate slip-on aid included
- Pack: single item

INDICATIONS

Chronic, post-traumatic or post-operative irritations of soft tissues in the area of lower ankle joint and the sole of the foot, e. g. arthrosis, ligament instability, foot problems due to static anomalies and defects in the shape of the foot due to dysfunction of the lower ankle joint e. g. skew and splayfoot

SIZES	Circ	umference of ankle
small medium large x-large		



- Ankle brace with laces for stabilisation of the lower and upper ankle joints
- Safe handling due to variable hook and eyelet elements
- Individually adjustable stabilisation level
- Plastic U-shaped stabilisation element
- Non-elastic breathable material

COLOUR: 🔳 grey

- THE MED-BENEFIT
- Stabilo[®] Strap for additional stabilisation and limitation of talus prolapse
- Can be worn with sports and walking shoes

INDICATIONS

Post-operative protection after ligament suture/reconstruction, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries

SIZES	Circumference of ankle		
b-measure	x-small up to 19 cm small 19 - 21 cm medium 21 - 23 cm large 23 - 25 cm x-large 25 - 27 cm		

right and left symmetrical 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Soft ankle brace for stabilisation of the lower and upper ankle joints
- Bilateral shells padded with soft memory foam
- Cross restraint featuring Dynamic-Force-Control (DFC) technology
- Elastic material that carries away moisture
- Brace can be opened fully from the rear for easy access
- Also for children, see page 35

COLOUR: 🗖 black

THE MED-BENEFIT

- Slim design

- Can be worn without difficulty with sports and walking shoes
- Limitation of talus prolapse with optimal stabilisation
- For positioning the ankle in pronation or supination

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries, supination/pronation trauma

SIZES	Circumference of ankle
b-measure	Children 17 – 19 cm small 19 – 21 cm medium 21 – 23 cm large 23 – 25 cm x-large 25 – 27 cm

please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Soft ankle joint brace for mobilisation and safe stabilisation of the ankle joint
- Modular structure, can be downgraded during the course of therapy
- Reduction of the talus prolapse
- Bilateral shells padded softly with memory foam
- Additional stabilisation by retrievable soleplate in the acute phase

COLOUR: I black/bordeaux

THE MED-BENEFIT

- Can be fully opened, therefore easy to put on
- Can be worn without difficulty with sports and walking shoes
- Easy to handle and stabilise without graduation using with BOA® Fit System

INDICATIONS

Post-operative, post-traumatic, e. g. conservative treatment of ankle distortions (level II-III) and ligament ruptures, immobilisation in the case of inflammatory conditions, post-operative protection after ligament suture/reconstruction, chronic instabilities of the upper and lower ankle joints, acute and chronic capsular ligament instabilities

SIZES

Circumference of ankle

b-measure	small medium large x-large	 19 - 21 - 23 - 25 -	23 25	cm cm

please indicate left or right 1 cm = 0.3937 inch

BORT MalleoXpress Active REF 100 580



PRODUCT-BENEFIT

- Soft ankle joint brace for stabilisation of the ankle joint
- Reduction of the talus prolapse by intelligent track control
- Bilateral shells padded softly with memory foam

COLOUR: black

THE MED-BENEFIT

- Can be fully opened, therefore easy to put on
- Can be worn without difficulty with sports and walking shoes
- Easy, safe handling
- Easy to handle and stabilise without graduation using with BOA® Fit System

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions and ligament ruptures, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, ligament injuries, supination/pronation trauma

Circumference of ankle

b-measure	small medium large x-large	_ 19 – 21 cm _ 21 – 23 cm _ 23 – 25 cm _ 25 – 27 cm
-----------	-------------------------------------	--

please indicate left or right 1 cm = 0.3937 inch



- Combined Air/Gel brace and foam pad to stabilise the ankle, while limiting pronation and supination
- Anatomically shaped plastic shells
- Two sets of padded variants:
- Pre-coolable Air/Gel pad to be used during the acute phase
- Memory foam padding for comfortable wear during the rehabilitation phase
- Also for children, see page 34

COLOUR: grey

THE MED-BENEFIT

- No impairment when walking
- Easy to put on thanks to circular Velcro fastening

Other product versions:





- With Air/Gel padding only

INDICATIONS

Post-operative protection after ligament suture, conservative treatment of ankle distortions of any severity, functional treatment of ligament rupturesankle fractures (Weber A), chronic ligament instability, shortening of time spent in plaster, ligament injuries, distortion of upper ankle joint, fibular ligament rupture, pronation and/or supination traumas

Ankle

SIZES	Fits height
Size 1/Children	up to 155 cm
Size 2	over 155 cm



PRODUCT-BENEFIT

- Ankle joint brace for stabilisation of the upper and lower ankle joint
- Stabilisation degree adjustable using restraints
- Very thin and slim construction
- Can be worn in shoes without any problems
- Effective stabilisation of the ankle joint and very comfortable to wear

COLOUR: grey

THE MED-BENEFIT

- Can be fully opened, therefore easy to put on and simplest handling
- High stabilisation combined with low self-weight thanks to ultra-modern material
- Lightweight and slim structure for good wearing comfort

INDICATIONS

Conservative treatment of ankle distortions and ligament ruptures, post-operative protection after ligament suture/ reconstruction, immobilisation in the case of inflammatory conditions, rheumatism, arthritis, tendinitis, supination/ pronation trauma

SIZES	h-measure
h-measure h-measure large	

Please indicate left or right 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Dynamic textile Tibialis Anterior Muscle brace
- 8-strap system
- Non-elastic lateral strapping
- Supports flexion of the foot in the upper ankle joint
- Stabilisation of the lower ankle joint
- Pronation effect

COLOUR: Dlack

THE MED-BENEFIT

- Easy to put on thanks to labelling of the bracing elements
- Can be applied using one hand
- Two support lengths towards the forefoot for optimal flexion

INDICATIONS

Peroneal paralysis of any severity, especially for mild to medium paralysis of the tibialis anterior muscle and the flexor hallucis longus (degree of severity 2-4)

SIZES

Circumference of ankle

b-measure	x-small small medium large x-large	 17 – 19 cm 19 – 21 cm 21 – 23 cm 23 – 25 cm 25 – 27 cm
	short long	shoe size 6 shoe size 6

Indicate right/left, short/long 1 cm = 0.3937 inch



- Polypropylene tibialis anterior muscle brace to flex the foot during the swing phase
- ¾ sole up to the metacarpophalangeal joints of the big toe



PRODUCT-BENEFIT

- Polypropylene tibialis anterior muscle brace to flex the foot during the swing phase
- Long sole design

COLOUR: White

THE MED-BENEFIT

- Thermoplastic reprocessing
- Can be used in closed shoes
- Highly flexible semi-shell in the calf area
- Velcro fastening for individual adjustment to suit different calf circumferences

INDICATIONS

Flaccid paralysis of the tibialis anterior muscle, reduced mobility of the patient, in the case of peroneal paralysis, lumbar nerve root damage, polio, stroke

COLOUR: White

THE MED-BENEFIT

- Thermoplastic reprocessing
- Can be used in closed shoes
- Velcro fastening for individual adjustment to suit different calf circumferences
- Light weight
- Easy to clean

INDICATIONS

Flaccid paralysis of the tibialis anterior muscle, reduced mobility of the patient, in the case of peroneal paralysis, lumbar nerve root damage, polio, stroke

SIZES	Shoe size UK
small	3 – 5.5
medium	6 – 8.5
large	9 – 13

SIZES	Sho	e size UK
small medium large x-large xx-large		2 - 4 4 - 5.5 5.5 - 7 7 - 9 9 - 11.5

please indicate left or right



PRODUCT-BENEFIT

- Orthopaedic shoe, closed design
- Anti-slip sole
- Breathable materials
- Pack: Single item

COLOUR: Dlack

THE MED-BENEFIT

- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual
- Can be fully opened
- Easy to put on
- Comfortable to wear

BORT Therapeutic Shoe, open design REF 930 180

PRODUCT-BENEFIT

- Orthopaedic shoe, open
- Anti-slip sole
- Breathable materials
- Pack: Single item

COLOUR: Dlack

THE MED-BENEFIT

- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual
- Can be fully opened
- Easy to put on
- Comfortable to wear

ACCESSORIES: REF 990 970 BORT Toe Cap

INDICATIONS

Post-operative and other wounds, even when treated with extensive and possibly padded bandaging, where the healing process takes place over an extended period of time and weight can be placed partially on the foot whilst walking, even outdoors. An additional insole is not required

Sh	oe size UK
	1.5 – 2
	3 – 4
	4.5 - 5.5
	6-7
	7.5 – 8.5
	9 – 10
	. 11 – 11.5

Shoe size UK

 1.5 – 2
 3-4
 4.5 - 5.5
 6-7
 7.5 – 8.5
 9 – 10
 11 – 11.5

right and left symmetrical

right and left symmetrical



- Orthopaedic shoe
- Extra-wide fit at the forefoot and toe area
- Zipper at the heel
- Flexible, anti-slip sole
- Pack: Pair or single item (REF) 930 320 R or 930 320 L)

COLOUR: Dlack

THE MED-BENEFIT

- Secure hold and comfort for sensitive feet
- Width can be individually adjusted thanks to Velcro fastening
- Can be fully opened, comfortable walking

INDICATIONS

Post-operative and other wounds, also being treated with large and possibly padded bandages whose healing process take place over an extended period of time and weight can be placed partially on the foot whilst walking, even outdoors. An additional insole is not required

	Sh	ioe siz	e UK
	. 3		7.5
	. 4		8.5
	4.5		9
	5.5		10
	. 6		11
	. 7		11.5
Supplied as pairs or single	item		

please indicate left or right

BORT Diabetic Orthopaedic Shoe REF 930 400 R/930 400 L



PRODUCT-BENEFIT

- Orthopaedic shoe, open with an optionally removable forefoot/toe cap
- Removable soft pad cushioning with a honeycomb structure
- Reinforced anti-slip sole with foot rolling motion aid
- Pack: Single item

COLOUR: 🗖 black

THE MED-BENEFIT

- Soft pad cushioning with individually removable honeycomb elements for local pressure relief
- Additional protection by the removable forefoot/toe cap
- No inner seams

INDICATIONS

Wounds in the foot and ankle region due to diabetes mellitus, peripheral artery occlusive disease (PAOD), further wounds in the foot and ankle region, post-traumatic, post-operative

Shoe size UK

			2 – 3
			3.5 – 4.5
 			5-6
 			6.5 – 7.5
			8 – 9
 			9.5 - 10.5
 	 	 	11 – 12
·	1.4		

please indicate left or right



PRODUCT-BENEFIT

- Forefoot relief shoe
- Anti-slip sole
- With additional protective toecap
- Breathable materials

- Pack: Single item

COLOUR: black

THE MED-BENEFIT

- Redistribution of pressure towards the middle or back of the foot (heel)
- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual



PRODUCT-BENEFIT

- Heel relief shoe
- Anti-slip sole
- Concave heel edge
- Breathable materials
- Pack: Single item

COLOUR: black

THE MED-BENEFIT

- Redistribution of pressure towards the middle or forefoot
- Improved stability during walking
- The foot bedding and EVA rubber sole can be easily adjusted and re-adjusted to suit any individual

INDICATIONS

Post-operative, wounds/injuries to the rear foot area at the sole of the foot. post-operative following corrective osteotomy in the forefoot area

INDICATIONS

Post-operative, wounds/injuries to the rear foot area and to the heel

Shoe size UK 1 - 23 - 4.555 - 77.5 - 910 - 11.5

Shoe size UK

 2-4
 4.5-6
 7 - 9
 9 - 11

right and left symmetrical

right and left symmetrical

SHOES/FOOT/ANKLE 175



PRODUCT-BENEFIT

- Elastic bandage support with pad for splayfoot
- Pad inserted in a slit pocket
- Support and lifting of the transverse arch of the metatarsal bones
- Pack: Pair

COLOUR: beige

THE MED-BENEFIT

- Very easy to use
- Removable pad



PRODUCT-BENEFIT

- Simple ankle support
- Anatomical design
- Closed heel
- Compression class II

COLOURS: beige blue

THE MED-BENEFIT

- Very thin, comfortable to wear
- Seamless circular knitted fabric

INDICATIONS

Chronic fallen arch disorders, fallen arches and splayfoot, metatarsalgia

INDICATIONS

Ankle joint injuries, after removal of plaster cast, tendon irritations

Circumference in cm around forefoot

16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26 SIZES

Circumference of ankle

b-measure	small _ medium _ large _ x-large _	 up to 2 21 – 2 23 – 2 25 – 2	3 cm 5 cm

Supplied as pairs

right and left symmetrical



PRODUCT-BENEFIT

- Textile Hallux Valgus Brace for corrective or postoperative stabilisation of the large toe and the metatarsophalangeal joint
- The side aluminium reinforcement splint can be replaced by a spiral spring rod during the course of therapy
- Plantar splayfoot pad
- Heel strap for securing positioning
- Can be opened fully using the Velcro fastener
- Additional toe straps for correction of the 2nd - 5th toes

COLOUR: Dlack

THE MED-BENEFIT

- Easy to adapt, also on bandages
- Very comfortable to wear due to soft padded material
- Low material thickness and individual adjustment options
- Can also be worn during the daytime and at night

BORT Valco® Hallux Valgus Splint REF 930 010

PRODUCT-BENEFIT

- Hallux Valgus Brace for physiological alignment of the metatarsophalangeal joint of the big toe
- Anatomical fit
- Three adjustable Velcro straps
- For use during periods of rest
- Pack: Single item

COLOUR: white/grey

THE MED-BENEFIT

- Simple and precise adjustment of the leverage force
- Inner surface of the splint entirely lined with foam
- Ventilation hole in the big toe area


INDICATIONS

Hallux valgus, after hallux surgery to ensure the success of treatment, postoperative capsule relief

SIZES Circumference in cm around forefoot small ______16 - 19 cm medium ______19 - 22 cm large _____22 - 25 cm

SIZES	Shoe size UK
small medium large	up to 4.5 4.5 – 7 over 7

please indicate left or right 1 cm = 0.3937 inch please indicate left or right



- Soft brace for corrective or postoperative stabilisation of the large toe and the metatarsophalangeal joint
- Stepless setting of the corrective pressure using the laterally positioned corrective tape
- Relief pad on the same level as the metatarsophalangeal joint

BORT Silicone Heel Cushion [REF] 930 050

PRODUCT-BENEFIT

- Silicone heel cushion for shock absorption
- Reduction of impact load on the ankle, knee and hip joints and the spine
- Shell-shaped heel
- Height: 6 mm (0.2 inches)
- Pack: Pair

COLOUR: 🗖 blue

THE GENERATION-BENEFIT

- Low material thickness and individual adjustment options
- Can also be worn during the daytime and at night
- Comfortable to wear due to thin, soft material
- Simple and easy to put on
- Attractive, functional design

INDICATIONS

Hallux valgus, after hallux surgery to ensure the success of treatment, post-operative capsule relief

SIZES	Shoe size UK
Size 0 Size 1	4-6 7-9
Size 2	10 – 11

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Achillodynia, Haglund's deformity, heel disorders, shock absorption in the case of arthrosis of the ankle, knee or hip joint

SIZES	Shoe size UK
small	up to 4.5
medium	5.5 – 7
large	7.5 – 9
x-large	10 - 11.5
Cupplied on pair	

Supplied as pairs

Foot nsoles



PRODUCT-BENEFIT

- Viscoelastic spur cushion with removable padding module
- Reduction of impact load on the ankle, knee and hip joints and at the spine
- Upper side is lined, under side is self-adhesive
- Height: 6 mm (0.2 inches)
- Pack: Pair

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

BORT Silicone Heel Spur Cushion with SoftSpot REF 930 070

PRODUCT-BENEFIT

- Silicone heel spur padding with Soft-Spot insole
- Reduction of impact load on the ankle, knee and hip joints and the spine
- Shell-shaped heel
- Soft anatomical insole for the heel spur
- Height: 6 mm (0.2 inches)
- Pack: Pair

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Heel spur, achillodynia, Haglund's deformity, axial shock absorption in the case of leg joint arthrosis

SIZES Shoe size UK A 0-2 B 3-6 C 7-9 D 10-13

INDICATIONS

Heel spur, achillodynia, Haglund's deformity, shock absorption in the case of arthrosis of the ankle, knee or hip joint

SIZES	Shoe size UK
small medium large x-large	up to 4.5 5.5 - 7 7.5 - 9 10 - 11.5

Supplied as pairs

Supplied as pairs



- Viscoelastic silicone heel lift to compensate differences in leg-length
- Upper side is lined, under side is self-adhesive
- Heights: 5 mm, 10 mm (0.2 and 0.4 inches)



PRODUCT-BENEFIT

- Silicone heel lift to compensate differences in leg-length
- Made of durable silicone
- Heights: 3 mm, 5 mm, 8 mm, 10 mm (0.1, 0.2, 0.3 and 0.4 inches)
- Also for children, see page 35

- Pack: Single item

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

- Pack: Single item

THE MED-BENEFIT

- Non-slip
- Fits any kind of shoe
- Washable
- Hygienic

INDICATIONS

Compensation of differences in leg-length, achilles tendon disorders (Achillodynia)

SIZES	Shoe size UK	SIZES	Shoe size UK
small large		Children small medium large	6 - 8.5
Heights	5 mm / 10 mm	Heights _3 mm/5	mm / 8 mm / 10 mm
right and left symmetrical		right and left symmetrical	

Quality brands for individual all-round care.

Products for rehabilitation



PRODUCT BENEFIT

- Simplifies putting on medical compression stockings
- Safe, independent handling
- Can also be used in case of reduced strength and restricted mobility

COLOUR: white/grey

THE AKTIVEN® BENEFIT

- Non-slip due to secure fixing
- Top product quality
- Attractive, functional design

GOLOON. La Winte/gr

INDICATIONS

Ailments or disabilities with considerable functional disorders of the finger, hand or arm joints, primarily the gripping function which make the correct putting on of compression stockings prescribed by a physician impossible without aids, e.g. in case of highly inflammatory or considerably degenerative joint diseases, paralysis, deforming diseases in the hand region, results of injury, malformations, amputations; special sock or stocking donning aids for compression stockings with grip extensions are indicated due to the impossibility to reach to the stockings with the arms, e.g. in cases of far-reaching spinal, hip and knee joint stiffness as well as obesity

SIZES

Size 1

Size 2

Circumference of calf

< 40 cm

>40 cm



- 100 % silicone therapy putty
- 6 different types of consistency
- Latex-free
- Box sizes: 80 g, 440 g

THE MED-BENEFIT

- Strengthening of muscles in the fingers, hand and forearm
- Comprehensive, therapeutic finger and hand exercises
- Exercise instructions with many examples available in each box

INDICATIONS

For regaining mobility of the hand

BORT Codyn[®] Seat Cushion REF 190 530



PRODUCT-BENEFIT

- Special ring cushion with recess in the coccyx area
- Inflatable, via valve
- Supplied with manual air pump
- Size: approx. 45 cm × 45 cm,
- (17.7 × 17.7 inches) – Inner diameter approx. 16 cm (6.3 inches)

COLOUR: blue

THE MED-BENEFIT

- Suitable for all types of seating
- Inflatable, therefore easy to carry and can be used anywhere

INDICATIONS

Coccygodynia, post coccyx contusion and fracture conditions, perianal relief of haemorrhoids, post low rectal surgery, periproctitis abscess, post episiotomy condition

HARDNESS

Colour

SIZE

extra-soft soft	 light blue pink
medium-soft medium	 light green dark blue
medium-hard hard	 violet dark green

One length only

Reha

182 BORT. Benefit at your side.



PRODUCT-BENEFIT

- High quality gripping aid for picking up large and small objects (bottles, needles, etc.)
- Special rubber grips allow calibrated and targeted grasping of objects
- Locking: Grasped objects are held even after the grip lever is released
- Also available in a collapsible version: REF 118 100, length 75 cm (29.5 inches), when closed 45 cm (17.7 inches)

THE MED-BENEFIT

- Two possible settings:
 - Short leverage: normal effort with unrestricted gripping ability of the hand
 - · Long leverage: less effort with restricted gripping ability of the hand
- Collapsible version, ideal for travel or wheelchair use

INDICATIONS

E.g paraplegia, restricted mobility



PRODUCT-BENEFIT

- Standard aluminium gripper
- Special rubber pads with anatomical handle
- Length: 75 cm (29.5 inches)

THE MED-BENEFIT

 For gripping even the smallest objects (needles etc.)

Lengths

 75 cm (REF 117 170)
 50 cm (REF 118 150)
 75 cm (REF 118 170)
 90 cm (REF 118 190)
 Collapsible 75/45 cm (REF 118 100)

1 cm = 0.3937 inch



- Self-adhesive padding for forearm crutches
- Fabric surface, can be cut easily
- Can be trimmed to size, therefore personal use
- Size: 11 × 7,5 cm (4.3 × 3.0 inches)
- Pack: Pair

COLOURS: blue black

BORT Grip Pads for Forearm Crutches [REF] 190 820

PRODUCT-BENEFIT

- Covered inside with durable neoprene
- Outside covered with skin-friendly fabric
- Easy to apply thanks to Velcro fastening
- High dimensional stability
- Size: 10 × 10 cm (3.9 × 3.9 inches)
- Pack: Pair

COLOUR: Dlack

APPLICATION

Gripping deficit caused by abnormal exertion in the hand or elbow area

BORT Special Detergent REF 599 070

- PRODUCT-BENEFIT
- Specifically for knitted medical compression fabrics
- Top purity when hand-washing
- For easy everyday use

Contents: 250 ml

THE OTC-BENEFIT

- Already dissolves skin fat at low temperatures
- Biodegradable
- The transparent bottle shows the filling level

APPLICATIONS

Long-term satisfaction when wearing knitted compression fabric, environmentally friendly to wash, high productivity

SIZE

One length only

Reha Jetergent

Quality brands for individual all-round care. OTC Brands



Practical aids for daily support.



Stabilising joint supports for free work and time.



Stabilising joint supports for free time and sport.



The thermal programme for feeling at ease.



Gentle silicone and gel foot protection.



Protection for sensitive feet – with fabric cover.

Practical aids for daily support.





PRODUCT-BENEFIT

- Clearly organised 7 day planner
- Braille lettering
- Size: 16 × 3 × 2 cm (6.3 × 1.2 × 0.8 inches)

COLOURS: Ctransparent blue



PRODUCT-BENEFIT

- Sharp stainless steel blade

COLOURS: 🗆 transparent 🔳 blue

THE OTC-BENEFIT

- Safe, controlled taking of medicines
- Easy to open

THE OTC-BENEFIT

- Makes cutting tablets and medicines easy and safe
- Easy to clean

Stabilising joint supports for work and free time.



THE ACTIVECOLOR*-BENEFIT

- Preventive application in the case of expected overstraining
- Effective stabilisation
- 🕒 Made in Germany



PRODUCT-BENEFIT

- Wrist support with adjustable Velcro fastening
- Compression can be individually adjusted

COLOURS: white beige blue blue black

THE ACTIVECOLOR®-BENEFIT

 Preventive application in the case of expected overstraining

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of wrist



medium ____ up to 19 cm large ____ 19 - 21 cm





- Stabilisation thanks to compressive knitted fabric

- Can be used for work or free time


INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES	Circumference of wrist
H	small up to 15 cm medium 15 - 17 cm large 17 - 19 cm x-large 19 - 21 cm

right and left symmetrical 1 cm = 0.3937 inch SIZES

Circumference of elbow



Stabilising joint supports for work and free time.



THE ACTIVECOLOR®-BENEFIT


Preventive application in the case of expected overstraining
 Stabilisation thanks to compressive knitted fabric

INDICATIONS

Back pain, lumbago, strains

SIZES

Circumference of waist



1 cm = 0.3937 inch

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of thigh

+- /	small	up to 43 cm
17	medium	43 – 50 cm
11	large	50 – 56 cm
17	x-large	56 – 61 cm



ActiveColor® Knee Support [REF] 1440	ActiveColor® Ankle Support REF 1450
PRODUCT-BENEFIT	PRODUCT-BENEFIT
 Round knitted knee support for prevention of overstraining Supporting effect Non-slip fit Compression class II 	 Circular knitted ankle brace to support the ankle joint Anatomical form, reinforced everted heel Can be worn with sports and walking shoes
COLOURS: 📕 beige 📕 blue 📕 black	- Compression class II COLOURS: — beige — blue — black

Can be used for work or free time

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES		Circumferer	nce of knee
e-{{	small medium large x-large xx-large		up to 32 cm 32 – 37 cm 37 – 42 cm 42 – 47 cm 47 – 52 cm

right and left symmetrical 1 cm = 0.3937 inch SIZES

Circumference of ankle

xx-large 27-29 cm

right and left symmetrical 1 cm = 0.3937 inch **ActiveColor**[®]

Stabilising joint supports for free time and sport.



THE ACTIVECOLOR® SPORT-BENEFIT

- Regulates moisture with comfort fibre: SENSIL
- Support from certain muscle groups with strong use
- 📀 Made in Germany



SENSIL.

PRODUCT-BENEFIT

- Compression wrist support for easy stabilisation of the wrist
- Anatomically contoured knit
- Hemmed thumb opening
- Regulates moisture with comfort fibre: SENSIL[®] Aquarius

COLOUR: I black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Preventive application in the case of expected overstraining
 Stabilisation thanks to compressive
- Stabilisation thanks to compressive knitted fabric

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of wrist

1/3	small		up to	15 cm
1 28	medium		15 —	17 cm
	large		17 –	19 cm
	x-large	_	19 —	21 cm





- Light neoprene support to stabilise the metacarpophalangeal joint and saddle joint of the thumb
- Supportive in load situations
- Additional stabilising elements at the thumb
- Variably adjustable using a Velcro fastener

COLOUR: I black/green

ActiveColor® Sport Wrist Support REF 1410 SP



PRODUCT-BENEFIT

- Light neoprene support for wrist stabilisation
- Supportive in load situations
- Variably adjustable using a Velcro fastener

COLOUR: black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Doesn't influence the gripping function
- Full support for sports activity
- Can be worn in water 👌

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of wrist

SIZES

small/ medium ____13 – 17 cm large/ x-large ____17 – 21 cm

please indicate left or right 1 cm = 0.3937 inch One size only

please indicate left or right 1 cm = 0.3937 inch ActiveColor[®] Sport

Stabilising joint supports for free time and sport.



PRODUCT-BENEFIT

- Round knitted compression elbow support
- Elastic in longitudinal and lateral directions
- Compression class II
- Regulates moisture with comfort fibre: SENSIL[®] Aquarius

COLOUR: 📕 black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Preventive application in the case of expected overstraining
- Stabilisation thanks to compressive knitted fabric

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SI	Z	ES
----	---	----

Circumference of elbow

1 -4	small	up to 26 cm
	medium	26 – 28 cm
	large	28 – 30 cm
(L)	x-large	30 – 32 cm

right and left symmetrical 1 cm = 0.3937 inch



PRODUCT-BENEFIT

- Light neoprene brace to relieve the elbow's muscle attachments
- Brace pressure individually adjustable using circumferential tape

&

COLOUR: 📕 black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Full support for sports activity
- Can be worn in water 🔬

INDICATIONS

Prevention of overstraining, irritations due to fatigue, 'tennis/golf elbow'

SIZES		erence of forearm s) below elbow joint)
	Size 1 Size 2 Size 3 Size 4	23 – 25 cm 25 – 27 cm 27 – 30 cm 30 – 34 cm





- Active thigh support
- Anatomically flat knit
- Compression class II
- Length: 21 cm (8.3 inches)
- Regulates moisture with comfort fibre: SENSIL® Aquarius

COLOUR: black/green

ActiveColor® Sport **Calf Support REF** 1320 SP

SENSIL

PRODUCT-BENEFIT

- Compression calf support
- Anatomically contoured knit
- Compression class II
- Length: 23 cm (9.06 inches)
- Regulates moisture with comfort fibre: SENSIL® Aquarius

COLOUR: black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Preventive application in the case of expected overstraining
- Stabilisation thanks to compressive knitted fabric Can be used for free time or when practicing sports

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

INDICATIONS

SIZES

After injuries to the calf, for treating sports injuries, e. g. ruptures of the muscle fibre, muscle iniury

SIZES

Circumference of thigh

17	small medium large x-large		up to 43 – 50 – 56 –	50 56	cm cm
----	-------------------------------------	--	-------------------------------	----------	----------

right and left symmetrical 1 cm = 0.3937 inch

medium _ large x-large

small

right and left symmetrical 1 cm = 0.3937 inch


up to 33 cm _ 33-38 cm 38–43 cm

43-47 cm

Circumference of calf

ActiveColor Sport

Stabilising joint supports for free time and sport.



PRODUCT-BENEFIT

- Round knitted knee support for prevention of overstraining
- With additional silicone grip top in all sizes, for anti-slip fit
- Optimal adhesion and pressure distribution on the edge of the thigh
- Compression class II
- Regulates moisture with comfort fibre: SENSIL® Aquarius

COLOUR: I black/green

THE ACTIVECOLOUR® SPORT-BENEFIT

- Preventive application in the case of expected overstraining
 Stabilization theory of the state o
- Stabilisation thanks to compressive knitted fabric

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES

Circumference of knee

	small	up to 32 cm
0 - (medium	 32 – 37 cm
	large	 37 – 42 cm
1 I	x-large	 42 – 47 cm
	xx-large	47 – 52 cm





- Circular knitted ankle support to relive the ankle joint
- Anatomical form, reinforced everted heel
- Can be worn with sports and walking shoes
- Compression class II
- Regulates moisture with comfort fibre: SENSIL[®] Aquarius

ActiveColor® Sport Compression Stocking REF 1310 SP



SENSIL.

PRODUCT-BENEFIT

- Sports compression stocking consisting of open-pored knitted mesh fabric
- Good climate management and high breathability
- Even compression pressure increases blood circulation and venous performance under load
- Reduces the period for muscle regeneration
- Regulates moisture with comfort fibre: SENSIL[®] Aquarius
- Pack: Pair

COLOUR: I black/green

COLOUR: I black/green

- Can be used for free time or when practicing sports

INDICATIONS

Prevention of overstraining, strains, sprains, irritations due to fatigue

SIZES	Circumference of ankle	
b-measure	small up to 21 cm medium 21 - 23 cm large 23 - 25 cm x-large 25 - 27 cm xx-large 27 - 29 cm	

	Size II Size III Size IV Size V	-

SIZES

Circumference of calf



The thermal programme for feeling at ease.



THE CLIMACARE®-BENEFIT

		Dress Size		
SIZES	Women	Men	SIZES	Circum
small medium large x-large	8/10 12/14 16/18 20/22	40		small medium large x-large (beige, silver, anthracite)

mference of hip up to 98 cm _ 98 – 108 cm 108 – 118 cm _ 118 – 128 cm


1 cm = 0.3937 inch







- Joint warmer for knee and elbow (only size small)
- Elastic knit, no restriction of movement
- Material composition:
 40 % angora, 30 % polyamide,
 25 % pure new wool, 5 % spandex
- Pack: Pair
- COLOURS: beige white



PRODUCT-BENEFIT

- Warms cold feet
- Non-slip pads provide greater safety when worn around the house
- Material composition:
 67 % pure new wool, 25 % angora,
 8 % polyamide
- Durable
- Pack: Pair

COLOURS: 🗌 white 🔳 anthracite

- Can be worn unnoticed under clothes

Very comfortable to wear thanks to the material composition

SIZES	Circumference of joint
small	up to 32 cm
medium	32 – 38 cm
large	38 – 44 cm
x-large (beige)	44 – 50 cm

SIZES	Shoe size UK
x-small (anthracite) .	0 - 2
small	3 - 4.5
medium	5.5 - 6
large	over 7

right and left symmetrical 1 cm = 0.3937 inch 1 cm = 0.3937 inch

Gentle silicone and gel foot protection.



PRODUCT-BENEFIT

- Slim silicone toe separator
- Separates dislocated and tight toes
- Available in three different lengths
- Contents: 2 pieces, includes talcum powder
- *length = length toe separator



PRODUCT-BENEFIT

- Slim gel toe separator with a ring for safe fixing
- Separates dislocated and tight toes
- Available in two different lengths
- Contents: 2 pieces

*length = length toe separator



PRODUCT-BENEFIT

- Silicone ring for relieving pressure on the fingers and toes
- Contents: 2 pieces, includes talcum powder



PRODUCT-BENEFIT

- Cover made of gel with a ventilation hole
- For relieving pressure on the fingers and toes
- Contents: 2 pieces
- *Ø = outside diameter of cap

*Ø = inner diameter ring

THE PEDISOFT®-BENEFI1

- Very comfortable to wear thanks to soft and smooth Silicone & Gel
- Optimal fit to any foot contour thanks to the "memory effect"
- The material returns to its original shape immediately after wearing





- Silicone toe pad
- For toe relief
- Fixation with a toe sleeve
- Contents: 1 piece, includes talcum powder



PRODUCT-BENEFIT

- Silicone forefoot insole
- For pressure relief in the forefoot area
- It is inserted directly into the shoe
- Long-lasting, durable and non-slip
- Contents: 1 pair, includes talcum powder



PRODUCT-BENEFIT

- Forefoot gel protection
- For relief of the ball or arch of the foot
- Fixed to the toe with a gel ring
- Contents: 1 pair



PRODUCT-BENEFIT

- Silicone forefoot pad
- Pressure relief in the forefoot area
- As with a toe spreader, the forefoot pad is fixed to the toe
- Contents: 2 pieces, includes talcum powder

Hygienic and long-lasting Easy to clean with soap or to disinfect, if required

Gentle silicone and gel foot protection.



PRODUCT-BENEFIT

- Wide, ergonomic silicone toe spreader
- Separates dislocated and tight toes
- Available in two different lengths
- Contents: 2 pieces, includes talcum powder

*length = length toe spreader



PRODUCT-BENEFIT

- Wide, ergonomic gel toe spreader with a ring for safe fixing
- Separates dislocated and tight toes
- Available in two different lengths
- Contents: 2 pieces

*length = length toe spreader



PRODUCT-BENEFIT

- Double gel toe spreader for safe fixing
- Separates dislocated and tight toes
- Available in two different lengths
- Contents: 2 pieces

*length = length toe spreader



PediSoft® Toe Spreader

PRODUCT-BENEFIT

- Ergonomic gel toe spreader with toe ball protection
- Separates dislocated and tight toes
- Contents: 1 piece

THE PEDISOFT®-BENEFIT

- Very comfortable to wear thanks to soft and smooth Silicone & Gel
- Optimal fit to any foot contour thanks to the "memory effect"
- The material returns to its original shape immediately after wearing





- Silicone protection of the ball of the foot for protection against pressure sores, e.g. in shoes
- Contents: 1 piece, includes talcum powder



PRODUCT-BENEFIT

 Little toe silicone shield for protection on the little toe

licone & Ge

 Contents: 1 piece, includes talcum powder



PRODUCT-BENEFIT

- For shock absorption and relief in case of heel spur
- For complaints in the ankle, knee and hip joints, in the heel area and in the arch of the foot
- Height: 6 mm (0.2 inches)
- Contents: 1 pair
- Designed for frequent use as it is durable
- Hygienic and long-lasting
- Easy to clean with soap or to disinfect, if required

Protection for sensitive feet with fabric cover.



Size 0 (shoe size 4-6), Size 1 (shoe size 7-9), Size 2 (shoe size 10-11) please indicate left or right

PRODUCT-BENEELT

- Soft splint with pad for metatarsophalangeal joint pressure relief in case of hallux valgus
- Stepless setting of the correction pressure using the laterally positioned corrective tape, Con. 1 pc

THE PEDISOFT®-TEXLINE-BENEFIT

- Pads provide optimal pressure distribution

- Optimal fit to any foot contour thanks to the "memory effect"

small (up to shoe size 9), large (from shoe size 10)

PRODUCT-BENEEIT

- For the simultaneous relief of pressure in the forefoot and ball area
- Surrounds the ball of the foot and the metacarpophalangeal joint softly
- Contents: 1 pair





- Toe correction pad with adjustable, elastic toe straps
- Corrects dislocation of hammer toes
- Up to three toes can be corrected at the same time
- Contents: 1 piece



PRODUCT-BENEFIT

- Silicone toe pad
- For toe relief
- Fixation with a toe sleeve
- Contents: 1 piece



PRODUCT-BENEFIT

- Foam tube for padding pressure sores
- Can be individually cut to the fitting length with scissors
- Length: 12 cm (4.7 inches)
- Packaging unit: 3 or 2 (large size)



PRODUCT-BENEFIT

- Tube bandage with round gel pad inside
- Partial pressure relief on the fingers and toes
- Length: approx. 3,5 cm (1.38 inches)
- Contents: 2 pieces
- The material returns to its original shape immediately after wearing
- Very hygienic
- Delicate wash 30 °C (with the exception of the Hammer Toe Correction Pad and the Pressure Protection)

PediSoft[®] TexLine

Protection for sensitive feet with fabric cover.



PRODUCT-BENEFIT

- Tube bandage coated with gel, can be trimmed to size, individually adaptable
- For relieving pressure on the fingers and toes
- Length: 15 cm (5.9 inches)
- Contents: 1 piece



PRODUCT-BENEFIT

- Closed toe cap, fully coated with gel
- For relieving pressure on the fingers and toes
- Contents: 1 piece, includes talcum powder

THE PEDISOFT®-TEXLINE-BENEFIT

- Pads provide optimal pressure distribution
- Protection from additional pressure from outside
- Comfortable to wear thanks to the fabric cover
- Optimal fit to any foot contour thanks to the "memory effect"





BORT GmbH P.O. Box 1330 | D-71367 Weinstadt, Germany Am Schweizerbach 1 | D-71384 Weinstadt Phone +49 (7151) 99200-0 Fax +49 (7151) 99200-50 www.bort.com | medical@bort.com



BORT. Benefit at your side.