



BORT SELLAFIX[®]

THUMB ORTHOSES



THE NEW RANGE:
MADE OF THERMOPLASTIC MATERIAL
IN A UNIQUE PERFORATED DESIGN

Made in Germany

WHICH THUMB JOINT DO YOU NEED TO HOLD?

THE CARPAL, METACARPAL OR PHALANGEAL JOINT?
- SELLAFIX® GIVES YOU THE CHOICE

BORT SellaFix® thumb orthoses are now even more comfortable to wear.

Up to 520 perforations allow plenty of air to reach the covered areas of the skin.

The edges are smoother and now no longer have seams. This means that the new SellaFix® orthoses cannot break when being fitted or if inadvertently squeezed.



Thanks to the thermoplastic material the orthoses can easily be reshaped and individually adapted.

Individually positionable terry padding strips provide exceptional comfort.

- > Localised adjustments can be made using a hair dryer
- > The thumb piece can be adjusted if briefly immersed in hot water (60°)
- > Minor reshaping operations can be performed using scissors

Indications:

Ulnar ligament lesion, rhizo-arthritis, fractures and tendon irritation
BORT SellaFix® metacarpal joint orthoses are available in three versions.



BORT Sella Fix® K

Art. no. 112 080

For immobilization of
> the carpal joint
while allowing full freedom of movement at the wrist in the case of strains and sprains



BORT Sella Fix® N

Art. no. 112 040

For immobilization of
> the carpal joint and
> the metacarpal joint




BORT Sella Fix® P

Art. no. 112 900

For immobilization of
> the carpal joint,
> the metacarpal joint
and
> the phalangeal joint

Measuring point:

Width of hand at
metacarpal
knuckle joint.



Sizes Hand width in cm

x-small	up to 7,5 cm
small	7,5 – 8,0 cm
medium	8,0 – 8,5 cm
large	8,5 – 9,5 cm
x-large	9,5 – 10,5 cm
xx-large	10,5 – 11,5 cm

Please indicate left or right.



BORT GmbH

P.O. Box 1330 · D-71367 Weinstadt
Ziegeleistr. 39-43 · D-71384 Weinstadt
medical@bort.de · www.bort.de

QUALITY
MANAGEMENT
certified acc. to
medical products guidelines
93/42/EWG



QUALITY
MANAGEMENT
we are certified
regular voluntary
control acc. to ISO 9001:2000

