

# NOVELTIES 2023



**BORT StabiloBasic**  
**special width**  
for the CARE of  
larger body shapes



## Generation

Shoulder brace which is particularly easy to put on

## Knee braces

Now in Med and Generation models: innovative BORT Polar curve joint

## Knee supports in Double Plus sizes

For particularly large thigh circumferences

**BORT. Benefit at your side.**



NEW

# BORT Generation Shoulder Brace

Shoulder brace which is particularly easy to put on



REF 215 700

## PRODUCT-BENEFIT

- Shoulder/arm adduction brace **with pouch** for immobilisation of the shoulder joint
- Waist and shoulder strap for safe fixing of the shoulder and arm

COLOUR: ■ blue



## THE GENERATION-BENEFIT

- Simple fitting process for the shoulder and waist strap
- Positioning aid for the back prevents twisting of the straps
- Quick-release fastener for simple and safe attachment of the arm pocket to the strap system



Suspension lock system  
for easy, comfortable handling



Specially constructed  
spine brace  
for a simple putting on  
process for the shoulder  
and waist straps

## SIZES

Body size

Size 1	up to 160 cm
Size 2	160 – 180 cm
Size 3	over 180 cm

please indicate left or right

Available during the first quarter of 2023

Generation





NEW

# BORT StabiloBasic special width padded Back Support

For the care of larger body shapes




SPECIAL WIDTH



REF 104 670

## PRODUCT-BENEFIT

- Back support for muscular stabilisation of the lumbar spine
- Support with anatomically-contoured knit
- Viscoelastic sacral pad, which can be freely positioned
- Very flat construction, therefor hardly visible under clothing
- Height: ventral 17 cm (6.7 inches), dorsal 23 cm (9.1 inches)


COLOUR:  silver

## THE BASIC-BENEFIT

- Additional elastic straps to lift the abdomen for relief of the lumbar spine **for larger sizes**
- Specially constructed fastening system simplifies putting on
- Breathable fabric

## SIZES

Circumference of waist

	Size 1	110 – 120 cm
	Size 2	120 – 130 cm
	Size 3	130 – 140 cm
	Size 4	140 – 150 cm





**NEW**

# BORT Plus and Double Plus sizes for knee supports

For the precise treatment of unusual leg proportions

With the introduction of Double Plus sizes, BORT extends its offer for particularly large thigh circumferences.

The Plus sizes offer for knee supports was also extended.

SIZES KNEE	d-measure	f-measure
Size 1	26 – 29 cm	32 – 36 cm
Size 2	29 – 33 cm	36 – 40 cm
Size 2+	29 – 33 cm	40 – 44 cm
Size 3	33 – 37 cm	40 – 44 cm
Size 3+	33 – 37 cm	44 – 48 cm
Size 4	37 – 41 cm	44 – 48 cm
Size 4+	37 – 41 cm	48 – 52 cm
Size 5	41 – 45 cm	52 – 56 cm
Size 5+	41 – 45 cm	56 – 60 cm
<b>Size 5++</b>	41 – 45 cm	<b>60 – 64 cm</b>
Size 6	45 – 49 cm	60 – 64 cm
Size 6+	45 – 49 cm	64 – 68 cm
<b>Size 6++</b>	45 – 49 cm	<b>68 – 72 cm</b>
Size 7	49 – 53 cm	68 – 72 cm
<b>Size 7+</b>	49 – 53 cm	<b>72 – 76 cm</b>



**NEW:**

**Double Plus sizes**

Extra wide circumference for an even broader care range



**BORT StabiloGen® Eco**  
**with a Grip Top**

**REF** 114 520 H

- Additional Double Plus sizes 5++ and 6++
- Additional sizes 7 and 7+



**BORT StabiloGen® Eco**  
**with a Grip Top Short**

**REF** 114 520 HK

- Additional Double Plus sizes 5++ and 6++
- Additional sizes 7 and 7+

## New sizes



**BORT Manu Basic®**

**REF** 112 920

- Additional sizes xx-large and xx-large plus



**BORT LumboXpress**  
**Back Support**  
**with Double Fastening**

**REF** 112 600

- Additional size 8

NEW DESIGN

# BORT Polar curve joint for knee braces

Based on the physiological rolling and gliding movement of the knee joint

Innovative: the BORT Polar curve joint has a wandering pivot axis.

Therefore, it can imitate the rolling-sliding motion sequence of the human knee joint.

Joint restriction is made without the use of tools.



## Less stress

Reduction of the stress load on the crucial ligaments



### BORT StabiloPro® Knee Brace

REF 182 150

- Now with an adjustable Polar curve joint



### BORT Generation Knee Brace

REF 215 400

- Now with an adjustable Polar curve joint

## Further knee braces with the BORT Polar curve joint



### BORT GenuXpress

REF 100 700

- Now with an adjustable Polar curve joint
- **NEW** Additional Plus sizes small plus, medium plus, large plus, x-large plus



### BORT GenuXpress ACL

REF 100 800

- Optimised joint
- **NEW** Additional Plus sizes small plus, medium plus, large plus, x-large plus

# Product modifications

## Simplified function



### BORT EpiBasic Sport

REF 122 600 SP

NEW DESIGN

- Positioning aid in the antecubital fossa
- Flexible stretching zone in the elbow area

## New colour



### BORT ManuStabilPro Open Style

REF 182 500

- **NEW**  
new colour grey



### BORT Vario Back Support with Pad

REF 112 610

- **NEW**  
new colour grey

All information without guarantee, error and changes reserved. Status: January 2023.

**BORT. Benefit at your side.**

BORT GmbH P.O. Box 1330 | D-71367 Weinstadt | Am Schweizerbach 1 | D-71384 Weinstadt  
Phone +49 (0) 7151 99200-0 | Fax +49 (0) 7151 99200-50 | [www.bort.com](http://www.bort.com) | [medical@bort.com](mailto:medical@bort.com)

M2047/2022-10/005 EN | Rev. 01

