

ActiveColor® Wrist Support

REF 1410

Function: Wrist complaints can arise both from a one-off burden (e.g. a strain) or also as a result of constant overstraining (carpal tunnel syndrome). The ActiveColor® Wrist Support with its variable Velcro fastening provides support and offers individually adjustable compression.

Colours: white, beige, blue, black

Indications: Prevention of overstraining, strains, sprains, irritations due to fatigue


Sizes: Circumference of wrist
medium up to 19 cm
large 19–21 cm

1 cm = 0,3937 inch




STABILISING JOINT SUPPORTS FOR FREE TIME AND SPORT

For overstrained joints and muscles



ActiveColor® Sport Thumb Splint

REF 1370 SP Can be worn in water



Function: Complaints due to overstress and wear of the thumb joints can have an impact due to pain. The ActiveColor® Sport thumb splint provides support in stress situations.

Sizes: Circumference of wrist
small/medium 13–17 cm
large/x-large 17–21 cm

Colour:
black/green

ActiveColor® Sport Wrist Support

REF 1410 SP Can be worn in water



Function: Wrist complaints are unpleasant during sports activities and can restrict motion sequences. The ActiveColor® Sport wrist support protects the joint and provides stability.


One size only

Colour:
black/green

1 cm = 0,3937 inch

ActiveColor® Sport Epi Brace

REF 1460 SP Can be worn in water



Function: Lightweight neoprene brace to relieve the elbow muscle insertions. The brace pressure can be set individually via the circumferential strap, providing full support for sports activities.

Sizes: Circumference of forearm
Size 1 23–25 cm, **Size 2** 25–27 cm,
Size 3 27–30 cm, **Size 4** 30–34 cm

Colour:
black/green

ActiveColor® Sport Calf Support

REF 1320 SP



Function: The compressing calf support relieves the calf muscles in case of heavy strain. It provides high stability via compression and supports you ideally for movement.

Sizes: Circumference of calf
small up to 33 cm, **medium** 33–38 cm,
large 38–43 cm, **x-large** 43–47 cm

Colour:
black/green

1 cm = 0,3937 inch

ActiveColor® Sport Knee Support

REF 1440 SP



Function: Overstress and unequal load support on the knee joints is one of the most frequent impairments. The supporting, anti-slip ActiveColor® Sport knee support **with an additional silicone grip top** provides good protection against overburdening.

Colour: black/green

Sizes: Circumference of knee
small up to 32 cm
medium 32–37 cm
large 37–42 cm
x-large 42–47 cm
xx-large 47–52 cm

1 cm = 0,3937 inch

ActiveColor® Sport Compression Stocking

REF 1310 SP



Function: The ActiveColor® Sport Compression Stocking was developed specifically for runners. Its even compression pressure increases blood circulation and venous performance under load and reduces the muscle regeneration period.

Sizes: Circumference of calf
Size II 25–31 cm, **Size III** 31–38 cm,
Size IV 38–44 cm, **Size V** 44–50 cm

Colour:
black/green

1 cm = 0,3937 inch

THE ACTIVECOLOR® SPORT-BENEFIT:

- Regulates moisture with comfort fibre: **AQUARIUS**
- Support from certain muscle groups with strong use
- Made in Germany

BORT GmbH | Am Schweizerbach 1
D-71384 Weinstadt-Benzach, Germany
medical@bort.com | www.bort.com



M2039/2024-06/002 EN | Rev.01



STABILISING JOINT SUPPORTS FOR FREE TIME AND WORK

For overstrained joints and muscles



BORT has developed the ActiveColor® Supports specifically for the simple, non-complex self-therapy of general joint complaints. The ActiveColor® Support adapts ideally to natural movement without restrictions. Additionally, the compression facilitates the joint in its physiological motion sequence.

Overview of the advantages:

- The ActiveColor® Supports facilitate the joint without disturbing the motion sequence.
- The ActiveColor® Supports are available in modern colours and can be combined with any clothing.
- Each support is available with a breathable storage bag for use when out (except REF 1490).
- During movement, there is practically no risk when using the ActiveColor® Support.



The support ought not be worn during periods of rest, e.g. overnight. In case of permanent discomfort, you must consult a physician.



**ActiveColor®
Thumb/Hand Support**



REF 1470

Function: The ActiveColor® Thumb/Hand Support serves to stabilise the wrist. The thumb opening ensures that the support is kept in place without slipping, also if the hand and arm are used intensely.

Colours: beige, blue, black

Indications: Prevention of overstraining, strains, sprains, irritations due to fatigue

Sizes: Circumference of wrist
small up to 15 cm
medium 15–17 cm
large 17–19 cm
x-large 19–21 cm



REF 1470 SP

With **AQUARIUS** thread for extra moisture transport.

1 cm = 0,3937 inch



**ActiveColor®
Elbow Support**



REF 1420

Function: The ActiveColor® Elbow Support with a high compression strength has a preventive effect in case of serious overstrain of the joints, tendons and the arm muscles. It also makes sense for all sports types with intense use of the arms.

Colours: beige, blue, black

Indications: Prevention of overstraining, strains, sprains, irritations due to fatigue

Sizes: Circumference of elbow
small up to 26 cm
medium 26–28 cm
large 28–30 cm
x-large 30–32 cm



REF 1420 SP

With **AQUARIUS** thread for extra moisture transport.

1 cm = 0,3937 inch



**ActiveColor®
Ankle Support**




REF 1450

Function: The ankle joint is constantly subject to injury by twisting, skew gait or stumbling. The ActiveColor® Ankle Support ensures a strong supportive effect and can be worn in sports shoes as well as in outdoor footwear.

Colours: beige, blue, black

Indications: Prevention of overstraining, strains, sprains, irritations due to fatigue

Sizes: Circumference of ankle
small up to 21 cm
medium 21–23 cm
large 23–25 cm
x-large 25–27 cm
xx-large 27–29 cm



REF 1450 SP

With **AQUARIUS** thread for extra moisture transport.

1 cm = 0,3937 inch



**ActiveColor®
Thigh Support**



REF 1430

Function: The ActiveColor® Thigh Support has a preventive effect when the thigh muscles are overstrained. Also for treating sports injuries of the thigh.

Colours: beige, blue

Indications: Prevention of overstraining, strains, sprains, irritations due to fatigue

Sizes: Circumference of thigh
small up to 43 cm
medium 43–50 cm
large 50–56 cm
x-large 56–61 cm



REF 1430 SP

With **AQUARIUS** thread for extra moisture transport.

1 cm = 0,3937 inch



**ActiveColor®
Knee Support**



REF 1440

Function: The overstraining and unequal load distribution on the knees is one of the commonest ailments. The supporting, anti-slip ActiveColor® Knee Support provides good protection against overstraining.

Colours: beige, blue, black

Indications: Prevention of overstraining, strains, sprains, irritations due to fatigue

Sizes: Circumference of knee
small up to 32 cm
medium 32–37 cm
large 37–42 cm
x-large 42–47 cm
xx-large 47–52 cm



REF 1440 SP

With **AQUARIUS** thread for extra moisture transport.

1 cm = 0,3937 inch



**ActiveColor®
Back Support**



REF 1490

Function: The ActiveColor® Back Support facilitates the spine and prevents dislocations in case of heavy burdens. Two dorsally integrated splints ensure additional back stability. With variable Velcro fastening for individual compression regulation.

Colours: beige, black

Indications: Back pain, lumbago, strains

Sizes: Circumference of waist
small 65–80 cm
medium 80–95 cm
large 95–110 cm
x-large 110–125 cm



REF 1490 SP

With **AQUARIUS** thread for extra moisture transport.

1 cm = 0,3937 inch

