

# Novelties 2022



- ⊕ GenuXpress ACL with an adjustable polar curve joint
- ⊕ New, extra sizes for various ankle supports
- ⊕ Colour change for various hand and back products



## Med

In the midst  
of healthy  
life



Range extension:  
BORT Corset Sport



Function and design revised:  
the Asymmetric® and  
Asymmetric® Plus knee supports



Range extension  
ActiveColor® Sport:  
with Epi Brace und Knee  
Support with a Grip Top



**BORT. Benefit at your side.**

## BORT GenuXpress ACL

Frame knee orthosis by the 4-point principle for guidance and stabilisation of the knee joint

### BRACE FRAME

- Stable aluminium frame construction using a 4-point principle
- Low self-weight
- Adjustable using tooth setting

### HARNESSES

- Padded and with individually adjustable length

### CONDYLE PADS

- Gentle knee joint guidance

### SELF-ADJUSTING TIBIA FITTING

- For full-surface, pressure-free fitting

### Innovative: BORT polar curve joint

#### ADJUSTABLE POLAR CURVE JOINT


- based on the physiological rolling and gliding movement of the knee joint
- Tool-free setting of extension and flexion

### LIGHTWEIGHT AND SLIM STRUCTURE

- for ideal wearing comfort

### EASYCLICK SYSTEM

- Easy to handle when putting the brace on and off



**BORT GenuXpress ACL**  
REF 100 800

easyClick-System

Self-adjusting tibia-fitting

BORT Polar curve joint

Med

**NEW**

PRODUCT-BENEFIT

- Frame knee orthosis by the 4-point principle for guidance and stabilisation of the knee joint
- Limitation of extension: 0°, 5°, 10°, 15°, 20° and 25°
- Limitation of flexion: 0° (180°), 30°, 45°, 60°, 75° and 90°

COLOUR: ■ grey

THE MED-BENEFIT

- High frame strength combined with low self-weight thanks to ultra-modern material
- Very good adaptability of brace frame
- Lightweight and slim structure for ideal wearing comfort
- Self-adjusting tibia fitting
- Tool-free setting of extension and flexion

BORT INNOVATION:

- Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint

BORT Polar curve joint

**NEW**

Available during the first quarter of 2022

Med

# BORT StabiloPro® Knee Support Open Style

Now with a new polar curve joint



**BORT StabiloPro® Knee Support Open Style**  
REF 182 300

BORT Polar curve joint

Med

**NEW DESIGN**

PRODUCT-BENEFIT

- Soft brace with movement limitation for guidance and relief of the knee joint
- Medial and lateral stabilisation
- Adjustable extension: 0°, 10°, 15°, 20°
- Adjustable flexion: 0°, 20°, 30°, 45°, 60°, 75°, 90°
- Joint splints cold-mouldable
- Breathable, elastic material for optimal comfort

COLOUR: ■ grey

THE MED-BENEFIT

- With four circular Velcro straps for particularly secure support
- It can be fully opened so is very easy to put on
- Plus-size system resulting in very good fitting
- Tool-free setting of extension and flexion

BORT INNOVATION:

- Adjustable polar curve joint, based on the physiological rolling and gliding movement of the knee joint

BORT Polar curve joint

**NEW**

**BORT. Benefit at your side.**

# BORT Corset Sport

The sporty side of BORT



## BACK



### PRODUCT-BENEFIT

- Corset for straightening and relief of the thoracic spine
- Straightening according to the rucksack principle
- Really easy to put on and off

COLOUR: ■ black/green

### THE SPORT-BENEFIT

- Free shoulder blade movement due to the slim design
- Conforms perfectly to the body contours, especially in the critical armpit area
- Modern, sporty design (with a reflecting label)

# BORT Asymmetric® BORT Asymmetric® Plus

Function and design revised



## KNEE



- ⊕ Change to a new design
- ⊕ SOFTflex knitted fabric for more elasticity, stability and very comfortable to wear
- ⊕ Woven silicone bonding zone in the thigh area



- ⊕ Change to a new design

# PRODUCT MODIFICATIONS



## FUNCTION AND DESIGN REVISED

**BORT ManuCarpal® Combi**  
REF 102 400

NEW DESIGN

Med

NEW

- Functionally shaped relief pad consisting of memory foam
- Relief of the carpal tunnel and the sensitive scar region (condition after surgery)
- Comfortable to wear

Available as of March 2022

**BORT Generation Back Brace**  
REF 215 800

NEW DESIGN

Generation

- + Grip straps with an improved gripping function

**BORT Generation Back Support**  
REF 215 850

NEW DESIGN

Generation

- + Grip straps with an improved gripping function

**BORT select® ManuZip® volar**  
REF 112 220

NEW DESIGN

select®

- + Relocation of the zip to the back of the hand

**BORT Stabilo Basic Lady Back Support with Pad**  
REF 104 680

NEW DESIGN

Med

- + Change to new Eco design
- + Easier to put on using the ergonomic latch plate

**BORT Stabilo Basic Back Support with Pad**  
REF 104 690

NEW DESIGN

Med

- + Change to new Eco design
- + Easier to put on using the ergonomic latch plate



# PRODUCT MODIFICATIONS



## NEW DESIGN



⊕ Change to new Eco design



⊕ Change to new Eco design

## NEW COLOUR



⊕ New colour grey



⊕ New colour grey



⊕ New colour grey



⊕ New colour grey




⊕ New colour grey



⊕ New colour grey



⊕ New colour grey



⊕ New colour grey




# PRODUCT MODIFICATIONS




## NEW SIZES




⊕ New size xx-large



⊕ New size xx-large



⊕ New size xx-large



⊕ New size xx-large



⊕ New size 47



⊕ New size xxx-large (only silver)

# OTC

## PRODUCT INNOVATION



## PEDISOFT®



### PRODUCT-BENEFIT

- Soft splint with pad for metatarsophalangeal joint pressure relief in case of hallux valgus
- Stepless setting of the correction pressure using the laterally positioned corrective tape
- Contents: 1 piece

COLOUR:  beige

### THE PEDISOFT®-BENEFIT

- Low material thickness and individual adjustment options
- Can also be worn during the daytime and at night
- Comfortable to wear due to thin, soft material
- Simple and easy to put on
- Attractive, functional design


**BORT. Benefit at your side.**

# OTC

# PRODUCT INNOVATION and MODIFICATIONS



## ACTIVECOLOR® SPORT




### PRODUCT-BENEFIT

- Light neoprene brace to relieve the elbow's muscle attachments
- Brace pressure individually adjustable using circumferential tape

COLOUR: ■ black/green

### THE ACTIVECOLOR® SPORT-BENEFIT

- Full support for sports activity
- Can be worn in water



- ⊕ Round knitted knee support for prevention of overstraining
- ⊕ With additional silicone grip top in all sizes, for anti-slip fit
- ⊕ Optimal adhesion and pressure distribution on the edge of the thigh

All information without guarantee, error and changes reserved. Status: January 2022.

**BORT. Benefit at your side.**

BORT GmbH P.O. Box 1330 | D-71367 Weinstadt | Am Schweizerbach 1 | D-71384 Weinstadt  
Phone +49 (0) 7151 99200-0 | Fax +49 (0) 7151 99200-50 | www.bort.com | medical@bort.com

