

NOVELTIES 2024



coxa pro

BORT CoxaPro Hip Brace

for securing the function and mobilisation
of the hip joint



Knee supports

BORT StabiloGen® Sport with a grip top. Also available in Double Plus and Plus sizes.

Generation

New wrist support, particularly easy to handle

BORT. Benefit at your side.

NEW

BORT CoxaPro Hip Brace

With simple restriction of movement for bending and stretching



REF 201 400

PRODUCT-BENEFIT

- Hip joint brace for securing the function and/or mobilisation of the hip joint
- Simple restriction of movement for bending and stretching
- Tool-free extension and flexion setting 0°, 15°, 30°, 45°, 60°, 75°, 90°
- Trochanter pad for intermittent massage
- Joint splint in the hip region and thigh can be adapted using cold-shaping (aluminium splint)

COLOUR: ■ grey

THE MED-BENEFIT

- Good functional security for the hip joint via two wide, opposing restraints on the thigh
- Light, slim design for very comfortable wearing
- High-quality functional cotton knitted fabric, particularly Velcro-resistant



With a trochanter pad for intermittent massage in case of tractus iliotibialis syndrome




INDICATIONS

All indications for which functional security and/or mobilisation of the hip joint is necessary, such as state after hip joint luxation, loosening hip endoprotheses, medium to severe coxarthrosis, muscular imbalance of the lumbar pelvic hip region, e.g. tractus iliotibialis syndrome, piriformis syndrome


SIZES

Circumference of iliac crest

	small	_____	85 – 100 cm
	medium	_____	100 – 115 cm
	large	_____	115 – 130 cm
	x-large	_____	130 – 145 cm

please indicate left or right

Accessory: Piriformis Pad for intermittent trigger point massage



Piriformis Pad

REF 991 126

PRODUCT-BENEFIT

- For trigger point massage of the piriformis muscle
- Can be used with the BORT CoxaPro Hip Brace
- Pack: single item

SIZES

One size only



NEW

BORT Generation Wrist Support

Easy and safe to handle due to practical lacing



REF 215 000

PRODUCT-BENEFIT

- Padded wrist support with palm reinforcement
- Elastic basic material
- Palm splint, anatomically pre-shaped, thermoplastic

COLOUR: ■ blue

THE GENERATION-BENEFIT

- Practical lacing enables easy, safe handling
- The specially padded zone in the palm area enables use combine with a walking aid/wheeled walker
- Attractive, functional design

SIZES

Circumference of wrist



small	_____	13 – 16 cm
medium	_____	16 – 19 cm
large	_____	19 – 23 cm

please indicate left or right

Our very first
BORT Generation
Wrist Support



Generation



NEW

BORT Talo Stabil® Eco Plus

With Stabilo® Strap for increased safety and stability



REF 054 650

PRODUCT-BENEFIT

- Active ankle joint support **with Stabilo® Strap**, partially elastic, freely positionable
- **SOFTflex knitted fabric** for more elasticity, stability and very comfortable to wear
- Treatment-supporting positioning of the foot in pronation, supination or neutral position
- Two silicone pads for intermittent massage
- Compression class II

COLOUR:  silver

THE ECO-BENEFIT

- Soft-velvety surface structure
- Pressure-reduced edge in the ankle area

SIZES

Circumference of ankle



small	_____	up to 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm
xx-large	_____	27 – 29 cm
xxx-large	_____	29 – 31 cm

right and left symmetrical

Med

NEU

BORT activemed Ankle Support Plus


With Stabilo® Strap for increased safety and stability



REF 220 650

PRODUCT-BENEFIT

- Anatomically contoured ankle support **with Stabilo® Strap** partially elastic, freely positionable
- Treatment-supporting positioning of the foot in pronation, supination or neutral position
- Two silicone pads for intermittent massage
- Microfibre knitted fabric, soft grip
- Compression class II

COLOUR:  graphite

THE MED-BENEFIT

- Pressure-reduced edge in the ankle area
- Anatomically knitted heel with a good fit

SIZES

Circumference of ankle



small	_____	up to 21 cm
medium	_____	21 – 23 cm
large	_____	23 – 25 cm
x-large	_____	25 – 27 cm

right and left symmetrical

Med

NEW

BORT StabiloGen® Sport with Grip Top

Sporty knee support, also available in Plus and Double Plus sizes



[REF] 114 520 SP-H


PRODUCT-BENEFIT

- Knee support for compression of soft tissue with a viscoelastic pad for the stabilisation and relief of the knee joint
- **With additional silicone grip top in all sizes**
- **SOFTflex 3D knitted fabric** with a velvety soft surface structure for more elasticity, stability and very comfortable to wear
- One lateral spiral stay on each side
- Highly elastic, soft flexion region in the hollow of the knee
- Compression class II
- **Double plus size system resulting in very good fitting**

COLOUR: ■ black/green

THE SPORT-PLUS

- **COOLMAX®**: wicks moisture away from the body, for dryness and comfort
- Ideal fit for sports activities
- **Latex-free fabric**
- **Anatomically and functionally optimised pad form** for safe guidance of the patella during intense activity
- Modern sporty design (with a reflecting label)
- Support length in size 3: 35 cm

SIZES		d-measure	f-measure
	Size 1	26 – 29 cm	32 – 36 cm
	Size 2	29 – 33 cm	36 – 40 cm
	Size 2+	29 – 33 cm	40 – 44 cm
	Size 3	33 – 37 cm	40 – 44 cm
	Size 3+	33 – 37 cm	44 – 48 cm
	Size 4	37 – 41 cm	44 – 48 cm
	Size 4+	37 – 41 cm	48 – 52 cm
	Size 5	41 – 45 cm	52 – 56 cm
	Size 5+	41 – 45 cm	56 – 60 cm
	Size 5++	41 – 45 cm	60 – 64 cm
	Size 6	45 – 49 cm	60 – 64 cm
	Size 6+	45 – 49 cm	64 – 68 cm
	Size 6++	45 – 49 cm	68 – 72 cm
	Size 7	49 – 53 cm	68 – 72 cm
	Size 7+	49 – 53 cm	72 – 76 cm

right and left symmetrical



The support for champions: Use under the toughest conditions

Jello Krahmer, current German Wrestling Champion 2023 (Greek-Roman style), wears the BORT StabiloGen® Sport with a Grip Top.

As a sponsor of the successful ASV Schorndorf wrestlers, BORT equips the team with various medical devices. Heavyweight Jello Krahmer was also involved in the product development of the new StabiloGen® Sport with a Grip Top.

Currently, Jello Krahmer, who competes in the highest weight category (130 kg), wears the support in Plus size 5+.

**Sport Support
in Double Plus
sizes**

Product developer Roman Freihaut when administering the BORT StabiloGen® Sport with a grip top in the size 5+.

Sport

NEW DESIGN

BORT Talo Stabil® Eco

Ankle Support now also available in wider and Plus sizes



BORT TaloStabil® Eco

REF 054 600

- Additional size xxx-large
- New Plus sizes xx-large+ and xxx-large+



SIZES	Circumference of ankle	b1-measure**
Children (blue)	17 – 19 cm	
small	up to 21 cm	
medium	21 – 23 cm	
large	23 – 25 cm	
x-large	25 – 27 cm	
xx-large	27 – 29 cm	32 – 34 cm
xx-large plus*	27 – 29 cm	34 – 36 cm
xxx-large*	29 – 31 cm	36 – 38 cm
xxx-large plus*	29 – 31 cm	38 – 40 cm

right and left symmetrical

* only in colour silver

** the b1 measurement is required for differentiation of the Plus size

BORT OmoSAT

Shoulder-Arm Adduction Brace with a simplified function



BORT OmoSAT

REF 121 980

- Quick-release fastener for simple and safe attachment of the arm pocket to the strap system
- Simple fitting process for the shoulder and waist strap

SIZES	Fits height
Size 1	up to 160 cm
Size 2	160 – 180 cm
Size 3	over 180 cm

right and left symmetrical

Med

BORT AktiVen® EXTRA: Foot length+

Phlebology/circular knitted fabric: BORT AktiVen® EXTRA medical compression stockings with a closed toe, also available immediately with foot length+.



AktiVen®

All information without guarantee, error and changes reserved. Status: January 2024.

BORT. Benefit at your side.

BORT GmbH P.O. Box 1330 | D-71367 Weinstadt | Am Schweizerbach 1 | D-71384 Weinstadt
Phone +49 (0) 7151 99200-0 | Fax +49 (0) 7151 99200-50 | www.bort.com | medical@bort.com