

# BORT Plus sizes

For the precise treatment of  
unusual leg and arm proportions



- ⊕ Supports and braces in Plus and Double Plus sizes
- ⊕ For more pronounced thigh and lower arm circumferences
- ⊕ Better fit, better function, more comfortable to wear

Double Plus sizes ++

Plus sizes +

Standard

NEW:

**Double Plus sizes**

Extra wide circumference for  
an even broader care range

5++ 6++

2+ 3+ 4+ 5+ 6+ 7+



**BORT. Benefit at your side.**

# Thinking one step

The number of means of orthopaedic leg care with exceptional proportions is constantly increasing. In particular, conical leg shapes with more pronounced thigh circumferences are common occurrences in medical stores nowadays – a challenge when selecting the fitting aid.

Numerous BORT knee supports and braces are therefore available in Plus and Double Plus sizes in addition to standard sizes, i.e. in intermediate sizes with a particularly loose-fitting cut for the thigh dimensions (f-measure). The lower leg dimension (d-measure) remains unchanged here and is equivalent to the respective standard size.

We also have a wrist support with Plus sizes for the C2-measure (lower arm dimension) in our programme for patients with special lower arm circumferences.

**With BORT supports and braces in Plus and Double Plus sizes, we can cover a wider care range. Fit, function and wearing comfort are improved.**

## NEW: Knee supports in Double Plus sizes

**BORT StabiloGen® Eco with a Grip Top**  
REF 114 520 H



Med

For stabilisation and relief of the knee joint  
**Double Plus sizes: 5++, 6++**  
Plus sizes: 2+, 3+, 4+, 5+, 6+, 7+


**BORT StabiloGen® Eco with a Grip Top Short**  
REF 114 520 HK



Med

For stabilisation and relief of the knee joint  
**Double Plus sizes: 5++, 6++**  
Plus sizes: 2+, 3+, 4+, 5+, 6+, 7+

SIZES KNEE	d-measure	f-measure
Size 1	26 – 29 cm	32 – 36 cm
Size 2	29 – 33 cm	36 – 40 cm
Size 2+	29 – 33 cm	40 – 44 cm
Size 3	33 – 37 cm	40 – 44 cm
Size 3+	33 – 37 cm	44 – 48 cm
Size 4	37 – 41 cm	44 – 48 cm
Size 4+	37 – 41 cm	48 – 52 cm
Size 5	41 – 45 cm	52 – 56 cm
Size 5+	41 – 45 cm	56 – 60 cm
<b>Size 5++</b>	41 – 45 cm	<b>60 – 64 cm</b>
Size 6	45 – 49 cm	60 – 64 cm
Size 6+	45 – 49 cm	64 – 68 cm
<b>Size 6++</b>	45 – 49 cm	<b>68 – 72 cm</b>
Size 7	49 – 53 cm	68 – 72 cm
<b>Size 7+</b>	49 – 53 cm	<b>72 – 76 cm</b>





# ahead: BORT Plus and



Trust the original: as the first company, BORT developed a Plus size system for knee supports in 2006 and is still one of the leading suppliers in this area today.



## Supports and braces in Plus sizes

### BORT select® StabiloGen®

[REF] 114 420



select®

For stabilisation and relief  
of the knee joint

Plus sizes: 2+, 3+, 4+, 5+, 6+

### BORT select® GenuZip®

[REF] 114 390



select®

For stabilisation and relief  
of the knee joint

Plus sizes: 2+, 3+, 4+, 5+, 6+

### BORT StabiloGen® Sport

[REF] 114 520 SP



Sport

For stabilisation and relief  
of the knee joint

Plus sizes: 2+, 3+, 4+, 5+, 6+

### BORT StabiloGen® Eco

[REF] 114 520



Med

For stabilisation and relief  
of the knee joint

Plus sizes: 2+, 3+, 4+, 5+, 6+

### BORT activemed Knee Support

[REF] 220 400



Med

For stabilisation and relief  
of the knee joint

Plus sizes: 2+, 3+, 4+, 5+, 6+

### BORT Knee Support with Patella Recess

[REF] 114 150



Med

For stabilisation and relief  
of the knee joint

Plus sizes: 2+, 3+, 4+, 5+, 6+





# special width for larger body shapes

For up to  
**165 cm**  
girth



BORT is your partner for exceptional medical care. And that not only for Plus sizes. Thus, we develop our special width products extra for patients with a large girth or knee circumference and excel with special properties:

- ⊕ For girths of up to 165 cm and knee circumferences up to 55 cm
- ⊕ Very comfortable to put on and off
- ⊕ Optimised wearing properties for larger body shapes
- ⊕ Robust design
- ⊕ Made in Germany – with the accustomed BORT quality

## BORT special width products

**BORT Stabilo Basic Special Width Padded Back Support**  
REF 104 670

**SPECIAL WIDTH**

For the muscular stabilisation of the lumbar spine  
**For girths of up to 150 cm**

**BORT PostOban® Special Width Thorax-Abdominal-Support**  
REF 104 180

**SPECIAL WIDTH**

For the circular stabilisation of the thoracic and abdominal region  
**For girths of up to 165 cm**

**BORT VarioPlus Back Support Special Width**  
REF 112 690

**SPECIAL WIDTH**

For stabilisation of the lumbar spine  
**For girths of up to 165 cm**

**BORT Stabilo® Knee Support Special Width**  
REF 114 490

**SPECIAL WIDTH**

For guidance and stabilisation of the knee joint  
**For knee circumferences up to 55 cm**

# From standard to made-to-measure: BORT is your size specialist




We have maintained a trustworthy partnership with our specialist shop customers for more than 40 years and convince with our selection, quality and experience.

A full range of medical devices from the areas of orthopaedics, phlebology and oedema therapy with proven therapeutic benefits and a known good value for money makes BORT a versatile, performant and reliable partner.



The wide range of supports, braces, medical compression stockings and accessories is characterised both by its high quality and its wide range of sizes – in particular for the precise care of unusual leg and arm proportions as well as larger body shapes.



For the precisely fitting care of your patients, place your trust in the wide range of sizes from BORT!