

So that it's fitting for everyone

BORT is your size specialist



- ⊕ BORT Plus and Double Plus sizes for conically shaped legs and arms
- ⊕ Specially wide products for particularly large body shapes



BORT. Benefit at your side.

BORT offers the full wide range of sizes



BORT has been a high-powered and trustworthy partner for its specialist retailer customers for more than 40 years.

A full range of medical devices from the areas of orthopaedics, phlebology and oedema therapy convinces with well-known good value for money.

Supports, braces and medical compression stockings from BORT excel both with their top quality and wide range of sizes.

Trust in BORT – your size specialist – for the precisely fitting care of your patients!





PLUS AND DOUBLE PLUS SIZES

For the care of unusual leg and arm shapes and sizes, BORT has incorporated supports and braces in Plus and Double Plus sizes in its range.

SPECIAL WIDTH

SPECIAL WIDTH


For particularly larger body shapes with a substantial girth or knee circumference, BORT special width products are first choice.



STANDARD SIZES

The normal BORT full range already includes a wide range of sizes and is constantly extended.

Forward-thinking: BORT Original Plus and Double Plus sizes



The needs for orthopaedic care for unusual body proportions is on the increase – in particular for conical leg shapes with large thigh circumferences.

Several BORT knee supports and braces are therefore also available in Plus and Double Plus sizes, i.e. in intermediate sizes with a particularly loose-fitting cut for the thigh dimensions (f-measure). The lower leg dimension (d-measure) remains identical here and is equivalent to the respective standard size. Plus sizes are also available for hand and ankle supports.

With BORT supports and braces in Plus and Double Plus sizes, we can cover a wider care range. Fit, function and wearing comfort are improved considerably.

For thigh circumferences up to **76 cm**



SIZES KNEE	d-measure	f-measure
Size 1	26 – 29 cm	32 – 36 cm
Size 2	29 – 33 cm	36 – 40 cm
Size 2+	29 – 33 cm	40 – 44 cm
Size 3	33 – 37 cm	40 – 44 cm
Size 3+	33 – 37 cm	44 – 48 cm
Size 4	37 – 41 cm	44 – 48 cm
Size 4+	37 – 41 cm	48 – 52 cm
Size 5	41 – 45 cm	52 – 56 cm
Size 5+	41 – 45 cm	56 – 60 cm
Size 5++	41 – 45 cm	60 – 64 cm
Size 6	45 – 49 cm	60 – 64 cm
Size 6+	45 – 49 cm	64 – 68 cm
Size 6++	45 – 49 cm	68 – 72 cm
Size 7	49 – 53 cm	68 – 72 cm
Size 7+	49 – 53 cm	72 – 76 cm

BORT StabiloGen® Eco
with a Grip Top
[REF] 114 520 H



Med

For stabilisation and relief of the knee joint

Double Plus sizes: 5++, 6++
Plus sizes: 2+, 3+, 4+, 5+, 6+, 7+

BORT StabiloGen® Eco
with a Grip Top Short
[REF] 114 520 HK



Med

For stabilisation and relief of the knee joint

Double Plus sizes: 5++, 6++
Plus sizes: 2+, 3+, 4+, 5+, 6+, 7+

BORT StabiloGen®
with a Grip Top
[REF] 114 520 SP-H

NEW



Sport

For stabilisation and relief of the knee joint

Double Plus sizes: 5++, 6++
Plus sizes: 2+, 3+, 4+, 5+, 6+, 7+

Knee support in Double Plus sizes



BORT select® StabiloGen®
[REF] 114 420



select®

For stabilisation and relief of the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT select® GenuZip®
[REF] 114 390



select®

For stabilisation and relief of the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT StabiloGen® Sport
[REF] 114 520 SP



Sport

For stabilisation and relief of the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT Knee Support
with Patella Recess
[REF] 114 150



Med

For stabilisation and relief of the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT activemed
Knee Support
[REF] 220 400



Med

For stabilisation and relief of the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT Osgood-Schlatter
Knee Support
[REF] 114 510



Med

For local decompression of tibial tuberosity
Plus sizes: 2+, 3+, 4+, 5+

Supports and braces in Plus sizes

**BORT Fillawnt
with Patella Ring**
REF 210 160



Med

For stabilisation and relief of
the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT Asymmetric®
REF 114 700



Med

For facilitate movement of
the patella
Plus sizes: 2+, 3+, 4+, 5+, 6+

BORT Asymmetric® Plus
REF 114 900



Med

For facilitate movement of
the patella
Plus sizes: 2+, 3+, 4+, 5+, 6+

**BORT Knee Support
with Adjustable Hinge**
REF 114 460



Med

For guidance of
the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

**BORT Stabilo® Knee Support
with Hinge**
REF 150 140



Med

For guidance of
the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

**BORT StabiloPro®
Knee Brace Open Style**
REF 182 300



Med

For stabilisation and relief of
the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

**BORT StabiloPro®
Knee Brace**
REF 182 150



Med

For stabilisation and relief of
the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+

**BORT Generation
Knee Brace**
REF 215 400



Generation

For guidance and relief of
the knee joint
Plus sizes: 2+, 3+, 4+, 5+, 6+, 7+

BORT GenuXpress
REF 100 700



Med

For guidance and stabilisation of
the knee joint
**Plus sizes: small plus,
medium plus, large plus, x-large plus**

BORT GenuXpress ACL
REF 100 800



Med


For guidance and stabilisation of
the knee joint
**Plus sizes: small plus,
medium plus, large plus, x-large plus**

Ankle support in Plus sizes
for patients with a large lower leg circumference

**PLUS-SIZES
ANKLE JOINT**

	b-measure	b1-measure*
small	up to 21 cm	
medium	21 – 23 cm	
large	23 – 25 cm	
x-large	25 – 27 cm	
xx-large	27 – 29 cm	32 – 34 cm
xx-large plus	27 – 29 cm	34 – 36 cm
xxx-large	29 – 31 cm	36 – 38 cm
xxx-large plus	29 – 31 cm	38 – 40 cm

Measuring Points:
Ankle joint



* The b1 measure is required for differentiation of the Plus size

BORT TaloStabil® Eco
REF 054 600



Med


For compression of soft tissues
**Plus sizes:
xx-large plus, xxx-large plus**

Wrist support in Plus sizes
for patients with a large lower arm circumference

**PLUS-SIZES
WRIST**

	C1-measure	C2-measure
x-small	up to 13 cm	14 – 16 cm
small	13 – 15 cm	16 – 18 cm
small plus	13 – 15 cm	18 – 20 cm
medium	15 – 17 cm	18 – 20 cm
medium plus	15 – 17 cm	20 – 22 cm
large	17 – 19 cm	20 – 22 cm
large plus	17 – 19 cm	22 – 24 cm
x-large	19 – 21 cm	22 – 24 cm
x-large plus	19 – 21 cm	24 – 26 cm
xx-large	21 – 23 cm	24 – 26 cm
xx-large plus	21 – 23 cm	26 – 28 cm

Measuring Point:
Circumference
of wrist



BORT ManuBasic®
REF 112 920



Med

To stabilise the wrist
**Plus sizes (in C2-measure):
small plus, medium plus, large plus,
x-large plus, xx-large plus**

In 2006, BORT was the first company which developed a Plus size system for knee supports and is one of the leading suppliers in this area today.



Provides excellent support when really needed: BORT StabuloGen® Sport with a Grip Top

BORT Sport supports are developed specifically for use during sports activity and combine maximum functionality with dynamic design.

The new BORT StabuloGen® Sport Knee Support with a Grip Top stabilises and relieves the knee joint. The latex-free SOFTflex 3D knitted fabric provides elasticity and high wearing comfort. Its optimised pad shape guarantees safe kneecap guidance.

The support model is also available in Double Plus and Plus sizes for particularly large thigh circumferences.



Our very first sports support with a **Double-Plus size**

SILICONE GRIP TOP

- For secure fixing
- Also integrated in all sizes

COOLMAX® HIGH-TECH FIBRE

- Wicks moisture away from the body
- For pleasant and dry wearing comfort



FUNCTIONAL COMFORT ZONE IN THE POPLITEAL FOSSA

- Highly elastic and breathable
- For high wearing comfort

Support length in size 3:
35 cm



The support for champions

Jello Krahmer, current German Wrestling Champion 2023 (Greek-Roman style), wears the BORT StabiloGen® Sport with a Grip Top.



Product developer Roman Freihaut when administering the BORT StabiloGen® Sport with a grip top in the size 5+

As a sponsor of the successful ASV Schorndorf wrestlers, BORT equips the team with various medical devices. Heavyweight Jello Krahmer was also involved in the product development of the new StabiloGen® Sport Knee Support with a Grip Top.

Use under the toughest conditions

Currently, Jello Krahmer, who competes in the highest weight category (130 kg), wears the sport support in Plus size 5+. His conclusion: "A strong knee support which provides me with stability and fits perfectly."

Our offer for larger body shapes: BORT special width

Millions of Germans suffer from obesity – and it's on the increase. The abnormally high proliferation of body fat is also a challenge for the medical store involved.

BORT is the size specialist for exceptional medical care. We develop our special width products extra for patients with a large girth or knee circumference and they excel with special properties:

- ⊕ For girths of up to 180 cm
- ⊕ Very comfortable to put on and off
- ⊕ Optimised wearing properties for larger body shapes
- ⊕ Robust design
- ⊕ Made in Germany – with the accustomed BORT quality

SPECIAL WIDTH

For up to
180 cm
girth




BORT Stabilo Basic
special Width Padded Back Support
 REF 104 670



SPECIAL WIDTH
Med

For muscular stabilisation of the lumbar spine
For a girth up to 160 cm

BORT VarioPlus
Back Support special Width
 REF 112 690



SPECIAL WIDTH
Med

For stabilisation of the lumbar spine
For a girth up to 165 cm

BORT PostOban® Special Width
Thorax-Abdominal-Support
 REF 104 180



SPECIAL WIDTH
Med

For circular stabilisation of thoracic and abdominal regions
For a girth up to 180 cm

