

BORT AktiVen®

Medical compression stockings



Size table in inches

| | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | |
|----|-----------|--------------------|----------------|----------------|----------------|-----------------|------------------|------------------|------------------|
| | cT | Pantyhose plus | 29.5–41.3 | 31.5–45.3 | 37.4–49.2 | 39.4–51.2 | 43.3–55.1 | 45.3–57.1 | 47.2–61.0 |
| | | Pantyhose standard | 23.6–35.4 | 25.6–37.4 | 29.5–43.3 | 31.5–45.3 | 33.5–47.2 | 35.4–49.2 | 37.4–51.2 |
| | cH | Pantyhose plus | 39.4–47.2 | 41.3–49.2 | 43.3–53.2 | 45.3–55.1 | 47.2–57.1 | 49.2–59.1 | 53.2–63.0 |
| | | Pantyhose standard | 33.5–41.3 | 35.4–43.3 | 37.4–47.2 | 39.4–49.2 | 41.3–51.2 | 43.3–53.2 | 45.3–55.1 |
| | cG | AT tights | 16.5–22.4 | 17.7–24.4 | 19.3–26.4 | 20.9–28.4 | 22.1–30.3 | 23.6–32.2 | 24.8–33.5 |
| | | Border wide | 19.3–22.4 | 20.9–24.4 | 22.4–26.4 | 24.4–28.4 | 26.0–30.3 | 27.6–32.2 | 28.7–33.5 |
| | | Border normal | 16.5–19.3 | 17.7–20.9 | 19.3–22.4 | 20.9–24.4 | 22.1–26.0 | 23.6–27.6 | 24.8–28.7 |
| | cF | | 15.4–20.5 | 16.1–22.1 | 17.3–23.6 | 18.5–25.6 | 19.7–27.2 | 20.9–28.7 | 22.4–30.7 |
| | cE | | 11.8–14.6 | 13.0–15.8 | 13.8–16.9 | 14.6–18.1 | 15.4–18.9 | 16.1–20.1 | 16.9–21.3 |
| | cD | | 10.6–13.0 | 11.4–13.8 | 12.6–15.0 | 13.4–15.8 | 14.2–17.3 | 15.0–18.9 | 15.8–19.7 |
| | cC | | 11.0–13.4 | 11.8–14.6 | 13.4–15.8 | 14.2–16.9 | 15.0–18.1 | 15.8–18.9 | 16.5–20.1 |
| | cB1 | | 9.1–10.6 | 9.5–11.4 | 10.2–12.2 | 11.4–13.4 | 12.2–14.6 | 13.0–15.4 | 13.8–16.1 |
| | cB | | 7.1–7.9 | 7.9–8.7 | 8.7–9.5 | 9.5–10.2 | 10.2–11.0 | 11.0–11.8 | 11.8–12.6 |
| | cY | | 10.6–12.2 | 11.4–13.0 | 11.8–13.8 | 12.2–14.6 | 12.6–15.0 | 13.0–15.8 | 13.8–16.9 |
| cA | | 6.7–8.7 | 7.5–9.5 | 8.3–10.2 | 9.1–11.4 | 9.8–12.6 | 10.6–13.4 | 11.0–13.8 | |

| | | |
|------------------------------|-----------------|------------------|
| AD length (D size) | short 13.4–15.4 | normal 15.8–18.1 |
| AG/AT length (G size) | short 24.0–27.6 | normal 28.0–31.9 |

Note on foot length: The closed foot toe for the qualities VITAL, TREND and 100 can be worn in shoe sizes 4–11. AktiVen® EXTRA stockings with a closed toe are also available in "Foot length+".

We recommend: Leg length "normal" + foot lengths $\ell Z > 10.2$ " (as of shoe size 8.5) = **normal/foot length+**
 Leg length "short" + foot lengths $\ell Z > 9.4$ " (as of shoe size 6) = **short/foot length+**

BORT. Benefit at your side.

BORT AktiVen®

Medical compression stockings



Size table in cm

| | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | |
|----|-----------|--------------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| | cT | Pantyhose plus | 75–105 | 80–115 | 95–125 | 100–130 | 110–140 | 115–145 | 120–155 |
| | | Pantyhose standard | 60–90 | 65–95 | 75–110 | 80–115 | 85–120 | 90–125 | 95–130 |
| | cH | Pantyhose plus | 100–120 | 105–125 | 110–135 | 115–140 | 120–145 | 125–150 | 135–160 |
| | | Pantyhose standard | 85–105 | 90–110 | 95–120 | 100–125 | 105–130 | 110–135 | 115–140 |
| | cG | AT tights | 42–57 | 45–62 | 49–67 | 53–72 | 56–77 | 60–82 | 63–85 |
| | | Border wide | 49–57 | 53–62 | 57–67 | 62–72 | 66–77 | 70–82 | 73–85 |
| | | Border normal | 42–49 | 45–53 | 49–57 | 53–62 | 56–66 | 60–70 | 63–73 |
| | cF | | 39–52 | 41–56 | 44–60 | 47–65 | 50–69 | 53–73 | 57–78 |
| | cE | | 30–37 | 33–40 | 35–43 | 37–46 | 39–48 | 41–51 | 43–54 |
| | cD | | 27–33 | 29–35 | 32–38 | 34–40 | 36–44 | 38–48 | 40–50 |
| | cC | | 28–34 | 30–37 | 34–40 | 36–43 | 38–46 | 40–48 | 42–51 |
| | cB1 | | 23–27 | 24–29 | 26–31 | 29–34 | 31–37 | 33–39 | 35–41 |
| | cB | | 18–20 | 20–22 | 22–24 | 24–26 | 26–28 | 28–30 | 30–32 |
| | cY | | 27–31 | 29–33 | 30–35 | 31–37 | 32–38 | 33–40 | 35–43 |
| cA | | 17–22 | 19–24 | 21–26 | 23–29 | 25–32 | 27–34 | 28–35 | |

| | | |
|-------------------------------|-------------|--------------|
| AD length (ℓD size) | short 34–39 | normal 40–46 |
| AG/AT length (ℓG size) | short 61–70 | normal 71–81 |

Note on foot length: The closed foot toe for the qualities VITAL, TREND and 100 can be worn in shoe sizes 4–11. AktiVen® EXTRA stockings with a closed toe are also available in “**Foot length+**”.

We recommend: Leg length “normal” + foot lengths ℓZ > 26 cm (as of shoe size 8.5) = **normal/foot length+**
 Leg length “short” + foot lengths ℓZ > 24 cm (as of shoe size 6) = **short/foot length+**

BORT. Benefit at your side.

Structure of the AktiVen® order numbers



The BORT AktiVen® product order numbers are based on clear classification.

They consist of:

- The 3-digit **REF** number (defines quality, design and toe)
- The 5-digit variant number (defines compression class, border/pantyhose, colour, length and size)

Thus, an 8-digit order number results which clearly defines each compression stocking model.

Note: For AD stockings with a grip top, the 3-digit REF number is extended with an H (example: 511H).

| REF number | | | Variant number | | | | | | | | |
|-------------------|------------------|---------------|----------------|-----|-----------------------------------|-----------------------|------------------|-------------------|------------------------------|------------------|----------------|
| Quality (digit 1) | Design (digit 2) | Toe (digit 3) | CCL (digit 4) | | Border/Pantyhose (digit 5) | | Colour (digit 6) | | Length (digit 7) | | Size (digit 8) |
| 5 = 100 | 1 = AD | 1 = open | 1 = CCL 1 | 100 | 0 = Simple border | AD | B = Beige | VITAL TREND | 1 = normal | | 1 |
| 6 = VITAL | 2 = AG | 2 = closed | 2 = CCL 2 | | 1 = Nubbed silicone border 5.0 cm | AD/AG VITAL EXTRA 100 | K = Caramel | VITAL TREND EXTRA | 2 = short | | 2 |
| 7 = EXTRA | 3 = AT | | | | 2 = Nubbed silicone border wide | | M = Marine | VITAL TREND | 3 = normal/ Foot length+* | EXTRA closed toe | 3 |
| 8 = TREND | | | | | 3 = Lace silicone border 5.0 cm | AG TREND | A = Anthracite | | 4 = short/ Foot length+* | EXTRA closed toe | 4 |
| | | | | | 4 = Lace silicone border wide | | S = Black | | | | 5 |
| | | | | | 5 = Pantyhose standard | AT | N = Skin | 100 | | | 6 |
| | | | | | 6 = Pantyhose plus (AT/U) | TREND 100 | | | | | 7 |
| | | | | | 7 = Nubbed silicone border 2.5 cm | AD | | | | | |

Sample order number: 512-20 N 27 100, AD, closed toe, CCL 2, simple border, colour: skin, length short, size 7

*Foot length+ = Foot length: /Z > 26 cm (length normal) / /Z > 24 cm (length short)

| | | | | | | | |
|---|---|----|---|---|---|---|---|
| 5 | 1 | 2- | 2 | 0 | N | 2 | 7 |
|---|---|----|---|---|---|---|---|

REF number

Variant number

BORT. Benefit at your side.