

BORT AktiVen®

Medical compression stockings



Size table

VITAL	TREND	100		1	2	3	4	5	6	7
			Pantyhose plus	75–105	80–115	95–125	100–130	110–140	115–145	120–155
	cT		Pantyhose standard	60–90	65–95	75–110	80–115	85–120	90–125	95–130
			Pantyhose plus	100–120	105–125	110–135	115–140	120–145	125–150	135–160
	cH		Pantyhose standard	85–105	90–110	95–120	100–125	105–130	110–135	115–140
			AT tights	42–57	45–62	49–67	53–72	56–77	60–82	63–85
			Border wide	49–57	53–62	57–67	62–72	66–77	70–82	73–85
	cG		Border normal	42–49	45–53	49–57	53–62	56–66	60–70	63–73
			cF	39–52	41–56	44–60	47–65	50–69	53–73	57–78
			cE	30–37	33–40	35–43	37–46	39–48	41–51	43–54
			cD	27–33	29–35	32–38	34–40	36–44	38–48	40–50
			cC	28–34	30–37	34–40	36–43	38–46	40–48	42–51
			cB ¹	23–27	24–29	26–31	29–34	31–37	33–39	35–41
			cB	18–20	20–22	22–24	24–26	26–28	28–30	30–32
			cY	27–31	29–33	30–35	31–37	32–38	33–40	35–43
		cA	17–22	19–24	21–26	23–29	25–32	27–34	28–35	

AD length (cD size)	short 34–39	normal 40–46
AG/AT length (cG size)	short 61–70	normal 71–81

BORT. Benefit at your side.

Structure of the AktiVen® order numbers



The BORT AktiVen® product order numbers are based on clear classification.

They consist of:

- The 3-digit **REF** number (defines quality, design and toe)
- The 5-digit variant number (defines compression class, border/pantyhose, colour, length and size)

Thus, an 8-digit order number results which clearly defines each compression stocking model.

REF number			Variant number							
Quality (digit 1)	Design (digit 2)	Toe (digit 3)	CCL* (digit 4)		Border/Pantyhose (digit 5)		Colour (digit 6)		Length (digit 7)	Size (digit 8)
5 = 100	1 = AD	1 = open	1 = CCL 1	100	0 = Simple border		AD	B = Beige	1 = normal	1
6 = VITAL	2 = AG	2 = closed	2 = CCL 2	VITAL TREND 100	1 = Nubbed silicone border		AD/AG VITAL 100	K = Caramel		2 = short
8 = TREND	3 = AT				2 = Nubbed silicone border wide			AG TREND	M = Marine	
					3 = Lace silicone border		AT VITAL TREND 100		A = Anthracite	
					4 = Lace silicone border wide			VITAL, TREND 100	S = Black	
					5 = Pantyhose standard		AT/U TREND		N = Skin	
					6 = Pantyhose plus					

Sample order number: 512-20 N 27 100, AD, closed toe, CCL 2, simple border, colour: skin, length short, size 7

*CCL = compression class

5	1	2-	2	0	N	2	7
REF number			Variant number				

BORT. Benefit at your side.