

Little heroes need good support.



FOR KIDS!

Med

In the midst
of healthy
life.



BORT. Benefit at your side.


Best protection so that little joints can grow healthily.

FOR KIDS!

Playing, frolicking, jostling – children simply enjoy movement. As an experienced manufacturer of supports suitable for children, BORT products meet all the requirements which enable the small patients' joints to become healthy again:

They relieve, protect and support. With a high therapeutic effectiveness and in first-class quality. Developed by our experts, who have adjusted each support specifically to the anatomical characteristics and provision requirements for children.


No. 127 260 KIDS



BORT Cervical Support

Anatomic cervical support to stabilise the lumbar spine. Can be worn during the day and at night. Use, for example, in case of pain in the lumbar spine region.


No. 121 350 KIDS



BORT OmoBasic® acc. to Gilchrist

Open shoulder-arm adduction brace for immobilisation of the shoulder joint. Flexible setting on the hand and upper arm. Use, e.g. in case of sprains.


No. 121 700 KIDS



BORT Clavicle Support

Backpack support. Adapts itself perfectly to the body contours, especially in the critical underarm area. Use e.g. for treatment of clavicle fractures.


No. 104 600 KIDS



BORT StabiloFix®

Corset for correcting the posture of the thoracic spine. Fastening in the chest area makes it possible to put on/remove the item independently. Use during, e.g. misalignment of thoracic spine.


No. 112 770 KIDS



BORT DigiSoft® Finger Brace

Finger brace for immobilisation of interphalangeal joints. Immobilisation of either one or two fingers. Use e.g. for positioning.

No. 112 050 KIDS



BORT SellaTex® Classic

Wrist brace with additional immobilisation of the metacarpophalangeal joint and saddle joint of the thumb. Adjustable thumb element. Slip-on aid. Use e.g. post-traumatically and post-operatively.

No. 112 730 KIDS

New
edition



Med

BORT ManuStabil® short

Rigid wrist brace with palm and dorsal reinforcement splints, individually adjustable. Can also be worn at night. Use e.g. in case of tenosynovitis.

No. 114 510 KIDS



Med

BORT Osgood-Schlatter Knee Support

Knee support for localised pressure relief of tibial tuberosity. Multi-zone silicone pad. Use e.g. in case of Osgood-Schlatter disease.

No. 114 520 KIDS



Med

BORT StabiloGen® Eco

Classical knee support for compression of soft tissues with pad for stabilisation and relief of the knee joint. Velvety soft surface structure. Use e.g. in case of ligament instability.

No. 150 140 KIDS



Med

BORT Stabilo® Knee Support with Articulated Joint

Knee support with side joint splints for guidance of the knee joint. Extension limit at 0°. Use e.g. in case of mild to moderate lateral ligament instability of the knee joint.

No. 145 000 KIDS



Med

BORT Immobilisation Splint with Patella Recess

Knee brace for immobilisation of the knee joint with a 0°–20° flexion. The support adapts to body shape thanks to 3-piece design. Use e.g. post-traumatically and post-operatively.

No. 054 600 KIDS



Med

BORT TaloStabil® Eco

Active ankle support for compression of soft tissue. Pressure relief at the base of the ankle. Use, e.g. in case of sprains and ligament instability.

No. 053 200 KIDS



Med

BORT Helix S Spiraldynamik® Lower Ankle Support

Ankle support for supporting torsion of the back of the foot and raising the arch of the foot. Stabilisation of the ankle and activation of foot muscles through compression. Use e.g. in case of ligament instability and defects in the shape of the foot due to dysfunction of the lower ankle joint.

No. 054 800 KIDS



Med

BORT TaloFX®

Soft ankle brace for stabilisation of the lower and upper ankle joints. Limitation of talus prolapse with optimal stabilisation. Use e.g. in case of conservative therapy for ankle joint sprains and ligament injuries.

No. 100 500 KIDS



Med

BORT MalleoStabil® Brace

Combined air/gel brace and foam pad to stabilise the ankle, while limiting pronation and supination. No impairment when walking. Use, e.g. in case of sprains of all degrees of severity, functional therapy of ligament injuries.

No. 930 200 KIDS



Med

BORT Heel Lift made of Silicone

Silicone heel lift to compensate differences in leg-length. Fits any kind of shoe. Use e.g. in case of compensation of differences in leg-length.

The KIDS-Benefit from BORT



- ⊕ Developed specially for children's anatomy
- ⊕ Comfortable to wear whilst providing maximum freedom of movement
- ⊕ Easy to put on, encourages daily wear
- ⊕ Colour design



Your competent partner for the provision
of care with the BORT KIDS range

(your dealer stamp)

Med



M2068|2019-11|002 EN



BORT. Benefit at your side.

BORT GmbH P.O. Box 1330 | D-71367 Weinstadt, Germany | Phone +49 (0) 7151 99200-0
Fax +49 (0) 7151 99200-50 | www.bort.com | medical@bort.com